



HD 2105.1, HD 2105.2 TEMPERATURE-pH METERS

The **HD2105.1** and **HD2105.2** are portable instruments with a large LCD display. They measure pH and redox potential (ORP) in mV. They measure temperature by using Pt100 or Pt1000 immersion, penetration or contact probes.

The electrode calibration can be carried out on one, two or three points and the calibration sequence can be chosen from a list of 13 buffers.

The temperature probes are equipped with an automatic recognition module and factory calibration data are stored inside.

The HD2105.2 is a **datalogger**; it stores up to 34,000 pH and temperature samples which can be transferred to a PC from the instrument connected via the RS232C and USB 2.0 serial ports. Storing interval, printing and baud rate can be configured by using the menu.

Both models are equipped with an RS232C serial port and can transfer the acquired measurements in real time to a PC or to a portable printer.

The **Max**, **Min** and **Avg** function calculate the maximum, minimum or average values. Other functions include: the relative measurement REL, the Auto-HOLD function, and the automatic turning off that can also be excluded.

The instruments have IP66 protection degree.

INSTRUMENT TECHNICAL CHARACTERISTICS

Measured quantities: pH, mV, °C, °F

Instrument

Dimensions (Length x Width x Height)	185x90x40mm
Weight	470g (complete with batteries)
Materials	ABS, rubber
Display	2x4½ digits plus symbols Visible area: 52x42mm

Operating conditions

Operating temperature	-5...50°C
Storage temperature	-25...65°C
Working relative humidity	0...90%RH without condensation

Protection degree

IP66

Power

Batteries	4 1.5V type AA batteries
Autonomy	200 hours with 1800mAh alkaline batteries
Power absorbed with instrument off	20µA
Mains-supply unit (SWD10)	Output mains adapter 12Vdc/1000mA

Security of stored data

Unlimited, independent of the state of charge of the batteries

Time

Date and time	Real time schedule
Accuracy	1min/month max drift

Measured values storage - model **HD2105.2**

Type	2000 pages containing 17 samples each
Quantity	34000 pairs of measurements composed of (pH or mV) and (°C or °F)
Storage interval	1s, 5s, 10s, 15s, 30s, 1min, 2min, 5min, 10min, 15min, 20min, 30min and 1h.

RS232C serial interface

Type	RS232C electrically isolated
Baud rate	It can be set from 1200 to 38400 bauds
Data bit	8
Parity	None
Stop bit	1
Flow Control	Xon/Xoff
Serial cable length	Max 15m
Print interval	Immediate or selectable between: 1s, 5s, 10s, 15s, 30s, 1min, 2min, 5min, 10min, 15min, 20min, 30min and 1h.

USB interface - model **HD2105.2**

Type	1.1 - 2.0 electrically isolated
------	---------------------------------

Connections

Input module for the temperature probes	8-pole male DIN45326 connector
pH/mV input	Female BNC
Serial interface RS232C	8-pole MiniDin connector
USB interface	MiniUSB B-type connector
Mains adapter	2-pole connector (positive at centre)



WA-2



HD2110CSNM



CP23

pH measurement by the instrument

Measurement range	-2.000...+19.999pH
Resolution	0.01 or 0.001pH selectable from menu
Accuracy	± 0.001pH ±1 digit
Input impedance	> 10 ¹² Ω
Calibration error @25°C	Offset > 20mV Slope < 50mV/pH or Slope < 63mV/pH Sensitivity < 85% or Sensitivity < 106.5%
Temperature compensation automatic/manual	-50...+150°C

mV measurement by the instrument

Measurement range	-1999.9...+1,999.9mV
Resolution	0.1mV
Accuracy	±0.1mV ±1 digit
Drift after 1 year	0.5mV/year

temperature measurement by the instrument

Pt100 measurement range	-200...+650°C
Pt1000 measurement range	-200...+650°C
Resolution	0.1°C
Accuracy	±0.1°C ±1 digit
Drift after 1 year	0.1°C/year

TECHNICAL DATA OF PROBES AND MODULES EQUIPPED WITH INSTRUMENT Temperature probes Pt100 sensor with SICRAM module

Model	Type	Application field	Accuracy
TP472I	Immersion	-196°C...+500°C	±0.25°C (-196°C...+300°C) ±0.5°C (+300°C...+500°C)
TP472I.0 1/3 DIN Thin Film	Immersion	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP473PI	Penetration	-50°C...+400°C	±0.25°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP473P.0 1/3 DIN Thin Film	Penetration	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP474C.I	Contact	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP474C.0 1/3 DIN Thin Film	Contact	-50°C...+300°C	±0.3°C (-50°C...+300°C)
TP475A.0 1/3 DIN Thin Film	Air	-50°C...+250°C	±0.3°C (-50°C...+250°C)
TP472I.5	Penetration	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP472I.10	Penetration	-50°C...+400°C	±0.30°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP49A.0 Class A Thin Film	Immersion	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AC.0 Class A Thin Film	Contact	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AP.0 Class A Thin Film	Penetration	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP875.I	Globe-thermometer Ø150mm	-30°C...+120°C	±0.25°C
TP876.I	Globe-thermometer Ø50mm	-30°C...+120°C	±0.25°C
TP87.0 1/3 DIN Thin Film	Immersion	-50°C...+200°C	±0.25°C
TP878.0 1/3 DIN Thin Film TP878.1.0 1/3 DIN Thin Film	Photovoltaic	+4°C...+85°C	±0.25°C
TP879.0 1/3 DIN Thin Film	Compost	-20°C...+120°C	±0.25°C

Common characteristics

Temperature drift @ 20°C	0.003%/°C
--------------------------	-----------

4 wires Pt100 and 2 wires Pt1000 Probes

Model	Type	Application field	Accuracy
TP47.100.0 1/3 DIN Thin Film	4 wires Pt100	-50...+250°C	1/3 DIN
TP47.1000.0 1/3 DIN Thin Film	2 wires Pt1000	-50...+250°C	1/3 DIN
TP87.100.0 1/3 DIN Thin Film	4 wires Pt100	-50...+200°C	1/3 DIN
TP87.1000.0 1/3 DIN Thin Film	2 wires Pt1000	-50...+200°C	1/3 DIN

Common features

Temperature drift @20°C

Pt100 0.003%/°C


Pt1000 0.005%/°C

A For the models of portable data logger series **HD21XX.2** has been implemented with a new serial port miniUSB type HID (Human Interface Device).

When making the connection to the PC by the USB cable Type A - Mini USB B-type coded CP23, **no USB driver installation is requested.**

B For the connection of the models **HD21XX.1** to the RS232 port of your PC, the USB/serial converter is available (**code C.206**). The converter is equipped with its own drivers that have to be installed before connecting the converter to the PC (please see the details in the CDRom supplied with the converter).

C The port with the MiniDIN connector which is present on every model is an RS232C type. By means of the cable coded HD2110CSNM, an RS232 port of a PC or the HD40.1. printer can be connected.



ORDER CODES

HD2105.1: The kit is supplied with: instrument HD2105.1, 4 1.5V alkaline batteries, operating manual, case and DeltaLog9 software

HD2105.2: The kit is supplied with: instrument **data logger** HD2105.1, 4 1.5V alkaline batteries, operating manual, case and DeltaLog9 software

Electrodes, temperature probes, calibration solutions, cable for the data download (from PC or printer) have to ordered separately.

HD2110CSNM: 8-pole connection cable MiniDin - Sub D 9-pole female for RS232C.

C.206: Cable for instruments of the series HD21...1 for direct connection to the USB input of a PC.

CP23: USB 2.0 connection cable type A - MiniUSB type B.

DeltaLog9: Software for download and management of the data on PC using Windows operating systems.

SWD10: Stabilized power supply at 100-240Vac/12Vdc-1A mains voltage.

HD40.1: The kit includes: 24-column portable thermal printer, serial interface RS232, 57mm paper width, four NiMh 1.2V rechargeable batteries, SWD10 power supply, instruction manual, 5 thermal paper rolls. It uses the optional cable HD2110 CSNM.

BAT-40: Spare battery pack for HD40.1 printer with built-in temperature sensor.

RCT: The kit includes 4 thermal paper rolls 57mm wide and 32mm diameter.

HD22.2: Laboratory electrode holder composed of base plate with built-in magnetic stirrer, shaft and replaceable electrode holder. Suitable diameter 12mm. Powered by power supplier SWD10 (optional).

HD22.3: Laboratory electrode holder composed of base plate. Flexible arm for free positioning. Suitable for electrodes with diameter 12mm.

pH Electrodes

KP 20: Gel pH combined electrode for general use, with S7 screw connector, EPOXY body.

KP 30: Gel pH combined electrode for general use, 1m cable with BNC, EPOXY body.

KP 50: Gel pH combined electrode, porous Teflon ring junction, suitable for emulsions, demineralised water and waste water with S7 screw connector, glass body.

KP 61: 3 diaphragm liquid filled pH combined electrode for wine, milk, cream, etc., S7 screw connector, liquid reference filling, glass body.

KP 62: 1 diaphragm gel pH combined electrode for general use, pure water, varnishes, gel filled, S7 screw connector, glass body.

KP 63: 1 liquid filled pH combined electrode for general use, varnishes, 1m cable with BNC, glass body.

KP 64: Liquid filled pH combined electrode, Teflon ring diaphragm, for wine, varnishes, emulsions, S7 screw connector, glass body.

KP 70: Pointed gel combined pH microelectrode diam. 6 x L=70 mm., with S7 screw connector, EPOXY body, glass tip, open junction for meat and cheese.

KP 80: Pointed gel pH combined electrode, with S7 screw connector, glass body, for cream, milk, viscous material, open junction.

KP100: Flat membrane gel combined pH electrode with S7 screw connector, glass body, for skin, leather, paper.

Characteristics and dimensions of the probes on page WA-76.

CP: 1.5m extension cable with BNC/S7 connector for electrode without cable, thread S7.

CP 5: 5m extension cable with BNC/S7 connector for electrode without cable, thread S7.

CP 10: 10m extension cable with BNC/S7 connector for electrode without cable, thread S7.

CP 15: 15m extension cable with BNC/S7 connector for electrode without cable, thread S7.

CE: S7 screw connector for pH electrode.

BNC: female BNC for extension cable

ORP Electrodes

KP 90: Redox Platinum liquid filled electrode with S7 screw connector, glass body.

KP 91: Gel Redox Platinum electrode, 1m cable with BNC, EPOXY body for general purpose light duty.

Characteristics and dimensions of the probes on page WA-76.

pH Buffer solutions

HD8642: Buffer solution 4.01pH - 200cc.

HD8672: Buffer solution 6.86pH - 200cc.

HD8692: Buffer solution 9.18pH - 200cc.

Redox Buffer solutions

HDR220: Redox buffer solution 220mV 0.5 l.

HDR468: Redox buffer solution 468mV 0.5 l.

Electrolyte solutions

KCL3M Ready to use solution for electrode refilling – 100 cc

Cleaning and maintenance

HD62PT: Diaphragm cleaning (tiourea in HCl) - 500ml.

HD62PP: Protein cleaning (pepsin in HCl) - 500ml.

HD62RF: Regeneration (fluorhydric acid) - 100ml.

HD62SC: Solution for electrode preservation - 500ml.

Temperature probes complete with SICRAM module

TP472I: Wire wound Pt100 sensor, immersion probe. Stem Ø 3 mm, length 300 mm. Cable length 2 m.

TP472I.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 230 mm. Cable length 2 m.

TP473PI: Wire wound Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP473P0: Thin film Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP474C.I: Wire wound Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP474C.0: Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP475A.0: Thin film Pt100 sensor, air probe. Stem Ø 4mm, length 230mm. Cable length 2 m.

TP472I.5: Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 500 mm. Cable length 2 m.

TP472I.10: Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 1000mm. Cable length 2 m.

TP49A.0: Thin film Pt100 sensor, immersion probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AC.0: Thin film Pt100 sensor, contact probe. Stem Ø 4 mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AP.0: Thin film Pt100 sensor, penetration probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP875.I: Wire wound Pt100 sensor, 150mm diameter globe-thermometer equipped with handle and SICRAM module. Cable length 2 m.

TP876.I: Wire wound Pt100 sensor, 50mm diameter globe-thermometer equipped with handle and SICRAM module. Cable length 2 m.

TP87.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 70 mm. Cable length 2 m.

TP878.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 2 m.

TP878.1.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 5 m.

TP879.0: Thin film Pt100 sensor, penetration probe for compost. Stem Ø 8 mm, length 1000 mm. Cable length 2 m.

Temperature probes without SICRAM module

TP47.100.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP47.1000.0: Thin film Pt1000 sensor, immersion probe. Probe's stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP47: Connector for Pt100 4-wire and Pt1000 2-wire probes without SICRAM module.

TP87.100.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 70mm. 4-wires connection cable with connector, length 1 m.

TP87.1000.0: Thin film Pt1000 sensor, immersion probe. Stem Ø 3mm, length 70mm. 2-wires connection cable with connector, length 1 m.



HD22.3



HD 2305.0 PORTABLE pH METER

The **HD2305.0** is a portable instrument with a large LCD display. It measures the pH and the redox potential (ORP) in mV. It measures the temperature by using Pt100 or Pt1000 immersion, penetration, contact or air probes. The electrode calibration can be carried out on one, two or three points at 4.01pH, 6.86pH and 9.18pH. The temperature probes are equipped with the automatic recognition SICRAM module which contains the factory calibration data stored inside.

The *Max*, *Min* and *Avg* function calculate the maximum, minimum or average values. Other functions include: the relative measurement REL, the Auto-HOLD function, and the automatic turning off that can also be excluded.

The instruments have IP67 protection degree.

INSTRUMENT TECHNICAL CHARACTERISTICS

Measured quantities: pH, mV, °C, °F

Instrument

Dimensions (Length x Width x Height)	140x88x38mm
Weight	160g (complete with batteries)
Materials	ABS
Display	2x4½ digits plus symbols Visible area: 52x42mm

Working conditions

Operating temperature	-5...50°C
Storage temperature	-25...65°C
Working relative humidity	0...90%RH without condensation

Protection degree **IP67**

Power

Batteries	3 1.5V type AA batteries
Autonomy	200 hours with 1800mAh alkaline batteries
Power absorbed with instrument off	<20µA

Connections

Input module for
temperature probes
pH/mV input

8-pole male DIN45326 connector
Female BNC

Measurement of pH by Instrument

Measurement range	-2.00...+19.99pH
Resolution	0.01
Accuracy	±0.01pH±1 digit
Input impedance	>10 ¹² Ω
Calibration error @25°C	Offset: >20mV Slope < 50mV/pH or Slope > 63mV/pH Sensitivity < 85% or Sensitivity > 106.5%
Temperature compensation automatic/manual	-50...+150°C

Measurement of mV by Instrument

Measurement range	-1999.9...+1999.9mV
Resolution	0.1mV
Accuracy	±0.1mV±1 digit
Drift after 1 year	0.5mV/year

Measurement of temperature by Instrument

Pt100 measurement range	-200...+650°C
Pt1000 measurement range	-200...+650°C
Resolution	0.1°C
Accuracy	±0.1°C±1 digit
Drift after 1 year	0.1°C/year

Temperature probes Pt100 sensor using SICRAM module

Model	Type	Application field	Accuracy
TP472I	Immersion	-196°C...+500°C	±0.25°C (-196°C...+300°C) ±0.5°C (+300°C...+500°C)
TP472I.0 1/3 DIN Thin Film	Immersion	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP473P.I	Penetration	-50°C...+400°C	±0.25°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP473P.0 1/3 DIN Thin Film	Penetration	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP474C.I	Contact	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP474C.0 1/3 DIN Thin Film	Contact	-50°C...+300°C	±0.3°C (-50°C...+300°C)
TP475A.0 1/3 DIN Thin Film	Air	-50°C...+250°C	±0.3°C (-50°C...+250°C)
TP472I.5	Penetration	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP472I.10	Penetration	-50°C...+400°C	±0.30°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP49A.0 Class A Thin Film	Immersion	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AC.0 Class A Thin Film	Contact	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AP.0 Class A Thin Film	Penetration	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP875.I	Globe-thermometer Ø150mm	-30°C...+120°C	±0.25°C
TP876.I	Globe-thermometer Ø50mm	-30°C...+120°C	±0.25°C
TP87.0 1/3 DIN Thin Film	Immersion	-50°C...+200°C	±0.25°C
TP878.0 1/3 DIN Thin Film TP878.1.0 1/3 DIN Thin Film	Photovoltaic	+4°C...+85°C	±0.25°C
TP879.0 1/3 DIN Thin Film	Compost	-20°C...+120°C	±0.25°C

Common characteristics

Temperature drift @ 20°C 0.003%/°C

4 wires Pt100 and 2 wires Pt1000 Probes

Model	Type	Application field	Accuracy
TP47.100.0 1/3 DIN Thin Film	4 wires Pt100	-50...+250°C	1/3 DIN
TP47.1000.0 1/3 DIN Thin Film	2 wires Pt1000	-50...+250°C	1/3 DIN
TP87.100.0 1/3 DIN Thin Film	4 wires Pt100	-50...+200°C	1/3 DIN
TP87.1000.0 1/3 DIN Thin Film	2 wires Pt1000	-50...+200°C	1/3 DIN

Common features

Temperature drift @20°C

Pt100 0.003%/°C

Pt1000 0.005%/°C

ORDERING CODES

HD2305.0: The kit is supplied with: instrument HD2305.0, 3 1.5V alkaline batteries, operating manual, case.

Electrodes, calibration solutions and temperature probes have to ordered separately.

HD22.2: Laboratory electrode holder composed of base plate with built-in magnetic stirrer, shaft and replaceable electrode holder. Suitable diameter 12mm. Powered by power supplier SWD10 (optional).

HD22.3: Laboratory electrode holder composed of base plate. Flexible arm for free positioning. Suitable for electrodes with diameter 12mm.

pH Electrodes

KP 20: Gel pH filled combined electrode for general use, with S7 screw connector, EPOXY body.

KP 30: Gel pH combined electrode for general use, 1m cable with BNC, EPOXY body.

KP 50: Gel pH combined electrode, porous Teflon ring junction, suitable for emulsions, demineralised water, waste water with S7 screw connector, glass body.

KP 61: 3 diaphragm liquid filled pH combined electrode for wine, milk, cream, etc., S7 screw connector, liquid reference filling, glass body.

KP 62: 1 diaphragm gel pH combined electrode for general use, pure water, varnishes, gel filled, S7 screw connector, glass body.

KP 63: liquid filled pH combined electrode for general use, varnishes, 1m cable with BNC, glass body.

KP 64: Liquid filled pH combined electrode, Teflon ring diaphragm, for wine, varnishes, emulsions, S7 screw connector, glass body.

KP 70: Pointed gel combined pH microelectrode diam. 6 x L=70 mm., with S7 screw connector, EPOXY body, glass tip, open junction for meat and cheese.

KP 80: Pointed gel pH combined electrode, with S7 screw connector, glass body, for cream, milk, viscous material, open junction.

KP100: Flat membrane gel combined pH electrode with S7 screw connector, glass body, for skin, leather, paper.

Characteristics and dimensions of the probes at page WA-76

CP: 1.5m extension cable with BNC/S7 connector for electrode without cable, thread S7.

CP5: 5m extension cable with BNC/S7 connector for electrode without cable, thread S7.

CP10: 10m extension cable with BNC/S7 connector for electrode without cable, thread S7.

CP15: 15m extension cable with BNC/S7 connector for electrode without cable, thread S7.

CE: S7 screw connector for pH electrode.

BNC: female BNC for extension cable

ORP Electrodes

KP 90: REDOX PLATINUM liquid filled electrode with S7 screw connector, glass body.

KP 91: Gel REDOX PLATINUM electrode, 1m cable with BNC, EPOXY body for general purpose light duty.

Characteristics and dimensions of the probes at page WA-76

pH Buffer solutions

HD8642: Buffer solution 4.01pH - 200cc.

HD8672: Buffer solution 6.86pH - 200cc.

HD8692: Buffer solution 9.18pH - 200cc.

Redox Buffer solutions

HDR220: Redox buffer solution 220mV 0.5 l.

HDR468: Redox buffer solution 468mV 0.5 l.

Electrolyte solutions

KCL3M Ready to use solution for electrode refilling – 100 cc



HD22.3



HD8642

HD8672

HD8692

Cleaning and maintenance

HD62PT: Diaphragm cleaning (tiourea in HCl) - 500ml.

HD62PP: Protein cleaning (pepsin in HCl) - 500ml.

HD62RF: Regeneration (fluorhydric acid) - 100ml.

HD62SC: Solution for electrode preservation - 200ml.

Temperature probes complete with SICRAM module

TP472I: Wire wound Pt100 sensor, immersion probe. Stem Ø 3 mm, length 300 mm. Cable length 2 m.

TP472I.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 230 mm. Cable length 2 m.

TP473PI: Wire wound Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP473PO: Thin film Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP474C.I: Wire wound Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP474C.0: Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP475A.0: Thin film Pt100 sensor, air probe. Stem Ø 4mm, length 230mm. Cable length 2 m.

TP472I.5: Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 500 mm. Cable length 2 m.

TP472I.10: Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 1000mm. Cable length 2 m.

TP49A.0: Thin film Pt100 sensor, immersion probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AC.0: Thin film Pt100 sensor, contact probe. Stem Ø 4 mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AP.0: Thin film Pt100 sensor, penetration probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP875.I: Wire wound Pt100 sensor, 150mm diameter globe-thermometer equipped with handle and SICRAM module. Cable length 2 m.

TP876.I: Wire wound Pt100 sensor, 50mm diameter globe-thermometer equipped with handle and SICRAM module. Cable length 2 m.

TP87.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 70 mm. Cable length 2 m.

TP878.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 2 m.

TP878.1.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 5 m.

TP879.0: Thin film Pt100 sensor, penetration probe for compost. Stem Ø 8 mm, length 1000 mm. Cable length 2 m.

Temperature probes without SICRAM module

TP47.100.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP47.1000.0: Thin film Pt1000 sensor, immersion probe. Probe's stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP47: Connector for Pt100 4-wire and Pt1000 2-wire probes without SICRAM module.

TP87.100.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 70mm. 4-wires connection cable with connector, length 1 m.

TP87.1000.0: Thin film Pt1000 sensor, immersion probe. Stem Ø 3mm, length 70mm. 2-wires connection cable with connector, length 1 m.





HD 2156.1, HD 2156.2 pH METER - CONDUCTIVITY METER - THERMOMETER

The **HD2156.1** and **HD2156.2** are portable instruments with a large LCD display. They measure pH, mV, redox potential (ORP), conductivity, liquid resistivity, total dissolved solids (TDS) and salinity using combined 4-ring and 2-ring conductivity/temperature probes. Temperature only is measured by Pt100 or Pt1000 immersion, penetration, contact or air probes.

The pH electrode calibration, as well as manual, can be carried out on one, two or three points and the calibration sequence can be chosen from a list of 13 buffers.

The calibration of the conductivity probe can be performed automatically in one or more of the 147 $\mu\text{S}/\text{cm}$, 1413 $\mu\text{S}/\text{cm}$, 12880 $\mu\text{S}/\text{cm}$ or 111800 $\mu\text{S}/\text{cm}$ solutions.

The HD2156.2 instrument is a **datalogger**. It stores up to 20,000 sets of three measurements composed of pH or mV, conductivity or resistivity or TDS or salinity and temperature: these data can be transferred to a PC from the instrument connected via the RS232C or USB 2.0 serial ports. The storing interval, printing, and baud rate can be configured by the menu.

Both models are fitted with an RS232C serial port and can transfer the acquired measurements to a PC or to a portable printer in real time.

The **Max**, **Min** and **Avg** function calculates the maximum, minimum or average values. Other functions include: the Auto-HOLD function and the automatic turning off which can also be excluded.

The instruments have IP66 protection degree.

INSTRUMENT TECHNICAL CHARACTERISTICS

Measured quantities: pH, mV, χ , Ω , TDS, NaCl, $^{\circ}\text{C}$, $^{\circ}\text{F}$

Instrument

Dimensions (Length x Width x Height)	185x90x40mm
Weight	470g (complete with batteries)
Materials	ABS, rubber
Display	2x4½ digits plus symbols Visible area: 52x42mm

Operating conditions

Working temperature	-5...50 $^{\circ}\text{C}$
Storage temperature	-25...65 $^{\circ}\text{C}$
Working relative humidity	0...90%RH without condensation
Protection degree	IP66

Power

Batteries	4 1.5V type AA batteries
Autonomy	200 hours with 1800mAh alkaline batteries
Power absorbed with instrument off	20 μA
Mains-supply unit	Output mains adapter 12Vdc/1000mA

Security of memorized data

Unlimited, independent of battery charge conditions

Time

Date and time	Real time schedule
Accuracy	1min/month max error

Measured values storage - model **HD2156.2**

Type	2000 pages containing 10 samples each
Quantity	20,000 sets of three measurements composed of pH or mV, χ , Ω or TDS or salinity and temperature.
Storage interval	1s, 5s, 10s, 15s, 30s, 1min, 2min, 5min, 10min, 15min, 20min, 30min and 1h.

Serial interface RS232C

Type	RS232C electrically isolated
Baud rate	Can be set from 1200 to 38400 bauds
Data bit	8
Parity	None
Stop bit	1
Flow Control	Xon/Xoff
Serial cable length	Max 15m
Print interval	Immediate or selectable between: 1s, 5s, 10s, 15s, 30s, 1min, 2min, 5min, 10min, 15min, 20min, 30min and 1h.

USB interface - model **HD2156.2**

Type	1.1 - 2.0 electrically isolated
------	---------------------------------

Connections

pH/mV input	Female BNC connector
Conductivity and Temperature input	8-pole male DIN45326 connector
Serial RS232C interface	8-pole MiniDin connector
USB interface	MiniUSB B-type connector
Mains adapter	2-pole connector (positive at centre)

Measurement of pH by Instrument

Measurement range	-2.000...+19.999pH
Resolution	0.01 or 0.001pH selectable from menu
Accuracy	$\pm 0.001\text{pH} \pm 1$ digit
Input impedance	$> 10^{12}\Omega$
Calibration error @25 $^{\circ}\text{C}$	Offset >20mV Slope < 50mV/pH or Slope > 63mV/pH Sensitivity < 85% or Sensitivity > 106.5%
Temperature compensation automatic/manual	-50...+150 $^{\circ}\text{C}$

Measurement of mV by Instrument

Measurement range	-1999.9...+1999.9mV
Resolution	0.1mV
Accuracy	$\pm 0.1\text{mV} \pm 1$ digit
Drift after 1 year	0.5mV/year

Measurement of conductivity			Resolution			Measurement of total dissolved solids (with coefficient χ /TDS=0.5)			Resolution		
Measuring range Kcell=0.01	0.000...1.999 μ S/cm	0.001 μ S/cm	Measuring range Kcell=0.01	0.000...1.999mg/l	0.005mg/l	Measuring range Kcell=0.01	0.000...1.999mg/l	0.005mg/l	Measuring range Kcell=0.01	0.000...1.999mg/l	0.005mg/l
Measuring range Kcell=0.1	0.00...19.99 μ S/cm	0.01 μ S/cm	Measuring range Kcell=0.1	0.00...19.99mg/l	0.05mg/l	Measuring range Kcell=0.1	0.00...19.99mg/l	0.05mg/l	Measuring range Kcell=0.1	0.00...19.99mg/l	0.05mg/l
Measuring range Kcell=1	0.0...199.9 μ S/cm	0.1 μ S/cm	Measuring range Kcell=1	0.0...199.9mg/l	0.5mg/l	Measuring range Kcell=1	0.0...199.9mg/l	0.5mg/l	Measuring range Kcell=1	0.0...199.9mg/l	0.5mg/l
	200...1999 μ S/cm	1 μ S/cm		200...1999mg/l	1mg/l		200...1999mg/l	1mg/l		200...1999mg/l	1mg/l
	2.00...19.99mS/cm	0.01mS/cm		2.00...19.99g/l	0.01g/l		2.00...19.99g/l	0.01g/l		2.00...19.99g/l	0.01g/l
	20.0...199.9mS/cm	0.1mS/cm		20.0...99.9g/l	0.1g/l		20.0...99.9g/l	0.1g/l		20.0...99.9g/l	0.1g/l
Measuring range Kcell=10	200...1999mS/cm	1mS/cm	Measuring range Kcell=10	100...999g/l	1g/l	Measuring range Kcell=10	100...999g/l	1g/l	Measuring range Kcell=10	100...999g/l	1g/l
Accuracy (conductivity)	$\pm 0.5\% \pm 1$ digit		Accuracy (total dissolved solids)	$\pm 0.5\% \pm 1$ digit		Accuracy (total dissolved solids)	$\pm 0.5\% \pm 1$ digit		Accuracy (total dissolved solids)	$\pm 0.5\% \pm 1$ digit	
Measurement of instrument's resistivity			Measurement of salinity			Measurement of salinity			Resolution		
Measuring range Kcell=0.01	up to 1G Ω ·cm (*)		Measurement range	0.000...1.999g/l	1mg/l	Measurement range	0.000...1.999g/l	1mg/l	Measurement range	0.000...1.999g/l	1mg/l
Measuring range Kcell=0.1	up to 100M Ω ·cm (*)			2.00...19.99g/l	10mg/l		2.00...19.99g/l	10mg/l		2.00...19.99g/l	10mg/l
Measuring range Kcell=1	5.0...199.9 Ω ·cm	0.1 Ω ·cm		20.0...199.9g/l	0.1g/l		20.0...199.9g/l	0.1g/l		20.0...199.9g/l	0.1g/l
	200...999 Ω ·cm	1 Ω ·cm	Accuracy (salinity)	$\pm 0.5\% \pm 1$ digit		Accuracy (salinity)	$\pm 0.5\% \pm 1$ digit		Accuracy (salinity)	$\pm 0.5\% \pm 1$ digit	
	1.00k...19.99k Ω ·cm	0.01k Ω ·cm	Temperature compensation	0...100°C with α_T selectable from 0.00 to 4.00%/°C		Temperature compensation	0...100°C with α_T selectable from 0.00 to 4.00%/°C		Temperature compensation	0...100°C with α_T selectable from 0.00 to 4.00%/°C	
	20.0k...99.9k Ω ·cm	0.1k Ω ·cm	automatic/manual			automatic/manual			automatic/manual		
	100k...999k Ω ·cm	1k Ω ·cm	Reference temperature	20°C or 25°C		Reference temperature	20°C or 25°C		Reference temperature	20°C or 25°C	
Measuring range Kcell=10	1...10M Ω ·cm	1M Ω ·cm	χ / TDS Conversion factor	0.4...0.8		χ / TDS Conversion factor	0.4...0.8		χ / TDS Conversion factor	0.4...0.8	
Accuracy (resistivity)	$\pm 0.5\% \pm 1$ digit		Cell constant K (cm ⁻¹)	0.01, 0.7, 1.0 and 10.0		Cell constant K (cm ⁻¹)	0.01, 0.7, 1.0 and 10.0		Cell constant K (cm ⁻¹)	0.01, 0.7, 1.0 and 10.0	

(*) The resistivity measurement is obtained from the reciprocal of the conductivity measurement: the indication of the resistivity, in the vicinity of the full scale, appears as in the following table

K cell = 0.01 cm ⁻¹		K cell = 0.1 cm ⁻¹	
Conductivity (μ S/cm)	Resistivity (M Ω ·cm)	Conductivity (μ S/cm)	Resistivity (M Ω ·cm)
0.001 μ S/cm	1000 M Ω ·cm	0.01 μ S/cm	100 M Ω ·cm
0.002 μ S/cm	500 M Ω ·cm	0.02 μ S/cm	50 M Ω ·cm
0.003 μ S/cm	333 M Ω ·cm	0.03 μ S/cm	33 M Ω ·cm
0.004 μ S/cm	250 M Ω ·cm	0.04 μ S/cm	25 M Ω ·cm

Standard solutions automatically detected @25°C

147 μ S/cm
1413 μ S/cm
12880 μ S/cm
111800 μ S/cm

TECHNICAL DATA OF PROBES EQUIPPED WITH INSTRUMENT

2 and 4 electrode conductivity probes

ORDER CODE	MEASUREMENT RANGE	DIMENSIONS
SP06T	K=0.7 5 μ S...200mS/cm 0...90°C 4-electrode cell in Pocaan/Platinum Max pressure 5bar	
SPT 401.001 not suitable for HD 2306.0	K=0.01 0.04...20 μ S/cm 0...120°C 2-electrode cell AISI 316 - Teflon Max pressure 5bar	
SPT01G	K=0.1 0.1 μ S...500 μ S/cm 0...80°C 2-electrode cell in Glass/Platinum Max pressure 5bar	
SPT1G	K=1 10 μ S...10mS/cm 0...80°C 2-electrode cell in Glass/Platinum Max pressure 5bar	
SPT10G	K=10 500 μ S...200mS/cm 0...80°C 2-electrode cell in Glass/Platinum Max pressure 5bar	

Temperature measurement of the instrument

Pt100 measuring range	-50...+200°C
Pt1000 measuring range	-50...+200°C
Resolution	0.1°C
Accuracy	±0.1°C ±1 digit
Drift after 1 year	0.1°C/year

TECHNICAL DATA OF PROBES AND MODULES EQUIPPED WITH INSTRUMENT Temperature probes Pt100 sensor with SICRAM module

Model	Type	Application field	Accuracy
TP472I	Immersion	-196°C...+500°C	±0.25°C (-196°C...+300°C) ±0.5°C (+300°C...+500°C)
TP472I.0 1/3 DIN Thin Film	Immersion	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP473PI	Penetration	-50°C...+400°C	±0.25°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP473P.0 1/3 DIN Thin Film	Penetration	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP474C.I	Contact	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP474C.0 1/3 DIN Thin Film	Contact	-50°C...+300°C	±0.3°C (-50°C...+300°C)
TP475A.0 1/3 DIN Thin Film	Air	-50°C...+250°C	±0.3°C (-50°C...+250°C)
TP472I.5	Penetration	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP472I.10	Penetration	-50°C...+400°C	±0.30°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP49A.0 Class A Thin Film	Immersion	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AC.0 Class A Thin Film	Contact	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AP.0 Class A Thin Film	Penetration	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP875.I	Globe-thermometer Ø150mm	-30°C...+120°C	±0.25°C
TP876.I	Globe-thermometer Ø50mm	-30°C...+120°C	±0.25°C
TP87.0 1/3 DIN Thin Film	Immersion	-50°C...+200°C	±0.25°C
TP878.0 1/3 DIN Thin Film	Photovoltaic	+4°C...+85°C	±0.25°C
TP878.1.0 1/3 DIN Thin Film			
TP879.0 1/3 DIN Thin Film	Compost	-20°C...+120°C	±0.25°C

Common characteristics

Temperature drift @ 20°C 0.003%/°C

4 wires Pt100 and 2 wires Pt1000 Probes

Model	Type	Application field	Accuracy
TP47.100.0 1/3 DIN Thin Film	4 wires Pt100	-50...+250°C	1/3 DIN
TP47.1000.0 1/3 DIN Thin Film	2 wires Pt1000	-50...+250°C	1/3 DIN
TP87.100.0 1/3 DIN Thin Film	4 wires Pt100	-50...+200°C	1/3 DIN
TP87.1000.0 1/3 DIN Thin Film	2 wires Pt1000	-50...+200°C	1/3 DIN

Common features

Temperature drift @20°C

Pt100 0.003%/°C

Pt1000 0.005%/°C

ORDER CODES

HD2156.1: The kit is supplied with: instrument HD2156.1, 4 1.5V alkaline batteries, operating manual, case and DeltaLog9 software.

HD2156.2: The kit is supplied with: instrument HD2156.2 **datalogger**, 4 1.5V alkaline batteries, operating manual, case and DeltaLog9 software.

pH/mV probes, conductivity probes, temperature probes, standard calibration solutions for various types of measures, connection cables for pH electrodes with S7 connector, cables for data transfer to PC or printer have to be ordered separately.

HD2110CSNM: 8-pole connection cable MiniDin - Sub D 9-pole female for RS232C.

C.206: Serial connection cable with USB connector for PC and 8-pole MiniDin male connector for the instrument HD2156.1.

CP23: USB 2.0 connection cable type A - MiniUSB type B (not suitable for HD2156.1).

DeltaLog9: Software for download and management of the data on PC using Windows operating systems.

SWD10: Stabilized power supply 100-240 Vac/12Vdc-1A mains voltage


HD40.1: 24-column portable thermal printer, serial interface, 57mm paper width, four NiMH 1.2V rechargeable batteries, SWD10 power supply, instruction manual, 5 thermal paper rolls. It uses the cable HD2110CSNM (optional).

RCT: The kit includes 4 thermal paper rolls 57mm wide and 32mm in diameter.

BAT-40: Spare battery pack for HD40.1 printer with built-in temperature sensor.

HD22.2: Laboratory electrode holder composed of base plate with built-in magnetic stirrer, shaft and replaceable electrode holder. Suitable diameter 12mm. Powered by bench-top meters of the series HD22...with cable HD22.2.1 (optional) or power supplier SWD10 (optional).

HD22.3: Laboratory electrode holder composed of base plate. Flexible arm for free positioning. Suitable for electrodes with diameter 12mm.



A For the models of portable data logger series **HD21XX.2** has been implemented with a new serial port miniUSB type HID (Human Interface Device).

When making the connection to the PC by the USB cable Type A - Mini USB B-type coded CP23, **no USB driver installation is requested.**

B For the connection of the models **HD21XX.1** to the RS232 port of your PC, the USB/serial converter is available (**code C.206**). The converter is equipped with its own drivers that have to be installed before connecting the converter to the PC (please see the details in the CDRom supplied with the converter).

C The port with the MiniDIN connector which is present on every model is an RS232C type. By means of the cable coded HD2110CSNM, an RS232 port of a PC or the HD40.1. printer can be connected.

pH Electrodes

- KP 20:** Gel pH filled combined electrode for general use, with S7 screw connector, EPOXY body.
- KP 30:** Gel pH combined electrode for general use, 1m cable with BNC, EPOXY body.
- KP 50:** Gel pH combined electrode, porous Teflon ring junction, suitable for emulsions, demineralised water, waste water with S7 screw connector, glass body.
- KP 61:** 3 diaphragm liquid filled pH combined electrode for wine, milk, cream, etc., S7 screw connector, liquid reference filling, glass body.
- KP 62:** 1 diaphragm gel pH combined electrode for general use, pure water, varnishes, gel filled, S7 screw connector, glass body.
- KP 63:** liquid filled pH combined electrode for general use, varnishes, 1m cable with BNC, glass body.
- KP 64:** Liquid filled pH combined electrode, Teflon ring diaphragm, for wine, varnishes, emulsions, S7 screw connector, glass body.
- KP 70:** Pointed gel combined pH microelectrode diam. 6 x L=70 mm., with S7 screw connector, EPOXY body, glass tip, open junction for meat and cheese.
- KP 80:** Pointed gel pH combined electrode, with S7 screw connector, glass body, for cream, milk, viscous material, open junction.
- KP100:** Flat membrane gel combined pH electrode with S7 screw connector, glass body, for skin, leather, paper.

Characteristics and dimensions of the probes at page WA-76

- CP:** 1.5m extension cable with BNC/S7 connector for electrode without cable, thread S7.
- CP 5:** 5m extension cable with BNC/S7 connector for electrode without cable, thread S7.
- CP 10:** 10m extension cable with BNC/S7 connector for electrode without cable, thread S7.
- CP 15:** 15m extension cable with BNC/S7 connector for electrode without cable, thread S7.
- CE :** S7 screw connector for pH electrode.
- BNC:** female BNC for extension cable

ORP Electrodes

- KP 90:** REDOX PLATINUM liquid filled electrode with S7 screw connector, glass body.
- KP 91:** Gel REDOX PLATINUM electrode, 1m cable with BNC, EPOXY body for general purpose light duty.

Characteristics and dimensions of the probes at page WA-76

pH Buffer solutions

- HD8642:** Buffer solution 4.01pH - 200cc.
- HD8672:** Buffer solution 6.86pH - 200cc.
- HD8692:** Buffer solution 9.18pH - 200cc.

Redox Buffer solutions

- HDR220:** Redox buffer solution 220mV 0.5 l.
- HDR468:** Redox buffer solution 468mV 0.5 l.

Electrolyte solutions

- KCL3M** Ready to use solution for electrode refilling – 100 cc

Cleaning and maintenance

- HD62PT:** Diaphragm cleaning (tiourea in HCl) - 500ml.
- HD62PP:** Protein cleaning (pepsin in HCl) - 500ml.
- HD62RF:** Regeneration (fluorhydric acid) - 100ml.
- HD62SC:** Solution for electrode preservation - 200ml.

Conductivity probes

Please refer to the purchasing codes reported on the table of the probes on page WA-9.

Conductivity buffer solutions

- HD 8747:** Calibration solution 0.001 mol/l corresponding to 147 µS/cm at 25°C, 200cc.
- HD 8712:** Calibration solution 0.1 mol/l corresponding to 12,880 µS/cm at 25°C, 200cc.
- HD 8714:** Calibration solution 0.01 mol/l corresponding to 1413 µS/cm at 25°C, 200cc.
- HD 87111:** Calibration solution 1 mol/l corresponding to 111800 µS/cm at 25°C, 200cc.

Temperature probes equipped with SICRAM module

- TP472I:** Wire wound Pt100 sensor, immersion probe. Stem Ø 3 mm, length 300 mm. Cable length 2 m.
- TP472I.0:** Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 230 mm. Cable length 2 m.
- TP473PI:** Wire wound Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

- TP473PO:** Thin film Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.
- TP474C.I:** Wire wound Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.
- TP474C.O:** Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.
- TP475A.O:** Thin film Pt100 sensor, air probe. Stem Ø 4mm, length 230mm. Cable length 2 m.
- TP472I.5:** Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 500 mm. Cable length 2 m.
- TP472I.10:** Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 1000mm. Cable length 2 m.
- TP49A.O:** Thin film Pt100 sensor, immersion probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle
- TP49AC.O:** Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 150mm. Cable length 2 m. Aluminium handle
- TP49AP.O:** Thin film Pt100 sensor, penetration probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle
- TP875.I:** Wire wound Pt100 sensor, 150mm diameter globe-thermometer equipped with handle. Cable length 2 m.
- TP876.I:** Wire wound Pt100 sensor, 50mm diameter globe-thermometer equipped with handle. Cable length 2 m.
- TP87.O:** Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 70 mm. Cable length 2 m.
- TP878.O:** Thin film Pt100 sensor, contact probe for solar panels. Cable length 2 m.
- TP878.1.O:** Thin film Pt100 sensor, contact probe for solar panels. Cable length 5 m.
- TP879.O:** Thin film Pt100 sensor, penetration probe for compost. Stem Ø 8 mm, length 1000 mm. Cable length 2 m.

Temperature probes without SICRAM module

- TP47.100.O:** Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.
- TP47.1000.O:** Thin film Pt1000 sensor, immersion probe. Probe's Stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.
- TP47:** Connector for Pt100 4-wire and Pt1000 2-wire probes without SICRAM module.
- TP87.100.O:** Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 70mm. 4-wires connection cable with connector, length 1 m.
- TP87.1000.O:** Thin film Pt1000 sensor, immersion probe. Stem Ø 3mm, length 70mm. 2-wires connection cable with connector, length 1 m.





TECHNICAL CHARACTERISTICS

pH simulation:	0 ÷ 14 pH
pH resolution:	0.1 pH
pH accuracy 20 ÷ 25°C:	0.002 pH
Thermal drift:	±0.0005 pH/°C from -5°C to 20°C and from 25°C to 50°
mV simulation:	±1999 mV
mV resolution:	1 mV
mV accuracy:	±100 µV
Thermal drift mV scale:	-199.9 ... +199.9: ±0.01 mV/°C from -5 to 20°C and from 25 to 50°C
Thermal drift mV scale:	-1999 ... +1999: ±0.05 mV/°C from -5 to 20°C and from 25 to 50°C
Noise 0 ÷ 10 Hz:	1 µV peak/peak
Simulation of temperature compensation:	-20 to 150°C (-4 to 302°F)
Output impedance:	100 kΩ 1%, 1 GΩ 5% (no leading load capacity)
Display:	LCD 2 lines, 3 ½ digits. Figure height approx. 12.5 mm.
Symbols:	pH, mV, °C, °F, HI imp., LO imp., 0.1 pH, 1 pH, 1 mV, 10 mV
Signals:	LOU, ER1, CAL
Working temperature:	-5 to 50°C (23 to 122°F)
Power supply:	9 Vdc alkaline battery. Low battery indication.
Consumption (instrument only):	5 mA lit up, 20 µA turned off
Autonomy:	about 200 hours
Dimensions:	187 x 72 x 38 mm.
Weight:	300 gr

ORDER CODES

HD 9609: The kit is made of: instrument HD 9609, adapter cables CP 9509BNC, CP 9509 T, carrying case

CP 9509BNC: Adapter cable L = 1 mt, male BNC connector on both ends

CP 9509 T: Adapter cable L = 1 mt, BNC connector on only one end

CP 9509 S7: Adapter cable L = 1 mt, BNC wall connector one end, S7 male connector on the other end.

HD 9609

pH AND mV SIMULATOR

GENERAL CHARACTERISTICS

The simulator **HD 9609** is a portable instrument for checking and calibrating pH and mV measuring instruments. The characteristics of this instrument satisfy any checking and calibrating requirements for both portable and panel-mounted instruments; it can be used in laboratories, in industry or for check out on field.

Despite its many functions, the instrument is simple to be used: a large display, with dual indication, and a series of symbols allow it to be used even by unskilled personnel.

The HD9609 outputs on channel A the simulation of signals of an electrode for measuring pH, ORP, ISFET, in the range:

- 0 to 14 pH, with resolution 0.10 pH;
- ±1999 mV, with resolution 1 mV.

It is possible to choose between two output impedance values:

- 100 kΩ, low impedance;
- 1 GΩ, high impedance.

The simulation of the temperature compensation electrode is programmed manually in the range from -20°C to +150°C, while the temperature is measured in Celsius or Fahrenheit degrees.

The pH simulation values can be manually set as desired, in steps of 0.1 or 1 pH. The mV simulation values can be manually set as desired, in steps of 1 or 10 mV.

The HD9609 is powered by an ordinary 9Vdc alkaline battery.

The electronics are housed in a sturdy ABS case with ergonomic lines.

In designing and making the instrument, each detail has been assessed and selected in order to provide an instrument with high performance and excellent long-term measurement stability.





DO 9403T-R1 pH/mV TRANSMITTER

The **DO 9403T-R1** pH transmitter converts the output of a pH electrode, with temperature compensation, into a 4÷20 mA signal. The pH or Redox electrode input circuit is galvanically isolated from the 4÷20 mA output signal.

An LCD indicator allows to display the value of the process signal and the various parameters. The accurate design and choice of the components make the instrument precise and reliable over time

The instrument works with a pH or Redox electrode and a temperature probe (Pt100 sensor, 100 Ω at 0°C).

Key functions

PRG Programming of the parameters is activated by pressing the PRG key. The Δ symbol lights up on the display and the message P1 appears to indicate that the parameter P1 is being programmed. Continuing pressing the PRG key, the messages P2, P3, P4, P5, P6, P7, P8, P9, P10 and the corresponding parameters are displayed in sequence. After P10 the instrument returns to normal function.

After displaying the parameter of interest, it is possible to view the value by pressing the OK button. In order to change the parameter, use the ▲ and ▼ buttons. Press the OK button again to confirm the value of the parameter.

SET Key for setting the relay intervention threshold. The Δ symbol and the REL symbol appear on the display, fixed or flashing, indicating the switching on or off threshold of relay A or of relay B.

°C/°F

- If this key is pressed it changes the temperature measuring unit to degrees Celsius or degrees Fahrenheit.
- When pressed together with the CAL key it activates the manual temperature setting function.
- If pressed during the pH calibration function it quits the calibration function without storing the calibration.




Fig.1 Active transmitter.




Fig.2 Active indicator.




Fig.3 Passive transmitter.




Fig.4

pH/mV - If this key is pressed it changes the measuring unit to mV or pH.

- When pressed together with the CAL key it activates the pH calibration function.

OK Confirms the programming parameters, or the relay SET values, and stores them.

CAL - When pressed together with the °C/°F key it activates the manual temperature setting function.

- When pressed together with the pH/mV key it activates the pH calibration function.

- Key used to confirm pH calibration and manual temperature setting.

▲ - Key for increasing the value displayed in the parameter programming phase.

- During the relay SET point programming phase.

- During the calibration phase.

▼ - Key for decreasing the value displayed in the parameter programming phase.

- During the relay SET point programming phase.

- During the calibration phase.

Setting the relay SET point

- Press the SET button; the Δ symbol appears on the display.
- The REL symbol and the letter A also light up on the display to indicate that the value shown corresponds to the switching on threshold of relay A.
- To change this value press the ▲ and ▼ buttons.
- Press SET; the REL symbol flashes and the letter A remains lit to indicate that the value shown corresponds to the switching off threshold of relay A.
- To change this value press the ▲ and ▼ buttons.
- Press the SET button; the REL symbol and the letter B light up to indicate that the value shown corresponds to the switching off threshold of relay B.
- To change this value press the ▲ and ▼ buttons.
- Press SET; the REL symbol flashes and the letter B remains lit to indicate that the value shown corresponds to the switching off threshold of relay B.
- To change this value press the ▲ and ▼ buttons.
- Press SET, the instrument stores the set parameters and returns to normal function. The REL and Δ symbols disappear.

NOTE: During the SET point setting phase (REL symbol lit or flashing) the instrument returns to normal function if no key is pressed for 2 minutes.

Manual temperature setting


If the temperature probe is not connected or if the probe is broken the measuring unit °C or °F flashes. In this case it is possible to set the temperature compensation value manually.

- Press the CAL key and the °C/°F key together; the Δ symbol appears and the manual temperature is shown with the measuring unit flashing.
- By using the ▲ and ▼ buttons, set the temperature value corresponding to the temperature of the liquid in which you wish to measure the pH.
- Press CAL to confirm this value. The Δ symbol switches off and the instrument returns to the previous display.

Calibration of the pH electrode

Calibration of the offset of the pH electrode:

- Immerse the electrode in the buffer solution used for calibrating the offset (6.86 pH).
- Press the CAL key and the pH/mV key together; the Δ symbol lights up on the display.
- By using the ▲ and ▼ buttons, adjust the pH value measured as a function of the liquid temperature.
- Press CAL to confirm this value. The Δ symbol switches off.



Dimensions

Calibration of the slope of the pH electrode:

- Immerse the electrode in the buffer solution used for calibrating the slope (4.01 or 9.18 pH).
- Press the CAL key and the pH/mV key together; the Δ symbol lights up on the display.
- By using the ▲ and ▼ buttons, adjust the pH value measured as a function of the liquid temperature.
- Press CAL to confirm this value. The Δ symbol switches off.

NOTE: If you want to quit without storing the new calibration, press the °C/°F key.

N.B.: The instrument can automatically recognize three standard calibration solutions: 4.01 pH, 6.86 pH and 9.18 pH.

Programming the parameters

- P1** Control of relay unit and analog output, pH or mV.
- P2** pH/mV value corresponding to 4 mA at output. May be set between -1.00 pH and 15.00 pH or between -1999 mV and +1999 mV.
- P3** pH/mV value corresponding to 20 mA at output. May be set between -1.00 pH and 15.00 pH or between -1999 mV and +1999 mV.
- P4** Delay time in the intervention of relay A. May be set between 0 and 255 seconds.
- P5** Delay time in the intervention of relay B. May be set between 0 and 255 seconds.
- P6** Calibration of Pt100 probe.
- P7** 4 mA output current calibration.
- P8** 20 mA output current calibration.
- P9** Input voltage calibration.
- P10** Display of the offset voltage value and of the slope value of the electrode.

To change one of these parameters press key PRG until the message corresponding to the parameter to be changed appears on the screen.

Press OK to show the value of the parameter.

By using the ▲ and ▼ buttons, bring the parameter displayed to the desired value.

Press OK again to confirm.

Parameter P10 cannot be altered, it can only be displayed.

NOTE **P6-P7-P8-P9: calibration procedure to be carried out at a laboratory by skilled workers.**

Combined electrode input	pH	-1.00 pH...15.00 pH (-500...+500 mV)
	ORP	-1999...+1999 mV
	Input impedance	>10 Tohm
	Cable length	<50 metres screened (about 5 nF)
	Accuracy	0.1% of reading ±1 digit ±0.01% of pH per °C of temperature drift
Temperature input	Pt100 2/4 wires	-50...199.9°C
	Transducer energizing	0.5 mA DC
	Cable length	<10 metres unscreened <20 metres screened (about 2 nF)
	Accuracy	0.2°C ±0.1% of reading ±2 digits ±0.01°C/°C
pH electrode temp. compensation	Automatic	According to Nernst
Current output	4.00...20.00 mA	Programmable and proportional to the pH or mV value
	Accuracy	0.5% of reading ±0.02 mA
	Insulation	2500 Vac 1 minute
R Load	Load resistance	$R_{Lmax} = \frac{V_{dc}-10}{0,022}$ $R_{Lmax} = 636 \Omega @ V_{dc} = 24 V_{dc}$
Relay output	A and B	Bistable, contact 3 A/230 Vac free potential
Power supply	Active	24 or 230 Vac -15/+10% 1 VA, 48...62 Hz, see fig. 1
	Passive	4÷20 mA, 2 wire configuration, 10÷35 V, see fig. 2
Temp.	Operation	0...50°C
	Storage	-20...70°C, no condensation
Case	External dimensions	120x80x56 mm
	Protection class	IP64

Pt100 probe calibration (100 Ω at 0°C)

(Calibration procedure to be carried out at a laboratory by skilled workers)

- Connect the Pt100 probe to the instrument. Press the PRG key until the message P6 appears on the display.
- Press the OK key; the currently measured temperature appears on the display.
- Immerse the Pt100 probe and a precision thermometer for reference in the zero calibration bath. Wait long enough for the reading to become stable.
- By using the ▲ and ▼ buttons, adjust the value of the temperature measured by the Pt100 probe so that it corresponds with the value on the reference thermometer.
- Immerse the Pt100 probe and a precision thermometer in the full scale calibration bath. Wait long enough for the reading to become stable.
- By using the ▲ and ▼ buttons, adjust the value of the temperature measured by the Pt100 probe so that it corresponds with the value on the reference thermometer.
- Press OK to confirm. To quit programming, press PRG repeatedly.

N.B.: If the temperature shown by the instrument is between $\pm 12^{\circ}\text{C}$, the instrument calibrates the probe offset, otherwise it calibrates the gain.

Calibrating the analog output

(Calibration procedure to be carried out at a laboratory by skilled workers)

- Connect a precision milliamperometer to the analog output.
- Press the PRG key until the message P7 appears on the display.
- Press OK; the message 4.0 appears on the display, indicating calibration at 4 mA.
- By using the ▲ and ▼ buttons, adjust the value of the output current so as to have an indication of 4,00 mA on the precision milliammeter.
- Press the PRG key until the message P8 appears on the display.
- Press OK; the message 20.0 appears on the display, indicating calibration at 20 mA.
- Using the ▲ and ▼ keys, adjust the value of the output current so as to have an indication of 20,00 mA on the precision milliammeter.
- Press OK to confirm. To quit programming, press PRG repeatedly.

Calibrating the voltage input

(Calibration procedure to be carried out at a laboratory by skilled workers)

- Press the PRG key until the message P9 appears on the display.
- Press OK; the mV value of the input appears on the display.
- Simulate a voltage of 0 mV at the input (if the value is between ± 25 mV the zero is calibrated, otherwise the full scale value is calibrated).
- Using the ▲ and ▼ keys, adjust the voltage value so as to have the correct voltage value on the display.
- Press the SET key, the REL symbol lights up on the display indicating that the instrument is measuring the voltage present at the input using the second measurement scale.
- Using the ▲ and ▼ keys, adjust the voltage value so as to have the correct voltage value on the display.
- Press the SET key, the REL symbol on the display switches off.
- Simulate a voltage of 450 mV at the input, corresponding to the full scale value of the first scale.
- Using the ▲ and ▼ keys, adjust the voltage value so as to have the correct voltage value on the display.
- Simulate a voltage of 1800 mV at the input, corresponding to the full scale value of the second scale.
- Using the ▲ and ▼ keys, adjust the voltage value so as to have the correct voltage value on the display.
- Press OK to confirm. To quit programming, press PRG repeatedly.

Display

Symbol	description
°C	indicates that the value shown is in °C.
°F	indicates that the value shown is in °F.
pH	indicates that the unit of the value shown is pH.
mV	Indicates that the unit of the value shown is milli Volts.
A	indicates that the relay A is in closed status.
B	indicates that the relay B is in closed status.
REL	<ul style="list-style-type: none">- indicates that the value shown corresponds to the closing thresholds of the contacts of relay A or B;- indicates that the offset of the second voltage measurement scale is being calibrated.

REL flashing indicates that the value shown corresponds to the opening thresholds of the contacts of relay A or B.

- Δ
- indicates that the instrument is in the parameter setting phase;
 - indicates that the closing and opening thresholds of relays A and B are being changed;
 - indicates that the manual compensation temperature is being changed;
 - indicates that the pH electrode is being calibrated.

Error signals

- OFL** - Warning which appears during measurement when the value to be displayed is out of scale.
- E1** - Error warning which appears during pH calibration to indicate that the offset value of the electrode is too high in absolute value.
- E2** - Error warning which appears during pH calibration to indicate that the difference between the mV readings given by the two buffer solutions used for calibration is too great.
- E3** - Error warning which appears during pH calibration to indicate that the mV readings given by the two buffer solutions used for calibration are too close (about 50 mV at 25°C).
- E4** - Reading error on the EEPROM.
- E5** - Error warning indicating that the slope calculation gives a value 20% lower than the nominal value or gives a negative value.
- E6** - Error warning indicating that the slope calculation gives a value 150% lower than the nominal value.

Order code

DO 9403T-R1: pH transmitter 4÷20 mA passive or active, power supply 24 Vac, 120x80x56 mm for use in the field.

HD 882 M100/300: Temperature probe with Pt100 sensor, miniature head, shaft Ø6x300 mm.

HD 882 M100/600: Temperature probe with Pt100 sensor, DIN B head, shaft Ø6x600 mm.

HD 8642: Buffer solution 4.01 pH.

HD 8672: Buffer solution 6.86 pH.

HD 8692: Buffer solution 9.18 pH.

HD R220: Buffer solution redox 220 mV

HD R468: Buffer solution redox 468 mV

HD62PT: Diaphragm cleaning (tiourea in HCl) - 500ml.

HD62PP: Protein cleaning (pepsin in HCl) - 500ml.

HD62RF: Regeneration (fluorhydric acid) - 100ml.

HD62SC: Solution for electrode preservation - 200ml

CP5T: Extension cable. Connector S7/wire - TERMINAL BOARD.

CP5/10T: Extension cable L=10m.Connector S7/wire.

KPI 10: Combined industrial electrode, S7 PG13,5 connector, glass body, Ag/AgCl sat KCl, Ø 12x120, temperature 0÷130°C, porous Teflon fitting.

KPI 11: Combined industrial electrode, S7 PG13,5 connector, Rytron body, Ag/AgCl sat KCl, temperature 0÷100°C, porous Teflon fitting.

KPI 12: Platinum electrode for Redox measurement, S7 PG13,5 connector, Ag/AgCl sat KCl, temperature 0÷130°C, pressure 6 bar.

KPI 13: Platinum electrode for Redox measurement, Rytron body, S7 PG13,5 connector, temperature 0÷100°C Ag/AgCl sat KCl.

Electrode dimensions at page WA-18




DO 9785T - DO 9765T pH or mV TRANSMITTERS

DO 9785T/DO 9765T pH transmitters convert the output of a pH electrode, with temperature compensation, into a 4÷20 mA signal. The pH electrode input circuit is galvanically insulated against the 4÷20 mA output signal.

An LCD indicator allows viewing of the process signal value and of the various parameters. The accurate design and choice of components make the instrument precise and reliable for a long working life.

The instrument works with a pH electrode or Redox and a temperature probe (Pt 100 sensor, 100 Ω at 0°C).



Technical characteristics

Combined electrode input	pH	-1,00 pH...15,00 pH (-500...+500 mV)
	ORP	-1999...+1999 mV
	Input impedance	>10 ¹² Ω
	Cable length	<50m screened (about 5 nF)
	Accuracy	0.1% of reading ±1 digit ±0.01% pH per °C of drift in temperature
Temperature input	Pt100 2/4 wires	-50÷200°C
	Transducer energizing	0,5 mA DC
	Cable length	<10 metres unscreened <50 metres screened (about 2 nF)
	Accuracy	0.2°C ±0.1% of reading ±2 digits ±0.01°C/°C
pH electrode compensation temperature	Automatic	According to Nernst
	Manual	-50÷200°C
Current output	4.00...20.00 mA	Programmable and proportional to the pH or mV value
	Accuracy	0.5% of reading ±0.02 mA
	Insulation	2500 Vac 1 minute
R _L Load	Load resistance	$R_{Lmax} = \frac{V_{dc}-10}{0,022}$ $R_{Lmax} = 636 \Omega @ V_{dc} = 24 V_{dc}$
Relay output	A and B	Bistable, contact 3A/230 Vac free potential
Power supply	Passive	4÷20 mA, 2 wire configuration, 10÷35 V, see fig. 2
	Active	24 or 230 Vac - 15/+10% 1 VA, 48...62 Hz, see fig. 1
DO 9765T case	External dimensions	120x122x56 mm wall mounting
	Protection class	IP64
DO 9785T case	External dimensions	96x96x126 mm panel mounting
	Protection class	IP54

Key functions

PRG Programming of the parameters is activated by pressing the PRG key plus the ▲ and ▼ keys. The message P1 appears on the display, indicating that the parameter P1 is being programmed. When the PRG key is pressed continuously, the messages P2, P3, P4, P5, P6, P7 and the corresponding parameters are displayed in sequence. After P7 the instrument returns to normal function.

SET Key for setting the relay intervention threshold. The ON or OFF symbol appears on the display, indicating the switching on or off threshold of relay A or of relay B.

°C/°F - If this key is pressed it changes the temperature measuring unit to degrees Celsius or degrees Fahrenheit.

- When pressed together with the CAL key it activates the manual temperature setting function.
- If pressed during the conductivity calibration function it quits the calibration function without storing the calibration.

pH/mV - If this key is pressed it changes the measuring unit to mV or pH.
- When pressed together with the CAL key it activates the pH calibration function.

OK Confirms the programming parameters, or the relay SET values, and stores them.

CAL - When pressed together with the °C/°F key it activates the manual temperature setting function.
- When pressed together with the pH/mV key it activates the pH calibration function.

- Key used to confirm pH calibration and manual temperature calibration.
▲ - Key for increasing the value displayed in the parameter programming phase.
- During the relay SET point programming phase.
- During the calibration phase.

▼ - Key for decreasing the value displayed in the parameter programming phase.
- During the relay SET point programming phase.
- During the calibration phase.

Setting the relay SET point

- Press the SET button; the ON symbol appears on the display with the letter A to indicate that the value shown corresponds to the switching on threshold of relay A.
- To change this value press the ▲ and ▼ keys.

- Press SET; the OFF symbol appears with the letter A to indicate that the switching off threshold of relay A is being displayed.
- To change this value press the ▲ and ▼ keys.
- Press the SET button; the ON symbol appears on the display with the letter B to indicate that the value shown corresponds to the switching on threshold of relay B.
- To change this value press the ▲ and ▼ keys.
- Press SET; the OFF symbol appears with the letter B to indicate that the switching off threshold of relay B is being displayed.
- To change this value press the ▲ and ▼ keys.
- Press SET, the instrument stores the values and returns to normal function.

NOTE: During the SET point setting phase (symbols ON or OFF lit) the instrument returns to normal function if no key is pressed for 2 minutes.

Temperature setting for manual compensation

If the temperature probe is not connected or if the probe is broken the measuring unit °C or °F flashes. In this case it is possible to set the temperature compensation value manually.

- Press the CAL key and the °C/°F key together; the message CAL appears at the bottom of the display.
- Using the ▲ and ▼ keys, set the temperature value corresponding to the temperature of the liquid in which you wish to measure the pH value.
- Press CAL to confirm this value. The message CAL disappears.

Calibration of the DO 9785T/DO 9765T with pH electrode

Calibration of the offset of the pH electrode:

- Immerse the electrode in the buffer solution used for calibrating the offset (6.86 pH).
- Press the CAL key and the pH/mV key together; the message CAL appears at the top of the display.
- Using the ▲ and ▼ keys, adjust the pH value measured as a function of the liquid temperature.
- Press CAL to confirm this value. The message CAL disappears.

Calibration of the slope of the pH electrode:

- Immerse the electrode in the buffer solution used for calibrating the slope (4.01 or 9.18 pH).
- Press the CAL key and the pH/mV key together; the message CAL appears at the top of the display.
- Using the ▲ and ▼ keys, adjust the pH value measured as a function of the liquid temperature.
- Press CAL to confirm this value. The message CAL disappears.

NOTE: If you want to quit without storing the new calibration, press the °C/°F key.

N.B.: The instrument can automatically recognize three standard calibration solutions: 4.01 pH, 6.86 pH and 9.18 pH.

Programming the parameters

- P1** Relay control unit and analog output, pH or mV.
- P2** pH/mV value corresponding to 4 mA at output. May be set between -1.00 pH and 15.00 pH or between -1999 mV and +1999 mV.
- P3** pH/mV value corresponding to 20 mA at output.. May be set between -1.00 pH and 15.00 pH or between -1999 mV and +1999 mV.
- P4** Delay time in the intervention of relay A. May be set between 0 and 255 seconds.
- P5** Delay time in the intervention of relay B. May be set between 0 and 255 seconds.
- P6** Calibration of Pt100 probe, calibration of output in current, calibration of output in voltage. (**Calibration procedure to be carried out at a laboratory by skilled personnel**).
- P7** Display of the offset voltage value and of the slope value of the electrode.

To change one of these parameters press key PRG until the message corresponding to the parameter to be changed appears on the screen. Using the ▲ and ▼ keys, bring the parameter displayed to the desired value. Press OK to confirm.

Parameter P7 cannot be altered.

Calibrating the voltage input (calibration procedure to be carried out at a laboratory by skilled workers)

- Press the PRG key until the message **P6** appears on the display.
- Press the **CAL key four times**; the message CAL appears at the top of the display and the mV value of the input at the bottom.
- Simulate a voltage of 0 mV at the input (if the value is between ±25 mV the zero is calibrated, otherwise the full scale value is calibrated).
- Using the ▲ and ▼ keys, adjust the voltage value so as to have the correct voltage value on the display.
- Press the **SET** button, the **ON** symbol appears on the display to indicate that the instruments is measuring the voltage of the input using the second scale of measurement.
- Using the ▲ and ▼ keys adjust the voltage value so as to have the correct voltage value on the display.
- Press the set button, on the display the symbol **ON** disappears
- Simulate a voltage of 450 mV at the input, corresponding to the full value of the first scale.
- Using the ▲ and ▼ keys, adjust the voltage value so as to have the correct voltage value on the display.




Fig.1 Active transmitter.




Fig.2 Passive transmitter.





Fig.3



Dimensions

- Simulate a voltage of 1800 mV at the input, corresponding to the full value of the second scale.
- Using the ▲ and ▼ keys, adjust the voltage value so as to have the correct voltage value on the display.
- Press OK to confirm.

Pt100 probe calibration (100 Ω at 0°C) (calibration procedure to be carried out at a laboratory by skilled workers)

- Connect the Pt100 probe to the instrument. Press the PRG key until the message **P6** appears on the display.
- Press the CAL key; the message CAL appears at the bottom of the display and the temperature is shown at the top.
- Immerse the Pt100 probe and a precision thermometer for reference in the zero calibration bath. Wait long enough for the reading to become stable.
- Using the ▲ and ▼ keys, adjust the value of the temperature measured by the Pt100 probe so that it corresponds with the value on the reference thermometer.
- Immerse the Pt100 probe and a precision thermometer in the full scale calibration bath. Wait long enough for the reading to become stable.
- Using the ▲ and ▼ keys, adjust the value of the temperature measured by the Pt100 probe so that it corresponds with the value on the reference thermometer.
- Press OK to confirm.

N.B.: If the temperature shown by the instrument is between $\pm 12^\circ\text{C}$, the instrument calibrates the probe offset, otherwise it calibrates the gain.

Calibrating the analog output (calibration procedure to be carried out at a laboratory by skilled workers)

- Press the PRG key until the message **P6** appears on the display.
- Connect a precision milliammeter to the analog output.
- Press the **CAL key twice**; the message CAL appears at the top of the display and the message 4.0 at the bottom, indicating calibration at 4 mA.
- Using the ▲ and ▼ keys, adjust the value of the output current so as to have an indication of 4.0 mA on the precision milliammeter.
- Press the CAL key; the message CAL appears at the top of the display and the message 20.0 at the bottom, indicating calibration at 20 mA.
- Using the ▲ and ▼ keys, adjust the value of the output current so as to have an indication of 20.0 mA on the precision milliammeter.
- Press OK to confirm.

Display

Symbol	description
°C	the value shown is in °C.
°F	the value shown is in °F.
pH	the unit of the value shown is pH.
mV	the unit of the value shown is milli Volts.
A	the relay A is in closed status.
B	the relay B is in closed status.
ON	the value shown corresponds to the closing thresholds of the contacts of relay A or B.
OFF	the value shown corresponds to the opening thresholds of the contacts of relay A or B.

Error signal

- OFL** - Warning which appears during measurement when the value to be displayed is out of scale.
- E1** - Error warning which appears during pH calibration to indicate that the offset value of the electrode is too high in absolute value.
- E2** - Error warning which appears during pH calibration to indicate that the difference between the mV readings given by the two buffer solutions used for calibration is too great.
- E3** - Error warning which appears during pH calibration to indicate that the mV readings given by the two buffer solutions used for calibration are too close (about 50 mV at 25°C).
- E4** - Reading error on the EEPROM.
- E5** - Error warning indicating that the slope calculation gives a value 20% lower than the nominal value or gives a negative value.
- E6** - Error warning indicating that the slope calculation gives a value 150% higher than the nominal value.

Order code

- DO 9785T**: pH transmitter 4÷20 mA passive or active, power supply 24 Vac with double display 96x96 mm, for panel mounting.
- DO 9765T**: pH transmitter 4÷20 mA passive or active, power supply 24 Vac with double display 122x120, for use on the field.
- HD 882 M100/300**: Temperature probe with Pt100 sensor, miniature head, shaft Ø6x300 mm.
- HD 882 M100/600**: Temperature probe with Pt100 sensor, miniature head, shaft Ø6x600 mm.
- HD 8642**: Buffer solution 4.01 pH.
- HD 8672**: Buffer solution 6.86 pH.
- HD 8692**: Buffer solution 9.18 pH.
- HD R220**: Buffer solution Redox 220 mV 0,5l.
- HD R468**: Buffer solution Redox 468 mV 0,5l.

WA-18

HD62PT: Diaphragm cleaning (tiourea in HCl) - 500ml.

HD62PP: Protein cleaning (pepsin in HCl) - 500ml.

HD62RF: Regeneration (fluorhydric acid) - 100ml.

HD62SC: Solution for electrode preservation - 200ml.

CP5: Extension cable 5m. Connector BNC/S7.

CP10: Extension cable 10m. Connector BNC/S7.

CP15: Extension cable 15m. Connector BNC/S7.

CP5T: Extension cable 5m. Connector S7/wire.


CP5/10T: Extension cable 10m. Connector S7/wire.

KPI 10: Combined industrial electrode, S7 PG13.5 connector, gel, glass body, Ag/AgCl sat KCl Ø12x120 mm, temperature 0÷130°C, porous Teflon fitting.


KPI 11: Combined industrial electrode, S7 brass 3/4" connector, Rytron body, Ag/AgCl sat KCl, temperature 0÷100°C, porous Teflon fitting.

KPI 12: Platinum electrode for Redox measurement, S7 PG13,5 connector, Ag/AgCl sat KCl, temperature 0÷130°C, pressure 6 bar.

KPI 13: Platinum electrode for Redox measurement, Rytron body, S7 PG13,5 connector, temperature 0÷100°C, Ag/AgCl sat KCl.



KPI 10 0...14 pH / **KPI 12 Redox** ± 1999 mV, 0...130°C



KPI 11 0...14 pH / **KPI 13 Redox** ± 1999 mV, 0...100°C





HD 2106.1, HD 2106.2 CONDUCTIVITY METERS - THERMOMETERS

The **HD2106.1** and **HD2106.2** are portable instruments with a large LCD display. They measure conductivity, liquid resistivity, total dissolved solids (TDS), and salinity using combined 4-ring and 2-ring conductivity/temperature probes. Temperature only is measured by Pt100 or Pt1000 immersion, penetration, contact or air probes. The probe calibration can be performed automatically in one or more than one of the 147 μ S, 1413 μ S, 12880 μ S or 111800 μ S/cm conductivity calibration solutions. The temperature probes are equipped with an automatic recognition module and factory calibration data are stored inside. The HD2106.2 is a **datalogger**. It memorizes up to 36,000 conductivity and temperature samples which can be transferred from the instrument connected to a PC via the RS232C and USB 2.0 serial ports. The storing interval, printing, and baud rate can be configured using the menu. Both models are fitted with an RS232C serial port and can transfer to a PC the acquired measurements or to a portable printer in real time. The *Max*, *Min* and *Avg* function calculates the maximum, minimum or average values. Other functions include: the relative measurement REL, the Auto-HOLD function, and the automatic turning off which can also be excluded.

The instruments have IP66 protection degree.



HD40.1



SWD10

INSTRUMENT TECHNICAL CHARACTERISTICS

Measured quantities: χ , Ω , TDS, NaCl, $^{\circ}\text{C}$, $^{\circ}\text{F}$

Instrument

Dimensions (Length x Width x Height)	185x90x40mm
Weight	470g (complete with batteries)
Materials	ABS, rubber
Display	2x4½ digits plus symbols Visible area: 52x42mm

Operating conditions

Working temperature	-5...50 $^{\circ}\text{C}$
Storage temperature	-25...65 $^{\circ}\text{C}$
Working relative humidity	0...90%RH without condensation
Protection degree	IP66

Power

Batteries	4 1.5V type AA batteries
Autonomy	200 hours with 1800mAh alkaline batteries
Power absorbed with instrument off	20 μ A
Mains (SWD10)	Output mains adapter 12Vdc / 1A

Security of memorized data

Unlimited, independent of battery charge conditions

Time

Date and time	In real time
Accuracy	1min/month max error

Measured values storage - model HD2106.2

Type	2000 pages containing 18 samples each
Quantity	36000 pairs of measurements [χ - $^{\circ}\text{C}$], [Ω - $^{\circ}\text{C}$], [TDS- $^{\circ}\text{C}$] or [Sal- $^{\circ}\text{C}$]
Selectable storage interval	1s, 5s, 10s, 15s, 30s, 1min, 2min, 5min, 10min, 15min, 20min, 30min and 1hour

Serial interface RS232C

Type	RS232C electrically isolated
Baud rate	Can be set from 1200 to 38400 baud
Data bit	8
Parity	None
Stop bit	1
Flow Control	Xon/Xoff
Serial cable length	Max 15m
Print interval	Immediate or selectable between: 1s, 5s, 10s, 15s, 30s, 1min, 2min, 5min, 10min, 15min, 20min, 30min and 1hour

USB interface - model HD2106.2

Type	1.1 - 2.0 electrically isolated
------	---------------------------------

Connections

Conductivity input	8-pole male DIN45326 connector
Input module for the temperature probes	8-pole male DIN45326 connector
Serial interface and USB	8-pole MiniUSB type B
Mains adapter	2-pole connector (positive at centre)

Measurement of conductivity

		Resolution
Measuring range Kcell=0.01	0.000...1.999 μ S/cm	0.001 μ S/cm
Measuring range Kcell=0.1	0.00...19.99 μ S/cm	0.01 μ S/cm
Measuring range Kcell=1	0.0...199.9 μ S/cm	0.1 μ S/cm
	200...1999 μ S/cm	1 μ S/cm
	2.00...19.99mS/cm	0.01mS/cm
	20.0...199.9mS/cm	0.1mS/cm
Measuring range Kcell=10	200...1999mS/cm	1mS/cm

Accuracy (conductivity)	$\pm 0.5\% \pm 1$ digit
-------------------------	-------------------------

Measurement of resistivity

		Resolution
Measuring range Kcell = 0.01	till 100G Ω -cm/(*)	
Measuring range Kcell=0.1	till 100M Ω -cm/(*)	
Measuring range Kcell	5.0...199.9 Ω -cm	0.1 Ω -cm
	200...999 Ω -cm	1 Ω -cm
	1.00k...19.99k Ω -cm	0.01k Ω -cm
	20.0k...99.9k Ω -cm	0.1k Ω -cm
	100k...999k Ω -cm	1k Ω -cm
	1...10M Ω -cm	1M Ω -cm
Measuring range Kcell=10	0.5...5.0 Ω -cm	0.1 Ω -cm

Accuracy (resistivity)	$\pm 0.5\% \pm 1$ digit
------------------------	-------------------------

Measurement of total dissolved solids (with coefficient χ /TDS=0.5)

Measuring range Kcell=0.01	0.000...19.999mg/l	0.005mg/l
Measuring range Kcell=0.1	0.00...19.99mg/l	0.05mg/l
Measuring range Kcell=1	0.0...199.9mg/l	0.5mg/l
	200...1999mg/l	1mg/l
	2.00...19.99g/l	0.01g/l
	20.0...99.9g/l	0.1g/l
Measuring range Kcell=10	100...999g/l	1g/l

Accuracy
(total dissolved solids) $\pm 0.5\% \pm 1$ digit

Measurement of salinity

Measurement range	0.000...1.999g/l	Resolution 1mg/l
	2.00...19.99g/l	10mg/l
	20.0...199.9g/l	0.1g/l

Accuracy (salinity) $\pm 0.5\% \pm 1$ digit

**Temperature compensation
automatic/manual**

0...100°C with α_T selectable from 0.00 to 4.00%/°C

Reference temperature 20°C or 25°C

χ / TDS Conversion factor 0.4...0.8

Preset cell constant values: K=0,01 - K=0,1 - K=0,7 - K=1 - K=10

**Standard solutions automatically
detected @25°C**

147µS/cm
1413µS/cm
12880µS/cm
111800µS/cm

Measurement of temperature

Pt100 measuring range	-50...+200°C
Pt1000 measuring range	-50...+200°C
Resolution	0.1°C
Accuracy	$\pm 0.5\% \pm 1$ digit
Drift after 1 year	0.1°C/year

(*) The resistivity measurement is obtained from the reciprocal of conductivity measurement. Close to the bottom of the scale, the indication of resistivity appears like reported in the table below:

K cell = 0.01 cm ⁻¹		K cell = 0.1 cm ⁻¹	
Conductivity (µS/cm)	Resistivity (MΩ·cm)	Conductivity (µS/cm)	Resistivity (MΩ·cm)
0.001 µS/cm	1000 MΩ·cm	0.01 µS/cm	100 MΩ·cm
0.002 µS/cm	500 MΩ·cm	0.02 µS/cm	50 MΩ·cm
0.003 µS/cm	333 MΩ·cm	0.03 µS/cm	33 MΩ·cm
0.004 µS/cm	250 MΩ·cm	0.04 µS/cm	25 MΩ·cm



TECHNICAL DATA OF PROBES AND MODULES EQUIPPED WITH INSTRUMENT **Temperature probes Pt100 sensor with SICRAM module**

Model	Type	Application field	Accuracy
TP472I	Immersion	-196°C...+500°C	±0.25°C (-196°C...+300°C) ±0.5°C (+300°C...+500°C)
TP472I.0 1/3 DIN Thin Film	Immersion	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP473PI	Penetration	-50°C...+400°C	±0.25°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP473P.0 1/3 DIN Thin Film	Penetration	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP474C.I	Contact	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP474C.0 1/3 DIN Thin Film	Contact	-50°C...+300°C	±0.3°C (-50°C...+300°C)
TP475A.0 1/3 DIN Thin Film	Air	-50°C...+250°C	±0.3°C (-50°C...+250°C)
TP472I.5	Penetration	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP472I.10	Penetration	-50°C...+400°C	±0.30°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP49A.0 Class A Thin Film	Immersion	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AC.0 Class A Thin Film	Contact	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AP.0 Class A Thin Film	Penetration	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP875.I	Globe-thermometer Ø150mm	-30°C...+120°C	±0.25°C
TP876.I	Globe-thermometer Ø50mm	-30°C...+120°C	±0.25°C
TP87.0 1/3 DIN Thin Film	Immersion	-50°C...+200°C	±0.25°C
TP878.0 1/3 DIN Thin Film TP878.1.0 1/3 DIN Thin Film	Photovoltaic	+4°C...+85°C	±0.25°C
TP879.0 1/3 DIN Thin Film	Compost	-20°C...+120°C	±0.25°C

Common characteristics

Temperature drift @ 20°C 0.003%/°C

4 wires Pt100 and 2 wires Pt1000 Probes

Model	Type	Application field	Accuracy
TP47.100.0 1/3 DIN Thin Film	4 wires Pt100	-50...+250°C	1/3 DIN
TP47.1000.0 1/3 DIN Thin Film	2 wires Pt1000	-50...+250°C	1/3 DIN
TP87.100.0 1/3 DIN Thin Film	4 wires Pt100	-50...+200°C	1/3 DIN
TP87.1000.0 1/3 DIN Thin Film	2 wires Pt1000	-50...+200°C	1/3 DIN

Common features

Temperature drift @20°C

Pt100 0.003%/°C

Pt1000 0.005%/°C

A For the models of portable data logger series **HD21XX.2** has been implemented with a new serial port miniUSB type HID (Human Interface Device).

When making the connection to the PC by the USB cable Type A - Mini USB B-type coded CP23, **no USB driver installation is requested.**

B For the connection of the models **HD21XX.1** to the RS232 port of your PC, the USB/serial converter is available (**code C.206**). The converter is equipped with its own drivers that have to be installed before connecting the converter to the PC (please see the details in the CD Rom supplied with the converter).

C The port with the MiniDIN connector which is present on every model is an RS232C type. By means of the cable coded HD2110CSNM, an RS232 port of a PC or the HD40.1. printer can be connected.

ORDER CODES

HD2106.1: The kit is composed of: instrument HD2106.1, 4 1.5V alkaline batteries, operating manual, case and DeltaLog9 software.

HD2106.2: The kit is composed of: instrument HD2106.2 **datalogger**, 4 1.5V alkaline batteries, operating manual, case and DeltaLog9 software.
Conductivity probes, temperature probes, standard calibration solutions, cables for data transfer to PC or printer have to be ordered separately.

HD2110CSNM: 8-pole connection cable MiniDin - Sub D 9-pole female for RS232C.

C.206: Serial connection cable for HD2106.1 instruments with USB connector for PC and 8-pole MiniDin male connector for the instrument.

CP23: Serial connection cable with USB connector type A - MiniUSB type B (not suitable for HD2106.1).

DeltaLog9: Software for download and management of the data on PC using Windows operating systems.

SWD10: Stabilized power supply 100-240 Vac/12Vdc-1A mains voltage

HD40.1: 24-column portable thermal printer, serial interface, 57mm paper width, four NiMH 1.2V rechargeable batteries, SWD10 power supply, instruction manual, 5 thermal paper rolls It uses the HD2110CSNM cable (optional).

RCT: The kit includes 4 thermal paper rolls 57mm wide and 32mm in diameter.

BAT-40: Spare battery pack for HD40.1 printer with built-in temperature sensor.

HD22.2: Laboratory electrode holder composed of base plate with built-in magnetic stirrer, shaft and replaceable electrode holder. Suitable diameter 12mm. Powered by bench-top meters of the series **HD22...**with cable HD22.2.1 (**optional**) or power supplier SWD10 (**optional**).

HD22.3: Laboratory electrode holder composed of base plate. Flexible arm for free positioning. Suitable for electrodes with diameter 12mm.

Conductivity probes

Please see the order codes reported in the probes' technical specifications.


Standard conductivity calibration solutions

HD8747: Standard calibration solution 0.001mol/l equal to 147µS/cm @25°C, 200cc.

HD8714: Standard calibration solution 0.01mol/l equal to 1413µS/cm @25°C, 200cc.

HD8712: Standard calibration solution 0.1mol/l equal to 12880µS/cm @25°C, 200cc.

HD87111: Standard calibration solution 1mol/l equal to 111800µS/cm @25°C, 200cc.



Temperature probes equipped with SICRAM module

- TP472I:** Wire wound Pt100 sensor, immersion probe. Stem Ø 3 mm, length 300 mm. Cable length 2 m.
- TP472I.0:** Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 230 mm. Cable length 2 m.
- TP473PI:** Wire wound Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.
- TP473PO:** Thin film Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.
- TP474C.I:** Wire wound Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.
- TP474C.O:** Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.
- TP475A.O:** Thin film Pt100 sensor, air probe. Stem Ø 4mm, length 230mm. Cable length 2 m.
- TP472I.5:** Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 500 mm. Cable length 2 m.
- TP472I.10:** Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 1000mm. Cable length 2 m.
- TP49A.O:** Thin film Pt100 sensor, immersion probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle
- TP49AC.O:** Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 150mm. Cable length 2 m. Aluminium handle
- TP49AP.O:** Thin film Pt100 sensor, penetration probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle
- TP875.I:** Wire wound Pt100 sensor, 150mm diameter globe-thermometer equipped with handle. Cable length 2 m.
- TP876.I:** Wire wound Pt100 sensor, 50mm diameter globe-thermometer equipped with handle. Cable length 2 m.
- TP87.O:** Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 70 mm. Cable length 2 m.
- TP878.O:** Thin film Pt100 sensor, contact probe for solar panels. Cable length 2 m.

TP878.1.O: Thin film Pt100 sensor, contact probe for solar panels. Cable length 5 m.

TP879.O: Thin film Pt100 sensor, penetration probe for compost. Stem Ø 8 mm, length 1000 mm. Cable length 2 m.

Temperature probes without SICRAM module

TP47.100.O: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP47.1000.O: Thin film Pt1000 sensor, immersion probe. Probe's Stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.






TP47: Connector for Pt100 4-wire and Pt1000 2-wire probes without SICRAM module.

TP87.100.O: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 70mm. 4-wires connection cable with connector, length 1 m.

TP87.1000.O: Thin film Pt1000 sensor, immersion probe. Stem Ø 3mm, length 70mm. 2-wires connection cable with connector, length 1 m.

TECHNICAL DATA OF PROBES AND MODULES EQUIPPED WITH INSTRUMENT

2 and 4 electrode conductivity probes

ORDER CODE	MEASUREMENT RANGE	DIMENSIONS
SP06T	K=0.7 5µS...200mS/cm 0...90°C 4-electrode cell in Pocan/Platinum Max pressure 5bar	
SPT 401.001 not suitable for HD 2306.0	K=0.01 0,04...20µS/cm 0...120°C 2-electrode cell AISI 316 - Teflon Max pressure 5bar	
SPT01G	K=0.1 0.1µS...500µS/cm 0...80°C 2-electrode cell in Glass/Platinum Max pressure 5bar	
SPT1G	K=1 10µS...10mS/cm 0...80°C 2-electrode cell in Glass/Platinum Max pressure 5bar	
SPT10G	K=10 500µS...200mS/cm 0...80°C 2-electrode cell in Glass/Platinum Max pressure 5bar	



HD 2306.0 CONDUCTIVITY METER - THERMOMETER

The **HD2306.0** is a portable instrument with a large LCD display. It measures conductivity, liquid resistivity, and total dissolved solids (TDS), using combined 4-ring and 2-ring conductivity/temperature probes. Temperature only is measured by Pt100 or Pt1000 immersion, penetration, contact or air probes. The probe calibration can be performed automatically in one or more than one of the 147µS, 1413µS, 12880µS/cm or 111800µS/cm conductivity calibration solutions. The temperature probes are equipped with an automatic recognition module and factory calibration data are stored inside.

The *Max*, *Min* and *Avg* function calculates the maximum, minimum or average values.

Other functions include: the relative measurement REL, the Auto-HOLD function, and the automatic turning off which can also be disabled.

The instrument has IP67 protection degree.

INSTRUMENT TECHNICAL CHARACTERISTICS

Measured quantities: χ , Ω , TDS, °C, °F

Instrument

Dimensions

(Length x Width x Height)

Weight

Materials

Display

140x88x38mm

160g (complete with batteries)

ABS

2x4½ digits plus symbols

Visible area: 52x42mm

Operating conditions

Working temperature

Storage temperature

Working relative humidity

Protection degree

-5...50°C

-25...65°C

0...90%RH without condensation

IP67

Power

Batteries

Autonomy

Power absorbed with instrument off

3 1.5V type AA batteries

200 hours with 1800mAh alkaline batteries

< 20µA

Connections

Conductivity input/temperature probes

8-pole male DIN45326 connector

Measurement of conductivity

Measuring range Kcell=0.1

Measuring range Kcell=1

0.00...19.99µS/cm

0.0...199.9µS/cm

200...1999µS/cm

2.00...19.99mS/cm

20.0...199.9mS/cm

200...1999mS/cm

Measuring range Kcell=10

Accuracy (conductivity)

±0.5%±1digit

Resolution

0.01µS/cm

0.1µS/cm

1µS/cm

0.01mS/cm

0.1mS/cm

1mS/cm

Measurement of resistivity

Measuring range Kcell=0.1

Measuring range Kcell=1

till 100MΩ·cm/(*)

5.0...199.9Ω·cm

200...999Ω·cm

1.00k...19.99kΩ·cm

20.0k...99.9kΩ·cm

100k...999kΩ·cm

1...10MΩ·cm

Measuring range Kcell=10

Accuracy (resistivity)

0.5...5.0Ω·cm

±0.5%±1digit

0.1Ω·cm

1Ω·cm

0.01kΩ·cm

0.1kΩ·cm

1kΩ·cm

1MΩ·cm

0.1Ω·cm

Measurement of total dissolved solids (with coefficient χ /TDS=0.5)

Measuring range Kcell=0.1

Measuring range Kcell=1

0.00...19.99mg/l

0.0...199.9mg/l

200...1999mg/l

2.00...19.99g/l

20.0...99.9g/l

Measuring range Kcell=10

Accuracy (conductivity)

100...999g/l

±0.5%±1digit

0.05mg/l

0.5mg/l

1mg/l

0.01g/l

0.1g/l

1g/l

Measurement of temperature

Pt100 measuring range

Pt1000 measuring range

Resolution

Accuracy

Drift after 1 year

-50...+200°C

-50...+200°C

0.1°C

±0.1°C ±1digit

0.1°C/year

Temperature compensation automatic/manual

0...100°C with α_t selectable from 0.00 to 4.00%/°C

Reference temperature

20°C or 25°C

χ / TDS Conversion factor

0.4...0.8

Preset cell constant values:

K=0.1 - K=0.7 - K=1 - K=10

Standard solutions automatically detected @25°C

147µS/cm

1413µS/cm

12880µS/cm

111800µS/cm



(*) The resistivity measurement is obtained from the reciprocal of conductivity measurement. Close to the bottom of the scale, the indication of resistivity appears like reported in the table below:

K cell = 0.1 cm ⁻¹	
Conductivity (μS/cm)	Resistivity (MΩ-cm)
0.01 μS/cm	100 MΩ-cm
0.02 μS/cm	50 MΩ-cm
0.03 μS/cm	33 MΩ-cm
0.04 μS/cm	25 MΩ-cm

TECHNICAL DATA OF PROBES AND MODULES EQUIPPED WITH INSTRUMENT Temperature probes Pt100 sensor with SICRAM module

Model	Type	Application field	Accuracy
TP472I	Immersion	-196°C...+500°C	±0.25°C (-196°C...+300°C) ±0.5°C (+300°C...+500°C)
TP472I.0 1/3 DIN Thin Film	Immersion	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP473PI	Penetration	-50°C...+400°C	±0.25°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP473P.0 1/3 DIN Thin Film	Penetration	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP474C.I	Contact	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP474C.0 1/3 DIN Thin Film	Contact	-50°C...+300°C	±0.3°C (-50°C...+300°C)
TP475A.0 1/3 DIN Thin Film	Air	-50°C...+250°C	±0.3°C (-50°C...+250°C)
TP472I.5	Penetration	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP472I.10	Penetration	-50°C...+400°C	±0.30°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP49A.0 Class A Thin Film	Immersion	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AC.0 Class A Thin Film	Contact	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AP.0 Class A Thin Film	Penetration	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP875.I	Globe-thermometer Ø150mm	-30°C...+120°C	±0.25°C
TP876.I	Globe-thermometer Ø50mm	-30°C...+120°C	±0.25°C
TP87.0 1/3 DIN Thin Film	Immersion	-50°C...+200°C	±0.25°C
TP878.0 1/3 DIN Thin Film TP878.1.0 1/3 DIN Thin Film	Photovoltaic	+4°C...+85°C	±0.25°C
TP879.0 1/3 DIN Thin Film	Compost	-20°C...+120°C	±0.25°C

Common characteristics

Temperature drift @ 20°C 0.003%/°C

4 wires Pt100 and 2 wires Pt1000 Probes

Model	Type	Application field	Accuracy
TP47.100.0 1/3 DIN Thin Film	4 wires Pt100	-50...+250°C	1/3 DIN
TP47.1000.0 1/3 DIN Thin Film	2 wires Pt1000	-50...+250°C	1/3 DIN
TP87.100.0 1/3 DIN Thin Film	4 wires Pt100	-50...+200°C	1/3 DIN
TP87.1000.0 1/3 DIN Thin Film	2 wires Pt1000	-50...+200°C	1/3 DIN

Common features

Temperature drift @20°C

Pt100 0.003%/°C

Pt1000 0.005%/°C

ORDERING CODES

HD2306.0: The kit is composed of: instrument HD2306.0, 3 1.5V alkaline batteries, operating manual, case. **Other conductivity probes, temperature probes, calibration solutions must be ordered separately.**

HD22.2: Laboratory electrode holder composed of base plate with built-in magnetic stirrer, shaft and replaceable electrode holder. Suitable diameter 12mm. Powered by bench-top meters of the series HD22...with cable HD22.2.1 (optional) or power supplier SWD10 (optional).

HD22.3: Laboratory electrode holder composed of base plate. Flexible arm for free positioning. Suitable for electrodes with diameter 12mm.

Conductivity probes

Please see the order codes reported in the probes' technical specifications.

Standard conductivity calibration solutions

HD8747: Standard calibration solution 0.001mol/l equal to 147μS/cm @25°C, 200cc.

HD8714: Standard calibration solution 0.01mol/l equal to 1413μS/cm @25°C, 200cc.

HD8712: Standard calibration solution 0.1mol/l equal to 12880μS/cm @25°C, 200cc.

HD87111: Standard calibration solution 1mol/l equal to 111800μS/cm @25°C, 200cc.

Temperature probes equipped with SICRAM module

TP472I: Wire wound Pt100 sensor, immersion probe. Stem Ø 3 mm, length 300 mm. Cable length 2 m.

TP472I.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 230 mm. Cable length 2 m.

TP473PI: Wire wound Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP473P.0: Thin film Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP474C.I: Wire wound Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP474C.0: Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP475A.0: Thin film Pt100 sensor, air probe. Stem Ø 4mm, length 230mm. Cable length 2 m.

TP472I.5: Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 500 mm. Cable length 2 m.

TP472I.10: Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 1000mm. Cable length 2 m.

TP49A.0: Thin film Pt100 sensor, immersion probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AC.0: Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AP.0: Thin film Pt100 sensor, penetration probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP875.I: Wire wound Pt100 sensor, 150mm diameter globe-thermometer equipped with handle. Cable length 2 m.

TP876.I: Wire wound Pt100 sensor, 50mm diameter globe-thermometer equipped with handle. Cable length 2 m.

TP87.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 70 mm. Cable length 2 m.

TP878.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 2 m.

TP878.1.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 5 m.

TP879.0: Thin film Pt100 sensor, penetration probe for compost. Stem Ø 8 mm, length 1000 mm. Cable length 2 m.



HD8747





HD8714

HD8712

HD8711

Temperature probes without SICRAM module

- TP47.100.0:** Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 230mm.
Connection cable 4 wires with connector, length 2 m.
- TP47.1000.0:** Thin film Pt1000 sensor, immersion probe. Probe's Stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.
- TP47:** Connector for Pt100 4-wire and Pt1000 2-wire probes without SICRAM module.
- TP87.100.0:** Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 70mm.
4-wires connection cable with connector, length 1 m.
- TP87.1000.0:** Thin film Pt1000 sensor, immersion probe. Stem Ø 3mm, length 70mm. 2-wires connection cable with connector, length 1 m.

TECHNICAL DATA OF PROBES AND MODULES EQUIPPED WITH INSTRUMENT		
2 and 4 electrode conductivity probes		
ORDER CODE	MEASUREMENT RANGE	DIMENSIONS
SP06T	K=0.7 5µS...200mS/cm 0...90°C 4-electrode cell in Pocer/Platinum Max pressure 5bar	
SPT01G	K=0.1 0.1µS...500µS/cm 0...80°C 2-electrode cell in Glass/Platinum Max pressure 5bar	
SPT1G	K=1 10µS...10mS/cm 0...80°C 2-electrode cell in Glass/Platinum Max pressure 5bar	
SPT10G	K=10 500µS...200mS/cm 0...80°C 2-electrode cell in Glass/Platinum Max pressure 5bar	



DO 9786T - R1 • DO 9766T - R1

CONDUCTIVITY TRANSMITTERS

DO 9786T/DO 9766T transmitters convert the output of a conductivity electrode with temperature compensation into a 4÷20 mA signal. The electrode input circuit is galvanically insulated against the 4÷20 mA output signal. An LCD indicator allows viewing of the process signal value and of the various parameters. The accurate design and choice of components make the instrument precise and reliable for a long working life. The instrument works in conjunction with a conductivity electrode and a temperature probe (Pt100 sensor, 100 Ω at 0°C).

Technical characteristics

Input conductivity	Measuring range	0.0...199.9 mS
	2/4 electrodes	Configurable cell constant 0.01...199.9 cm ⁻¹
	Transducer energizing	Square wave 10...1000 mV, depending on conductivity, 200...1600 Hz, depending on conductivity
	Input impedance	>100 MΩ
	Cable length	<10 metres unscreened <50 metres screened (about 2 nF)
	Accuracy	0.5% of reading ±2 digits ±0.01% per °C of drift in temperature
Input temperature	Pt100 2/4 wires	-50...199.9°C
	Transducer energizing	0.5 mA dc
	Cable length	<10 metres unscreened <50 metres screened (about 5 nF)
	Accuracy	0.2°C ±0.1% of reading ±0.01°C/°C of drift in temperature
Compensation temperature	None	
	manual	Linear 0.00...4.00%/°C -50...+200°C
	automatic	Linear 0.00...4.00%/°C -50...+200°C
	Reference temperature	20 or 25°C Configurable
Current output	4.00...20.00 mA	Programmable and proportional to conductivity
	Accuracy	0.5% of reading ±0.02 mA
	Insulation	2500 Vac 1 minute
R Load	Load resistance	$R_{Lmax} = \frac{V_{dc}-10}{0,022}$ $R_{Lmax} = 636 \Omega @ V_{dc} = 24 V_{dc}$
Relay output	A and B	Bistable, contact 3A/230 Vac free potential
Power supply	Passive	4÷20 mA, 2 wire configuration, 10÷35 V, see fig. 2
	Active	24/230 Vac - 15/+10% 1 VA, 48...62 Hz, see fig. 1
Case DO 9766T	External dimensions	120x122x56 mm wall mounting
	Protection class	IP64
Case DO 9786T	External dimensions	96x96x126 mm panel mounting
	Protection class	IP44


Key functions

- PRG** Programming of the parameters is activated by pressing the PRG key plus the ▲ and ▼ keys. The message P1 appears on the display, indicating that the parameter P1 is being programmed. When the PRG key is pressed continuously, the messages P2, P3, P4, P5, P6, P7, P8 and the corresponding parameters are displayed in sequence. After P8 the instrument returns to normal function.
- SET** Key for setting the relay intervention threshold. The ON or OFF symbol appears on the display, indicating the switching on or off threshold of relay A or of relay B.
- °C/°F** - If this key is pressed it changes the temperature measuring unit to degrees Celsius or degrees Fahrenheit.
 - When pressed together with the CAL key it activates the manual temperature setting function.
 - If pressed during the conductivity calibration function it quits the calibration function without storing the calibration.
- X** When pressed together with the CAL key it activates the conductivity calibration function.
- OK** Confirms the programming parameters, or the relay SET values, and stores them.
- CAL** - When pressed together with the °C/°F key it activates the manual temperature setting function.
 - When pressed together with the X key it activates the conductivity calibration function.
 - Key used to confirm conductivity calibration and manual temperature calibration.
- ▲ - Key for increasing the value displayed in the parameter programming phase.
 - During the relay SET point programming phase.
 - During the calibration phase.
- ▼ - Key for decreasing the value displayed in the parameter programming phase.
 - During the relay SET point programming phase.
 - During the calibration phase.

Setting the relay SET point

- Press the SET button; the ON symbol appears on the display with the letter A to indicate that the value shown corresponds to the switching on threshold of relay A.
- To change this value press the ▲ and ▼ keys.
- Press SET; the OFF symbol appears with the letter A to indicate that the switching off threshold of relay A is being displayed.
- To change this value press the ▲ and ▼ keys.
- Press the SET button; the ON symbol appears on the display with the letter B to indicate that the value shown corresponds to the switching on threshold of relay B.
- To change this value press the ▲ and ▼ keys.
- Press SET; the OFF symbol appears with the letter B to indicate that the switching off threshold of relay B is being displayed.
- To change this value press the ▲ and ▼ keys.
- Press SET, the instrument stores the values and returns to normal function.

NOTE: During the SET point setting phase (symbols ON or OFF lit) the instrument returns to normal function if no key is pressed for 2 minutes.



Manual temperature setting

If the temperature probe is not connected or if the probe is broken the measuring unit °C or °F flashes. In this case it is possible to set the temperature compensation value manually.

- Press the CAL key and the °C/°F key together; the message CAL appears at the bottom of the display.
- Using the ▲ and ▼ keys, set the temperature value corresponding to the temperature of the liquid in which you wish to measure conductivity.
- Press CAL to confirm this value. The message CAL disappears.

Calibration of the DO 9786T-R1 / DO 9766T-R1 with conductivity probes

Calibration of the DO 9786T-R1 / DO 9766T-R1 transmitters with conductivity probes:

- Immerse the probe in the buffer solution used for calibration.
- Press the CAL key and the X key together; the message CAL appears at the top of the display.
- The instrument can automatically recognize two standard calibration solutions: a 0.1 molar solution of KCl and a 0.01 molar solution of KCl. The instrument proposes the conductivity value as a function of the measured temperature if the temperature probe is connected, or the manually set temperature.
- Using the ▲ and ▼ keys, adjust the conductivity value measured as a function of the liquid temperature.
- Press CAL to confirm this value. The message CAL disappears.

NOTE: If you want to quit without storing the new calibration, press the °C/°F key.

N.B.: Before calibrating the probe set a cell constant close to the cell constant of the probe that you wish to calibrate with key PRG, function P2. If the message E1 appears during calibration, the instrument is indicating that the probe gain is too high; quit programming (°C/°F button) and increase the value of the cell constant. Likewise, if E2 appears, the instrument is indicating that the probe gain is too low; quit calibration and decrease the cell constant. Repeat the calibration operation.

Programming the parameters

- P1 Temperature coefficient. May be set between 0 and 4.00%/°C (0 and 2.22%/°F).
- P2 Cell constant. May be set between 0.01 and 199.9.
- P3 Conductivity value corresponding to 4 mA at output. May be set between 0 and 199.9 mS.
- P4 Conductivity value corresponding to 20 mA at output. May be set between 0 and 199.9 mS.
- P5 Delay time in the intervention of relay A. May be set between 0 and 250 seconds.
- P6 Delay time in the intervention of relay B. May be set between 0 and 250 seconds.
- P7 Reference temperature of the conductivity measurement. May be set between the values 20.0 or 25.0°C.
- P8 Calibration of Pt100 probe and calibration of analog output in current (see Pt100 probe calibration and analog output calibration).

To change one of these parameters (except P8) press key PRG until the message corresponding to the parameter to be changed appears on the screen. Using the ▲ and ▼ keys, bring the parameter displayed to the desired value. Press OK to confirm.

Pt100 probe calibration

(Factory calibration to be performed by skilled workers)

- Connect the Pt100 probe to the instrument. Press the PRG key until the message P8 appears on the display.
- Press the CAL key; the message CAL appears at the bottom of the display and the temperature is shown at the top.
- Immerse the Pt100 probe and a precision thermometer for reference in the zero calibration bath. Wait long enough for the reading to become stable.
- Using the ▲ and ▼ keys, adjust the value of the temperature measured by the Pt100 probe so that it corresponds with the value on the reference thermometer.
- Immerse the Pt100 probe and a precision thermometer in the full scale calibration bath. Wait long enough for the reading to become stable.
- Using the ▲ and ▼ keys, adjust the value of the temperature measured by the Pt100 probe so that it corresponds with the value on the reference thermometer.
- Press OK to confirm.

N.B.: If the temperature shown by the instrument is between +12°C, the instrument calibrates the probe offset, otherwise it calibrates the gain.

Calibrating the analog output

(Factory calibration to be performed by skilled workers)

- Press the PRG key until the message P8 appears on the display.
- Connect a precision milliammeter to the analog output.
- Press the **CAL key twice**; the message CAL appears at the top of the display and the message 4.0 at the bottom, indicating calibration at 4 mA.
- Using the ▲ and ▼ keys, adjust the value of the output current so as to have an indication of 4,000 mA on the precision milliammeter.
- Press the CAL key; the message CAL appears at the top of the display and the message 20.0 at the bottom, indicating calibration at 20 mA.
- Using the ▲ and ▼ keys, adjust the value of the output current so as to have an indication of 20,000 mA on the precision milliammeter.
- Press OK to confirm.

Display

Symbol Description

°C indicates that the value shown is in °C.

°F indicates that the value shown is in °F.

µS indicates that the unit of the value shown is micro Siemens.

mS indicates that the unit of the value shown is milli Siemens.

A indicates that the relay A is in closed status.

B indicates that the relay B is in closed status.

ON indicates that the value shown corresponds to the closing thresholds of the contacts of relay A or B.

OFF indicates that the value shown corresponds to the opening thresholds of the contacts of relay A or B.




Fig.1 Active transmitter.




Fig.2 Passive transmitter.





Fig.3 Calibration analog output.



Dimensions

Error signals

- OFL** - Warning which appears during measurement when the value to be displayed is out of scale.
- E1** - Error warning which appears during conductivity calibration to indicate that the probe gain is too high. Press P2 to increase the cell constant value.
- E2** - Error warning which appears during conductivity calibration to indicate that the probe gain is too low. Press P2 to decrease the cell constant value.
- E3** - Error warning which appears to indicate that the instrument is unable to recognize the buffer solution used for automatic calibration. Press the ▲ or ▼ key to make this indication disappear.
- E4** - Reading error on the EEPROM.

APPENDIX

Table of compatibility between range and sensor

Conductivity Range	Nominal cell constant			
	0.01÷0.2	0.2÷2	2÷20	20÷199.9
0÷19.99 µS	√			
0÷199.9 µS	√	√		
0÷1999 µS	√	√	√	
0÷199.9 mS	√	√	√	√
0÷19.99 mS		√	√	√
0÷199.9 mS			√	√
0÷1999 mS				√

Temperature sensor

Temperature	Pt100	Temperature	Pt100
-50°C	80.25 Ω	100°C	138.50 Ω
-25°C	90.15 Ω	125°C	147.94 Ω
0°C	100.00 Ω	150°C	157.32 Ω
25°C	109.73 Ω	175°C	166.62 Ω
50°C	119.40 Ω	199°C	175.47 Ω
75°C	128.98 Ω		

Calculating the temperature coefficient of a solution

If the temperature coefficient of the solution is not known, it may be determined using the DO 9786T/DO 9766T.

- Set the temperature coefficient at 0.0%/°C (parameter P1).

The following measurements should be taken as close as possible to the work point, between 5°C and 70°C, for greater accuracy.

- Immerse the probe in the testing liquid. Wait for the measurement to become stable.
- Take note of the temperature and of the conductivity.
- Increase the solution temperature by at least 10°C.
- Take note of the temperature and of the conductivity.
- Calculate the temperature coefficient using the following formula:

$$\alpha = \frac{(G_x - G_y) \times 100\%}{G_y(T_x - 20) - G_x(T_y - 20)} \quad (\text{reference temperature } 20^\circ\text{C})$$

Where:

G_x conductivity at temperature T_x
G_y conductivity at temperature T_y

N.B.: if the reference temperature is 25°C, replace 20 with 25.

- Set the temperature coefficient with the value calculated as above (parameter P1).

Calibration of the instrument for measuring conductivity

The conductivity measurement depends strongly on the temperature of the liquid that is to be measured; this relationship must be considered during calibration.

Calibration of the instrument alone by means of a precision resistance

This is a sure and accurate method for calibrating the instrument alone, but it does not allow for the variations of the cell constant that may occur, nor of the state of efficiency and cleanness of the cell.

The precision resistance used for calibration will be chosen according to the scale that you want to calibrate. Typical values are the following:

Conductivity	Resistance
100,0 µS	10.000 Ω
500,0 µS	2.000 Ω
1000 µS	1.000 Ω
5000 µS	200 Ω
10,00 mS	100 Ω
50,00 mS	20 Ω
100,0 mS	10 Ω
500,0 mS	2 Ω
1000 mS	1 Ω

The precision resistance will be connected to the end of the cable that connects the probe to the instrument. This ensures greater accuracy of calibration. Disable the temperature compensation α_T during the calibration of the instrument with the precision resistances.

Calibration with standard solutions

In this case too, for the calibration of the instrument, cable and measuring probes in a standard solution, the greatest attention must be paid to the temperature of the solutions and the cleanness of the measuring cell. It is advised not to carry out calibration below 500 µS/cm. Solutions with low conductivity must be kept closed in their containers. Contact with the air increases their value due to the absorption of CO₂.

The regulations for the preparation of standard solutions with a base of KCl dissolved in water with a high degree of purity supply the method and percentages of KCl and water to be mixed. DELTA OHM supplies four solutions for calibration:

HD8747: Standard calibration solution 0.001 mol/l equal to 147 µS/cm @25°C, 200cc.

HD8714: Standard calibration solution 0.01 mol/l equal to 1413 µS/cm @25°C, 200cc.

HD8712: Standard calibration solution 0.1 mol/l equal to 12880 µS/cm @25°C, 200cc.

HD87111: Standard calibration solution 1 mol/l equal to 111800 µS/cm @25°C, 200cc.

Care and maintenance of the conductivity cell

In conductivity measurement systems in industrial plants, if the installation is correctly made, readings are generally reliable for a long time. The important thing is to carry out correct, programmed maintenance of the measuring cell.

Abrasion of the cable due to continued swinging movements must be avoided, as must the formation of deposits and scale on the cell which can change its geometrical structure.

The cell must always be immersed in the liquid that is to be measured. In the industrial field, measurements may range from highly pure water to sewage or water contaminated by corrosive substances.

It is good practice to check the compatibility of the materials of which the cell and the connecting cable are made with the liquid in which the measurement is to be taken. Check that there are no floating bodies, suspended granules that may be more or less conductive, or which could get stuck inside the cell, thus leading to incorrect measurements.

To clean the cell use detergents or substances suitable for the material of which the cell is made.

Selecting the cell constant and installation

The measurement range of the liquid to be examined determines the choice of the cell constant to be used.

Installation of the cell will vary according to the application. On the whole, the following points must be considered:

- Choose the correct cell and cell constant, suitable for the application.
- Use suitable materials, cable, cell, supports, so as to resist corrosion and the influence of atmospheric agents.
- The sensor/cell must be firmly fixed, in a place where they are easily accessible for maintenance.
- The liquid in which the sensor is immersed must be a representative part of the whole that is to be measured.
- There must be a moderate flow of liquid so that an updated sample of liquid arrives at the electrodes. Excessive movement or flow causes turbulence and air bubbles between the electrodes. As an air bubble is not conductive, it modifies the volume of the cell and changes the constant.
- Install the sensor in such a way that sludge or particles of material cannot be deposited inside it.
- If installed in containers where high currents are circulating, the conductivity cell may present measuring problems.
- The maintenance and cleaning interval depends on the quality of the liquid in which the cell is installed.

ORDER CODE

DO 9786T-R1: Conductivity transmitter 4÷20 mA passive or active, power supply 24 Vac with double display 96x96 mm, **for panel mounting.**

DO 9766T-R1: Conductivity transmitter 4÷20 mA passive or active, power supply 24 Vac with double display 122x120 mm, **for use on the field.**

SPT 86: Combined industrial conductivity and temperature probe in POCAN with 4 platinum electrodes, cell constant K = 0.7, 1.5 meters cable, Pt100 with 4 wires. Temperature 0÷90°C.

SPTKI 10: Combined industrial conductivity probe in Glass with 2 platinum electrodes, cell constant K = 1, **S7/PG13** screw-joint, 2 wires output: eurostandard S7. Temperature 0÷100°C.

SPTKI 11: Combined industrial conductivity and temperature probe in Rytron with 2 electrodes, cell constant K = 1, 5 meters cable, Pt100 with four wires. Temperature 0÷80°C.

SPTKI 12: Combined industrial conductivity and temperature probe in Rytron with 2 platinum electrodes, cell constant K = 01, 5 meters cable, Pt100 with four wires. Temperature 0÷80°C.

SPTKI 13: Combined industrial conductivity and temperature probe in Rytron with 2 platinum electrodes, cell constant K = 10, 5 meters cable, Pt100 with four wires. Temperature 0÷80°C.

SPT 401.001: Combined conductivity and temperature probe with 2-electrode in AISI316 stainless steel, cell constant 0.01. Measuring range 0.04 ... 20 µS, 0 ... 120 °C.

HD 882 M100/300: Temperature probe with Pt100 sensor, miniature head, shaft Ø6x300 mm.

HD 882 M100/600: Temperature probe with Pt100 sensor, miniature head, shaft Ø6x600 mm.

HD 8747: Calibration solution 0.001 mol/l corresponding to 147 µS/cm at 25°C, 200cc.

HD 8712: Calibration solution 0.1 mol/l corresponding to 12880 µS/cm at 25°C, 200cc.

HD 8714: Calibration solution 0.01 mol/l corresponding to 1413 µS/cm at 25°C, 200cc.


HD 87111: Calibration solution 1 mol/l corresponding to 111800 µS/cm at 25°C, 200cc.

CP 5T: Extension cable to connect the probe SPT KI 10 to the transmitter, L=5m (S7 on one side, wires on the other side)


CP 5/10T: Extension cable to connect the probe SPT KI 10 to the transmitter, L=10m (S7 on one side, wires on the other side)

CP 5/20T: Extension cable to connect the probe SPT KI 10 to the transmitter, L=20m (S7 on one side, wires on the other side)


SPT 401.001
cell constant K=0.01




SPT 86
cell constant K=0.7




SPTKI 10
cell constant K=1




SPTKI 11
cell constant K=1



SPTKI 12
cell constant K=0.1



SPTKI 13
cell constant K=10



	Cell constant	Measuring range	Temperature range	Material	Electrodes	Temperature sensor	Max. pressure	Connection
SPT 401.001	K=0.01	0.04÷20µS/cm	0÷120°C	AISI 316 - PTFE	2 AISI 316	-	12bar	4-pole connector
SPT 86	K=0.7	5µS÷20mS	0÷90°C	Pocan	4 platinum	Pt100 4 wire	6bar	1.5 m cable
SPTKI 10	K=1	100µS÷200mS	0÷100°C	Glass	2 platinum	-	6bar	S7
SPTKI 11	K=1	100µS÷10mS	0÷80°C	Rytron	2 platinum	Pt100 4 wire	6bar	5 m cable
SPTKI 12	K=0.1	1µS÷1mS	0÷80°C	Rytron	2 platinum	Pt100 4 wire	6bar	5 m cable
SPTKI 13	K=10	10µS÷200mS	0÷80°C	Rytron	2 platinum	Pt100 4 wire	6bar	5 m cable



HD 2109.1, HD 2109.2 DISSOLVED OXYGEN - TEMPERATURE METERS

The **HD2109.1** and **HD2109.2** are portable instruments with a large LCD display. They measure the concentration (in mg/l) of dissolved Oxygen in liquids, the saturation index (in %) and the temperature using SICRAM combined probes of polarographic and galvanic type with two or three electrodes and integrated temperature sensor. Temperature only is measured by Pt100-SICRAM or direct 4 wire-immersion, penetration, contact or air probes. Thanks to an internal pressure sensor, the instruments automatically compensate for barometric pressure. The instrument anticipates automatic compensation of the Oxygen probe membrane permeability and of the salinity of the liquid being examined. The dissolved Oxygen probe's quick calibration function guarantees timely correctness of the performed measurements. The dissolved Oxygen and the temperature probes are equipped with an automatic recognition module and factory calibration data are stored inside. The HD2109.2 is a **datalogger**. It stores up to 18,000 dissolved Oxygen concentration, saturation index measurements, barometric pressure and temperature samples which can be transferred from the instrument connected to a PC via the RS232C and USB 2.0 serial ports. The storing interval, printing, and baud rate can be configured using the menu. Both models are fitted with an RS232C serial port and can transfer to a PC the acquired measurements or to a portable printer in real time. *The Max, Min and Avg* function calculates the maximum, minimum or average values. Other functions include: the relative measurement REL, the Auto-HOLD function, and the automatic turning off which can also be excluded.

The instruments have IP66 protection degree.



WA-34



INSTRUMENT TECHNICAL CHARACTERISTICS

Measured quantities: mg/l O₂, sat.% O₂, mbar, °C, °F

Instrument

Dimensions (Length x Width x Height)	185x90x40mm
Weight	470g (complete with batteries)
Materials	ABS, rubber
Display	2x4½ digits plus symbols Visible area: 52x42mm

Operating conditions

Working temperature	-5...50°C
Storage temperature	-25...65°C
Working relative humidity	0...90%RH without condensation
Protection degree	IP66

Power

Batteries	4 1.5V type AA batteries
Autonomy	200 hours with 1800mAh alkaline batteries
Power absorbed with instrument off	20µA
With dissolved oxygen probe	40µA
Mains (SWD10)	Output mains adapter 12Vdc / 1A

Security of memorized data

Unlimited, independent of battery charge conditions

Time

Date and time	Schedule in real time
Accuracy	1min/month max error

Measured values storage - model HD2109.2

Type	2000 pages containing 9 samples each
Quantity	18,000 samples composed of 4 parameters: mg/l O ₂ - %O ₂ - mbar - (°C or °F)
Selectable storage interval	1s, 5s, 10s, 15s, 30s, 1min, 2min, 5min, 10min, 15min, 20min, 30min and 1hour

Serial interface RS232C

Type	RS232C electrically isolated
Baud rate	Can be set from 1200 to 38400 baud
Data bit	8
Parity	None
Stop bit	1
Flow Control	Xon/Xoff
Serial cable length	Max 15m
Print interval	Immediate or selectable between: 1s, 5s, 10s, 15s, 30s, 1min, 2min, 5min, 10min, 15min, 20min, 30min and 1hour

USB interface - model HD2109.2

Type	1.1 - 2.0 electrically isolated
------	---------------------------------

Connections

Input for Oxygen and temperature probes	8-pole male DIN45326 connector
RS232C serial interface	8-pole MiniDin connector
USB interface	MiniUSB type B
Mains adapter	2-pole connector (positive at centre)

Measurement of the concentration of dissolved Oxygen

Measurement range	0.00...90.00mg/l
Resolution	0.01mg/l
Accuracy (0...90% 1013mbar, 20...25°C)	±0.03mg/l±1digit



HD40.1



SWD10

Measurement of the saturation index of dissolved Oxygen

Measurement range	0.0...600.0%
Resolution	0.1%
Accuracy	±0.3% ±1digit (in the range 0.0...199.9%) ±1% ±1digit (in the range 200.0...600.0%)

Measurement of barometric pressure

Measurement range	0.0...1100.0mbar
Resolution	0.1mbar
Accuracy	±2mbar±1digit between 18 and 25°C ±(2mbar+0.1mbar/°C) in the remaining range

Setting the salinity

Setting range	0.0...70.0g/l
Resolution	0.1g/l

Temperature measurement with the sensor inside the dissolved Oxygen probe

Measurement range	0...+45°C
Resolution	0.1°C
Accuracy	±0.1°C ±1digit
Drift after 1 year	0.1°C/year

Temperature measurement by Instrument with Pt100 probe

Pt100 measurement range	-200...+650°C
Resolution	0.1°C
Accuracy	±0.1°C ±1 digit
Drift after 1 year	0.1°C/year

Temperature compensation

Automatic	0...50°C
-----------	----------

TECHNICAL DATA OF PROBES AND MODULES EQUIPPED WITH INSTRUMENT

Temperature probes Pt100 sensor with SICRAM module

Model	Type	Application field	Accuracy
TP472I	Immersion	-196°C...+500°C	±0.25°C (-196°C...+300°C) ±0.5°C (+300°C...+500°C)
TP472I.0 1/3 DIN Thin Film	Immersion	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP473P.I	Penetration	-50°C...+400°C	±0.25°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP473P.0 1/3 DIN Thin Film	Penetration	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP474C.I	Contact	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP474C.0 1/3 DIN Thin Film	Contact	-50°C...+300°C	±0.3°C (-50°C...+300°C)
TP475A.0 1/3 DIN Thin Film	Air	-50°C...+250°C	±0.3°C (-50°C...+250°C)
TP472I.5	Penetration	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP472I.10	Penetration	-50°C...+400°C	±0.30°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP49A.0 Class A Thin Film	Immersion	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AC.0 Class A Thin Film	Contact	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AP.0 Class A Thin Film	Penetration	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP875.I	Globe-thermometer Ø150mm	-30°C...+120°C	±0.25°C
TP876.I	Globe-thermometer Ø50mm	-30°C...+120°C	±0.25°C
TP87.0 1/3 DIN Thin Film	Immersion	-50°C...+200°C	±0.25°C
TP878.0 1/3 DIN Thin Film	Photovoltaic	+4°C...+85°C	±0.25°C
TP878.1.0 1/3 DIN Thin Film			
TP879.0 1/3 DIN Thin Film	Compost	-20°C...+120°C	±0.25°C

Common characteristics

Temperature drift @ 20°C	0.003%/°C
--------------------------	-----------


4 wires Pt100 Probes

Model	Type	Application field	Accuracy
TP47.100.0 1/3 DIN Thin Film	4 wires Pt100	-50...+250°C	1/3 DIN
TP87.100.0 1/3 DIN Thin Film	4 wires Pt100	-50...+200°C	1/3 DIN

Common features

Temperature drift @20°C
Pt100

0.003%/°C



A For the models of portable data logger series **HD21XX.2** has been implemented with a new serial port miniUSB type HID (Human Interface Device).

When making the connection to the PC by the USB cable Type A - Mini USB B-type coded CP23, **no USB driver installation is requested.**

B For the connection of the models **HD21XX.1** to the RS232 port of your PC, the USB/serial converter is available (**code C.206**). The converter is equipped with its own drivers that have to be installed **before** connecting the converter to the PC (please see the details in the CDROM supplied with the converter).

C The port with the MiniDIN connector which is present on every model is an RS232C type. By means of the cable coded HD2110CSNM, an RS232 port of a PC or the HD40.1. printer can be connected.




DO9700

DO9701

Oxygen probe – dimensions and characteristics

Model	D09709SS	D09709SS.5	D09709SS.1	D09709SS.5.1
Type	Polarographic probe, Silver anode, Platinum cathode		Galvanic probe, Zinc anode, Silver cathode	
O ₂ measuring range	0.00...60.00mg/l		0.00...20.00mg/l	
Functioning temperature	0...45°C		0...50°C	
Accuracy instrument with probe	±1% f.s.		±2% f.s.	
Membrane	Replaceable		Replaceable	
Cable length	2m	5m(*)	2m	5m(*)

(*) Cable with connector



ORDER CODES

HD2109.1: The kit is composed of: instrument HD2109.1, calibrator D09709/20 (for polarographic probe) or D09709/21 (for galvanic probe), 4 1.5V alkaline batteries, operating manual, case and DeltaLog9 software.

The probes and data transfer cable must be ordered separately.

HD2109.2: The kit is composed of: instrument HD2109.2 **data logger**, calibrator D09709/20 (for polarographic probe) or D09709/21 (for galvanic probe), 4 1.5V alkaline batteries, operating manual, case and DeltaLog9 software.

The probes and data transfer cable must be ordered separately.

HD2110CSNM: 8-pole connection cable MiniDin - Sub D 9-pole female for RS232C.

C.206: Cable for instruments of the series HD21...1 for direct connection to the USB input of a PC.

CP23: USB 2.0 connection cable type A - MiniUSB type B (not suitable for HD2109.1).

DeltaLog9: Software for download and management of the data on PC using Windows operating systems.

SWD10: Stabilized power supply 100-240 Vac/12Vdc-1A mains voltage

HD40.1: 24-column portable thermal printer, serial interface, 57mm paper width, four NiMH 1.2V rechargeable batteries, SWD10 power supply, instruction manual, 5 thermal paper rolls.

RCT: The kit includes 4 thermal paper rolls 57mm wide and 32mm in diameter.

BAT-40: Spare battery pack for HD40.1 printer with built-in temperature sensor.

HD22.2: Laboratory electrode holder composed of base plate with built-in magnetic stirrer, shaft and replaceable electrode holder. Suitable diameter 12mm. Powered by bench-top meters of the series HD22...with cable HD22.2.1 (**optional**) or power supplier SWD10 (**optional**).

HD22.3: Laboratory electrode holder composed of base plate. Flexible arm for free positioning. Suitable for electrodes with diameter 12mm.

Solutions

D09700: zero oxygen solution.

D09701: electrolyte solution for polarographic probes D09709 SS and D09709 SS.5.

D09701.1: electrolyte solution for galvanic probes D09709 SS.1 and D09709 SS.5.1.

Combined dissolved Oxygen/temperature probes

DO 9709 SS Polarographic combined oxygen and temperature probe with possibility of membrane replacement. Ø12mm x 120mm. 2m cable. The code includes: probe, 2 membranes, electrolyte solution and zero point solution.

DO 9709 SS.5 Polarographic combined oxygen and temperature probe with possibility of membrane replacement. Ø12mm x 120mm. 5m cable. The code includes: probe, 2 membranes, electrolyte solution and zero point solution.

DO 9709 SS.1 Galvanic combined galvanic oxygen and temperature probe with possibility of membrane replacement. Ø12mm x 76mm. Ø16mm tip with membrane. 2m cable. The code includes: probe, 2 membranes in total, electrolyte solution and zero point solution.

DO 9709 SS.5.1 Galvanic combined galvanic oxygen and temperature probe with possibility of membrane replacement. Ø12mm x 76mm. Ø16mm tip with membrane. 5m cable. The code includes: probe, 2 membranes in total, electrolyte solution and zero point solution.

Accessories

DO 9709/20: Calibrator for polarographic probes DO 9709SS and DO 9709SS.5

DO 9709/21: Calibrator for galvanic probes DO 9709SS.1 and DO 9709SS.5.1

DO 9709 SSK: Kit of accessories for probes DO 9709SS and DO 9709SS.5: 3 membranes, zero point solution and electrolyte.

DO 9709/21K: Kit of accessories for probes DO 9709SS.1 and DO 9709SS.5.1: 3 membranes, zero point solution and electrolyte.

Temperature probes equipped with SICRAM module

TP472I: Wire wound Pt100 sensor, immersion probe. Stem Ø 3 mm, length 300 mm. Cable length 2 m.

TP472I.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 230 mm. Cable length 2 m.

TP473P.I: Wire wound Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP473P.O: Thin film Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP474C.I: Wire wound Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP474C.O: Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP475A.O: Thin film Pt100 sensor, air probe. Stem Ø 4mm, length 230mm. Cable length 2 m.

TP472I.5: Thin film Pt100 sensor, immersion probe. Stem Ø 6mm, length 500 mm. Cable length 2 m.

TP472I.10: Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 1000mm. Cable length 2 m.

TP49A.O: Thin film Pt100 sensor, penetration probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AC.O: Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AP.O: Thin film Pt100 sensor, penetration probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP875.I: Wire wound Pt100 sensor, 150mm diameter globe-thermometer equipped with handle. Cable length 2 m.

TP876.I: Wire wound Pt100 sensor, 50mm diameter globe-thermometer equipped with handle. Cable length 2 m.

TP87.O: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 70 mm. Cable length 2 m.

TP878.O: Thin film Pt100 sensor, contact probe for solar panels. Cable length 2 m.

TP878.1.O: Thin film Pt100 sensor, contact probe for solar panels. Cable length 5 m.

TP879.O: Thin film Pt100 sensor, penetration probe for compost. Stem Ø 8 mm, length 1000 mm. Cable length 2 m.

Temperature probes without SICRAM module

TP47.100.O: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP47: Connector for Pt100 4-wire probes without SICRAM module.

TP87.100.O: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 70mm. 4-wires connection cable with connector, length 1 m.



HD 98569 MULTIPARAMETER INSTRUMENT: pH - CONDUCTIVITY DISSOLVED OXYGEN – TEMPERATURE

The **HD 98569** is a portable multi-parameter data logger for electrochemical measures: **pH**, **conductivity**, **dissolved oxygen** and **temperature**. It is fitted with a large back-lighted LCD display.

The instrument measures:

- **pH, mV, redox potential (ORP)** with pH, redox or combined pH/temperature electrodes complete with **SICRAM module**;
 - **conductivity, resistivity** in liquids, **total dissolved solids (TDS)**, and **salinity** with combined 4-ring and 2-ring conductivity and temperature probes **with SICRAM module**.
 - **Concentration of dissolved oxygen** in liquids (in mg/l), **saturation index** (in %) using **SICRAM combined probes** of polarographic type with two or three electrodes or galvanic type, with integrated temperature sensor.
- The instrument is fitted with input for the measurement of **temperature** with Pt100 immersion, penetration, contact or air probes with SICRAM module.
- The pH electrode calibration can be carried out up to five points and the calibration sequence can be chosen from a list of 8 buffers. Temperature compensation can be automatic or manual.
 - The conductivity probe calibration can be performed with automatically detected conductivity calibration solutions: 147µS/cm, 1413µS/cm, 12880µS/cm, 111800µS/cm or manually with calibration solutions having different values.
 - The dissolved oxygen probe's quick calibration function guarantees long-term correctness of the performed measurements.
 - pH, conductivity dissolved oxygen and temperature probes fitted with SICRAM module can store factory and calibration data inside.

The HD 98569 is a **data logger**, it stores up to 200 single screens (labels) and up to 9000 samples in continuous storage mode: pH or mV, conductivity or resistivity or TDS or salinity, concentration of dissolved oxygen and saturation index and temperature.

The data can be transferred from the instrument connected to a PC via the multi-standard RS232C serial port and USB 2.0-1.1.

The instruments equipped with **HD22BT** Bluetooth option can transfer the data without any connection to a PC fitted with USB/Bluetooth converter HD USBKL1, or to the printer **HD40.2** with Bluetooth interface or to a PC with Bluetooth input.

The serial connection RS232C can be used for direct printing of labels with a 24 column printer (**HD40.1** or **HD40.2**).

The software **DeltaLog11** (vers. 2.0 and subsequent ones) allows instrument management and configuration, and data processing on PC.

Technical characteristics of HD 98569

Measured values

pH - mV
 χ - Ω - TDS - NaCl
 mg/l O₂ - %O₂
 °C - °F

Instrument

Dimensions

(LengthxWidthxHeight)

250x100x50mm

Weight

640g (complete with batteries)

Materials

ABS, rubber

Display

Graphic, back lighted LCD, 56x38mm.
 128x64 points

Operating conditions

Working temperature

-5 ... 50°C

Storage temperature

-25 ... 65°C

Working relative humidity

0 ... 90% RH without condensate

Protection degree

IP66

Power

Batteries

4 batteries 1.5V type AA

Autonomy (with probes connected)

25 hours with 1800mAh alkaline batteries

Mains (cod. SWD10)

12Vdc/1A (positive at centre)

Security of memorized data

Unlimited

Time

Date and hour

Schedule in real time

Accuracy

1min/month max. drift

Continuous storage (LOG key)

Quantity

9000 samples of the three inputs

Type

organised in 1800 pages containing 5 samples each

Storage interval

1s ... 999s

Storage on command (MEM key)

Quantity

200 samples of the three inputs

Type

organised in 200 pages containing 1 sample each



- ① Only conductivity probes with SICRAM module.
- ② Input for O₂ and temperature probes or for only temperature probes with SICRAM module.
- ③ Input for pH, mV, pH and temperature probes or for only temperature probes with SICRAM module.
- ④ External Power supply.
- ⑤ RS232 or USB interface.

Calibration storage pH and Dissolved Oxygen	Last 8 pH and dissolved oxygen calibrations. The last 2 are saved in the SICRAM memory of the probe as well.
Conductivity	Last calibration is saved in the SICRAM memory of the probe.
RS232C serial interface	
Type	RS232C electrically isolated
Baud rate	Can be set from 1200 to 38400 baud
Data bit	8
Parity	None
Stop bit	1
Flow control	Xon/Xoff
Length of serial cable	Max 15m

USB interface Type	1.1 - 2.0 electrically isolated
Bluetooth interface	optional for PCs fitted with Bluetooth input. The interface can be installed in Delta Ohm only.

Connections Enabled inputs for temperature probes with SICRAM module	pH/mV and O ₂ inputs.
Input for pH/temperature with SICRAM module	8-pole male DIN45326 connector
Input for conductivity/temperature with SICRAM module	8-pole male DIN45326 connector
Input for dissolved oxygen/temperature with SICRAM module	8-pole male DIN45326 connector
RS232C / USB interface	8-pole MiniDin female connector
Bluetooth	Optional
Mains adapter	2-pole (Ø5.5mm- Ø2.1mm). Positive at centre (e.g. SWD10).

■ Measurement of pH by instrument

Measuring range	-9.999...+19.999pH
Resolution	0.01 o 0.001pH selectable from menu
Accuracy	±0.001pH ±1 digit
Input impedance	>10 ¹² Ω
Calibration error @25°C	IOffset1 > 20mV Slope > 63mV/pH or Slope < 50mV/pH Sensitivity > 106.5% or Sensitivity < 85%
Calibration points	Up to 5 points from a list of 8 automatically detected buffers
Temperature compensation	-50...150°C
Automatically detected standard solutions @25°C	1.679pH - 4.000pH - 4.010pH 6.860pH - 7.000pH - 7.648pH 9.180pH - 10.010pH

Measurement of mV by instrument

Measuring range	-1999.9...+1999.9mV
Resolution	0.1mV
Accuracy	±0.1mV ±1 digit
Drift after 1 year	0.5mV/year

■ Measurement of conductivity by instrument

Measurement range (K cell=0.01)	0.000...1.999µS/cm	Resolution 0.001µS/cm
Measurement range (K cell=0.1)	0.00...19.99µS/cm	0.01µS/cm
Measurement range (K cell=1)	0.0...199.9µS/cm	0.1µS/cm
	200...1999µS/cm	1µS/cm
	2.00...19.99mS/cm	0.01mS/cm
	20.0...199.9mS/cm	0.1mS/cm
Measurement range (K cell=10)	200...1999mS/cm	1mS/cm
Accuracy (conductivity) instrument	±0.5% ±1 digit	

Measurement of resistivity by instrument		Resolution (*)
Measurement range (K cell=0.01)	Up to 1GΩcm	(*)
Measurement range (K cell=0.1)	Up to 100MΩ·cm	0.1Ωcm
Measurement range (K cell=1)	5.0...199.9Ω·cm	1Ω·cm
	200...999Ω·cm	0.01kΩcm
	1.00k...19.99kΩ·cm	0.1kΩcm
	20.0k...99.9kΩ·cm	1kΩcm
	100k...999kΩ·cm	1MΩ·cm
	1...10MΩ·cm	0.1Ωcm
Measurement range (K cell=10)	0.5...5.0Ω·cm	
Accuracy (resistivity) instrument	±0.5% ±1 digit	

(*) The resistivity measurement is obtained from the reciprocal of conductivity measurement. Close to the full scale, the indication of resistivity appears like reported in the table below:

K cell = 0.01 cm ⁻¹		K cell = 0.1 cm ⁻¹	
Conductivity (µS/cm)	Resistivity (MΩ·cm)	Conductivity (µS/cm)	Resistivity (MΩ·cm)
0.001 µS/cm	1000 MΩ·cm	0.01 µS/cm	100 MΩ·cm
0.002 µS/cm	500 MΩ·cm	0.02 µS/cm	50 MΩ·cm
0.003 µS/cm	333 MΩ·cm	0.03 µS/cm	33 MΩ·cm
0.004 µS/cm	250 MΩ·cm	0.04 µS/cm	25 MΩ·cm

Measurement of total dissolved solids (with coefficient χ /TDS=0.5)

Measurement range (K cell=0.01)	0.00...1.999mg/l	Resolution 0.005mg/l
Measurement range (K cell=0.1)	0.00...19.99mg/l	0.05mg/l
Measurement range (K cell=1)	0.0...199.9 mg/l	0.5 mg/l
	200...1999 mg/l	1 mg/l
	2.00...19.99 g/l	0.01 g/l
	20.0...199.9 g/l	0.1 g/l
	100...999 g/l	1 g/l
Measurement range (K cell=10)		
Accuracy (total dissolved solids) instrument	±0.5% ±1 digit	

Measurement of salinity

Measurement range	0.000...1.999g/l	Resolution 1mg/l
	2.00...19.99g/l	10mg/l
	20.0...199.9 g/l	0.1 g/l
Accuracy (salinity) instrument	±0.5% ±1 digit	

Automatic/manual temperature compensation

	0...100°C with $\alpha_T = 0.00...4.00\%/^{\circ}\text{C}$
--	--

Reference temperature

	0...50°C (Default values 20°C or 25°C)
--	--

Conversion factor χ / TDS

	0.4...0.8
--	-----------

Admitted cell constants K (cm⁻¹)

	0.01...20.00
--	--------------

Automatically detected standard solutions (@25°C)

	147µS/cm
	1413µS/cm
	12880µS/cm
	111800µS/cm

■ Measurement of concentration of dissolved oxygen

Measurement range	0.00...90.00mg/l
Resolution	0.01mg/l
Accuracy instrument	±0.03mg/l ±1 digit (0...90%, 1013mbar, 20...25°C)

Measurement of saturation index of dissolved oxygen

Measurement range	0.0...600.0%
Resolution	0.1%
Accuracy instrument	±0.3% ±1 digit (in the range 0.0...199.9%) ±1% ±1 digit (in the range 200.0...600.0%)

Salinity setting

Setting	directly from menu or automatically by conductivity measurement
Setting range	0.0...70.0g/l
Resolution	0.1g/l

Temperature measurement with the sensor inside the O₂ probe

Measurement range	0.0...50.0°C
Resolution	0.1°C
Accuracy instrument	±0.1°C ±1 digit
Drift after 1 year	0.1°C/year
Automatic temperature compensation	0...50°C

■ Measurement of temperature by instrument

Pt100 Measurement range	-50...+150°C
Resolution	0.1°C
Accuracy instrument	±0.1°C ±1 digit
Drift after 1 year	0.1°C/year

TECHNICAL DATA OF PROBES AND MODULES EQUIPPED WITH INSTRUMENT Temperature probes Pt100 sensor with SICRAM module

Model	Type	Application field	Accuracy
TP472I	Immersion	-196°C...+500°C	±0.25°C (-196°C...+300°C) ±0.5°C (+300°C...+500°C)
TP472I.0 1/3 DIN Thin Film	Immersion	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP473PI	Penetration	-50°C...+400°C	±0.25°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP473P.0 1/3 DIN Thin Film	Penetration	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP474C.I	Contact	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP474C.0 1/3 DIN Thin Film	Contact	-50°C...+300°C	±0.3°C (-50°C...+300°C)
TP475A.0 1/3 DIN Thin Film	Air	-50°C...+250°C	±0.3°C (-50°C...+250°C)
TP472I.5	Penetration	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP472I.10	Penetration	-50°C...+400°C	±0.30°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP49A.0 Class A Thin Film	Immersion	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AC.0 Class A Thin Film	Contact	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AP.0 Class A Thin Film	Penetration	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP875.I	Globe-thermometer Ø150mm	-30°C...+120°C	±0.25°C
TP876.I	Globe-thermometer Ø50mm	-30°C...+120°C	±0.25°C
TP87.0 1/3 DIN Thin Film	Immersion	-50°C...+200°C	±0.25°C
TP878.0 1/3 DIN Thin Film TP878.1.0 1/3 DIN Thin Film	Photovoltaic	+4°C...+85°C	±0.25°C
TP879.0 1/3 DIN Thin Film	Compost	-20°C...+120°C	±0.25°C

Common characteristics

Temperature drift @ 20°C	0.003%/°C
--------------------------	-----------

24 column printing example

HD 98569
pH / chi / Oxy / temperature
Ser num=12345678

2007 - 01 - 31 12:00:00

LAB POSITION #1

Operator = Amministratore

SAMPLE ID = 00000001

pH EL sernum = 01234567
pH = 7.010
pH out of calibration !

O₂ EL sernum = 76543210
mg/l O₂ = 5.59

chi EL sernum = 98756410
mS = 2.177

Temp = 25.0°C ATC

ORDERING CODES

HD 98569: The kit is composed of: instrument **data logger** HD 98569 for measurement of pH - redox - conductivity - resistivity - TDS - salinity - concentration of dissolved oxygen-saturation index - temperature, 4 1.5V batteries type AA, instructions manual, software DeltaLog11 (vers. 2.0 and subsequent ones), carrying case and SICRAM module pH471.1 (cable 1 meter).

The pH/mV electrodes, conductivity probes, dissolved oxygen probes, temperature probes, standard reference solutions for different measurement types, connection cables for data download to PC or printer have to be ordered separately.

HD2110CSNM: 8-pole connection cable Mini Din - Sub D 9-pole female for RS232C, for connection to PC with RS232C USB input.

HD2101/USB: Connection cable USB 2.0 connector type A - 8-pole Mini Din for connection to PC with USB input.

DeltaLog11: Further unit of software (vers. 2.0 and subsequent ones) for data download and management on PC using Windows operating systems.

SWD10: Stabilized power supply at 100-240Vac/12Vdc-1A mains voltage.

HD40.1: 24-column portable thermal printer, serial interface, 57mm paper width, four NiMH 1.2V rechargeable batteries, SWD10 power supply, instruction manual, 5 thermal paper rolls It uses HD2110 CSNM cable (optional).

HD40.2: 24-column portable thermal printer, **Bluetooth and serial interface**, 57mm paper width, four NiMH 1.2V rechargeable batteries, SWD10 power supply, instruction manual, 5 thermal paper rolls. Requires the module HD22BT (optional) or the cable HD 2110 CSNM (optional).

RCT: The kit includes 4 thermal paper rolls 57mm wide and 32mm in diameter.

BAT-40: Spare battery pack for HD40.1 printer with built-in temperature sensor.

HD22.2: Laboratory electrode holder composed of basis plate with incorporated magnetic stirrer, staff and replaceable electrode holder. Height max. 380mm. For Ø12mm electrodes. Powered by bench top meters of series HD22... with cable HD22.2.1 (optional) or supplier SWD10 (optional)

HD22.3: Laboratory electrode holder with metal basis plate. Flexible electrode holder for free positioning. For Ø 12mm probes.


HD22BT: Bluetooth module for wireless data transmission from instrument PC. **The fitting of the module into the instrument is made exclusively by Delta Ohm, at the time of placing the order.**

SICRAM Modules with S7 input for pH electrodes

pH 471.1: SICRAM module for pH electrodes with S7 standard connection, cable L=1m.


pH 471.2: SICRAM module for pH electrodes with S7 standard connection, cable L=2m.

pH 471.5: SICRAM module for pH electrodes with S7 standard connection, cable L=5m.




pH Electrodes to be connected to pH471... SICRAM module


KP20: Combined pH electrode for general use, GEL-filled, with screw connector S7, body in Epoxy,




KP 50: Combined pH electrode for heavy pollutants, varnishes, emulsions, gel-filled, with S7 screw connector, body in glass.




KP 61: Combined pH electrode, 3 diaphragms for milk, cream, etc. reference filling solution KCl 3M, with screw connector S7, body in glass.




KP 62: Combined pH electrode, 1 diaphragm for pure water, paints, etc. GEL-filled, with screw connector S7, body in glass




KP 64: Combined pH electrode for water, varnishes, emulsions, etc. reference filling solution KCl 3M, with S7 screw connector, body in glass.




KP 70: Combined pH electrode, micro diam. 6.5mm, open junction, GEL-filled, for paste, bread, cheese, etc, with S7 connector, body in epoxy, glass tip.



KP 80: Combined pointed pH electrode, gel-filled, for cream, milk and viscous substances, with screw connector S7, body in glass.




KP100: Flat membrane gel combined pH electrode with S7 screw connector, glass body, for skin, leather, paper.




ORP electrodes to be connected to pH471... SICRAM module

KP90: REDOX PLATINUM electrode, with screw connector S7, reference filling solution KCl 3M, body in glass.



pH electrodes with SICRAM module

KP63TS: Combined pH/temperature electrode, Pt100 sensor, reference filling solution KCl 3M, with SICRAM module, body in glass, Ag/AgCl sat KCl, single diaphragm, for general purpose, 1 m cable length.



pH buffer solutions

HD8642: Buffer solution 4.01pH - 200cc.

HD8672: Buffer solution 6.86pH - 200cc.

HD8692: Buffer solution 9.18pH - 200cc.

Redox buffer solutions

HDR220: Redox buffer solution 220mV 500cc.

HDR468: Redox buffer solution 468mV 500cc.

Electrolyte solutions

KCL 3M: 100cc ready for use solution for refilling of electrodes.

Cleaning and maintenance

HD62PT: Diaphragm cleaning (tiourea in HCl) - 500cc.

HD62PP: Protein cleaning (pepsin in HCl) - 500cc.

HD62RF: Regeneration (fluorhydric acid) - 100cc.


HD62SC: Solution for electrode preservation - 500cc.

Combined conductivity and temperature probes with SICRAM module


SP06TS: Combined conductivity and temperature 4-electrode cell, body in Pocan.

Cell constant K=0.7.


Measurement range 5µS/cm ... 200mS/cm, 0...90°C, max pressure 5 bar.




SPT401.001S: Combined conductivity and temperature 2-electrode cell in stainless steel AISI 316. Cell constant $K=0.01$. Cable 2m.
Measurement range $0.04\mu\text{S}/\text{cm} \dots 20\mu\text{S}/\text{cm}$, $0 \dots 120^\circ\text{C}$. Measurement in closed-cell, max pressure 5 bar.




SPT016GS: Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant $K=0.1$.
Measurement range $0.1\mu\text{S}/\text{cm} \dots 500\mu\text{S}/\text{cm}$, $0 \dots 80^\circ\text{C}$, max pressure 5 bar.



SPT16GS: Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant $K=1$.
Measurement range $10\mu\text{S}/\text{cm} \dots 10\text{mS}/\text{cm}$, $0 \dots 80^\circ\text{C}$, max pressure 5 bar.



SPT10GS: Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant $K=10$.
Measurement range $500\mu\text{S}/\text{cm} \dots 200\text{mS}/\text{cm}$, $0 \dots 80^\circ\text{C}$, max pressure 5 bar.




Standard calibration solutions


HD8747: Standard calibration solution $0.001\text{mol}/\text{l}$ equal to $147\mu\text{S}/\text{cm}$ @ 25°C - 200cc.
HD8714: Standard calibration solution $0.01\text{mol}/\text{l}$ equal to $1413\mu\text{S}/\text{cm}$ @ 25°C - 200cc.
HD8712: Standard calibration solution $0.1\text{mol}/\text{l}$ equal to $12880\mu\text{S}/\text{cm}$ @ 25°C - 200cc.
HD87111: Standard calibration solution $1\text{mol}/\text{l}$ equal to $111800\mu\text{S}/\text{cm}$ @ 25°C - 200cc.

Combined dissolved Oxygen/temperature probes


DO 9709 SS Polarographic combined oxygen and temperature probe with possibility of membrane replacement. $\varnothing 12\text{mm} \times 120\text{mm}$. 2m cable. The code includes: probe, 2 membranes, electrolyte solution and zero point solution.




DO 9709 SS.5 Polarographic combined oxygen and temperature probe with possibility of membrane replacement. $\varnothing 12\text{mm} \times 120\text{mm}$. 5m cable. The code includes: probe, 2 membranes, electrolyte solution and zero point solution.



DO 9709 SS.1 Galvanic combined galvanic oxygen and temperature probe with possibility of membrane replacement. $\varnothing 12\text{mm} \times 76\text{mm}$. $\varnothing 16\text{mm}$ tip with membrane. 2m cable. The code includes: probe, 2 membranes, electrolyte solution and zero point solution.



DO 9709 SS.5.1 Galvanic combined galvanic oxygen and temperature probe with possibility of membrane replacement. $\varnothing 12\text{mm} \times 76\text{mm}$. $\varnothing 16\text{mm}$ tip with membrane. 5m cable. The code includes: probe, 2 membranes, electrolyte solution and zero point solution.



Accessories

DO 9709/20: Calibrator for polarographic probes DO 9709SS and DO 9709SS.5
DO 9709/21: Calibrator for galvanic probes DO 9709SS.1 and DO 9709SS.5.1
DO 9709 SSK: Kit of accessories for polarographic probes DO 9709SS and DO 9709SS.5: 3 membranes, zero point solution and electrolyte.
DO 9709/21K: Kit of accessories for galvanic probes DO 9709SS.1 and DO 9709SS.5.1: 3 membranes, zero point solution and electrolyte.
DO9700: zero oxygen solution.
DO9701: electrolyte solution for polarographic probes DO9709 SS and DO9709 SS.5.
DO9701.1: electrolyte solution for galvanic probes DO9709 SS.1 and DO9709 SS.5.1.

Temperature probes equipped with SICRAM module

TP472I: Wire wound Pt100 sensor, immersion probe. Stem $\varnothing 3\text{mm}$, length 300 mm. Cable length 2 m.
TP472I.0: Thin film Pt100 sensor, immersion probe. Stem $\varnothing 3\text{mm}$, length 230 mm. Cable length 2 m.
TP473P.I: Wire wound Pt100 sensor, penetration probe. Stem $\varnothing 4\text{mm}$, length 150 mm. Cable length 2 m.
TP473P.0: Thin film Pt100 sensor, penetration probe. Stem $\varnothing 4\text{mm}$, length 150 mm. Cable length 2 m.
TP474C.I: Wire wound Pt100 sensor, contact probe. Stem $\varnothing 4\text{mm}$, length 230mm, contact surface $\varnothing 5\text{mm}$. Cable length 2 m.
TP474C.0: Thin film Pt100 sensor, contact probe. Stem $\varnothing 4\text{mm}$, length 230mm, contact surface $\varnothing 5\text{mm}$. Cable length 2 m.
TP475A.0: Thin film Pt100 sensor, air probe. Stem $\varnothing 4\text{mm}$, length 230mm. Cable length 2 m.
TP472I.5: Thin film Pt100 sensor, penetration probe. Stem $\varnothing 6\text{mm}$, length 500 mm. Cable length 2 m.
TP472I.10: Thin film Pt100 sensor, penetration probe. Stem $\varnothing 6\text{mm}$, length 1000mm. Cable length 2 m.
TP49A.0: Thin film Pt100 sensor, immersion probe. Stem $\varnothing 2,7\text{mm}$, length 150mm. Cable length 2 m. Aluminium handle
TP49AC.0: Thin film Pt100 sensor, contact probe. Stem $\varnothing 4\text{mm}$, length 150mm. Cable length 2 m. Aluminium handle
TP49AP.0: Thin film Pt100 sensor, penetration probe. Stem $\varnothing 2,7\text{mm}$, length 150mm. Cable length 2 m. Aluminium handle
TP875.I: Wire wound Pt100 sensor, 150mm diameter globe-thermometer equipped with handle. Cable length 2 m.
TP876.I: Wire wound Pt100 sensor, 50mm diameter globe-thermometer equipped with handle. Cable length 2 m.
TP87.0: Thin film Pt100 sensor, immersion probe. Stem $\varnothing 3\text{mm}$, length 70 mm. Cable length 2 m.
TP878.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 2 m.
TP878.1.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 5 m.
TP879.0: Thin film Pt100 sensor, penetration probe for compost. Stem $\varnothing 8\text{mm}$, length 1000 mm. Cable length 2 m.



HD 3405.2 BENCH-TOP pH METER

The **HD3405.2** is a bench top instrument for electrochemical measures: **pH, and temperature**. The displayed data can be stored (**datalogger**) and can be transferred to PC or serial printer thanks to the multi-standard serial port RS232C and USB2.0 and software DeltaLog9 (Vers.2.0 and subsequent ones). The storing and printing parameters can be set from menu.

The **HD3405.2** measures **pH, redox potential (ORP)** in mV. It measures **temperature** with Pt100 or Pt1000 immersion, penetration or contact probes. The pH electrode calibration can be carried out on one, two or three points and the calibration sequence can be chosen from a list of 13 buffers. The display shows continually the temperature in °C or °F and one of the parameters according to the connected probe type. Printing and storage always include the temperature in °C or °F and one selectable parameter for each probe type.

Other functions include: Max, Min and Avg function, the Auto-HOLD function, the automatic turning off which can also be excluded.

The instruments have IP66 protection degree.

Technical characteristics HD3405.2

pH - mV - °C/°F measurement

Instrument

Dimensions (Length x Width x Height)	220x120x55mm
Weight	460g (complete with batteries)
Materials	ABS, rubber
Display	2x4½ characters plus symbols visible area: 52x42mm

Operating conditions

Working temperature	-5 ... 50°C
Storage temperature	-25 ... 65°C
Working relative humidity	0 ... 90% RH without condensation
Protection degree	IP66

Power

Batteries	3 batteries 1.5V type AA
Autonomy (only batteries)	100 hours with 1800mAh alkaline batteries
Mains (cod. SWD10)	Output mains adapter 100-240Vac/ 12Vdc-1A

Security of memorized data

Unlimited

Storage of measured values

Type	2000 pages of 17 samples each
Quantity	34,000 sets of measures made up of [pH or mV] and [°C or °F].

Selectable storage interval

1s, 5s, 10s, 15s, 30s, 1min, 2min, 5min, 10min, 15min, 20min, 30min and 1hour

Time

Date and hour	Schedule in real time
Accuracy	1min/month max departure

Serial interface RS232C

Type	RS232C electrically isolated
Baud rate	Can be set from 1200 to 38400 baud
Data bit	8
Parity	None
Stop bit	1
Flow Control	Xon/Xoff
Serial cable length	Max 15m
Selectable print interval	immediate or 1s, 5s, 10s, 15s, 30s, 1min, 2min, 5min, 10min, 15min, 20min, 30min and 1hour

USB Interface



Type	1.1 - 2.0 electrically isolated
------	---------------------------------

Connections

Serial interface and USB	8-pole MiniDin connector
Mains adapter (cod. SWD10)	2-pole connector (positive at centre) 12Vdc/1A

Measurement connections

Temperature probe input with SICRAM module or TP47 module	8-pole male DIN45326 connector
pH/mV input	female BNC



pH Measurement

Measurement range	-2.000...+19.999pH
Resolution	0.01 or 0.001pH selectable from menu
Accuracy	±0.001pH ±1 digit
Input impedance	>10 ¹² Ω
Calibration error @25°C	IOffset1 > 20mV Slope < 50mV/pH or Slope > 63mV/pH Sensitivity < 85% or Sensitivity > 106.5%
Automatic / manual temperature compensation	-50...+150°C

mV Measurement

Measurement range	-1999.9...+1999.9mV
Resolution	0.1mV
Accuracy	±0.1mV ±1 digit
Drift after 1 year	0.5mV/year

Automatically detected pH standard solutions (@25°C)

1.679pH - 2.000pH - 4.000pH - 4.008pH
4.010pH - 6.860pH - 6.865pH - 7.000pH
7.413pH - 7.648pH - 9.180pH - 9.210pH
10.010pH

Temperature Measurement

Pt100 measurement range	-200...+650°C
Pt1000 measurement range	-200...+650°C
Resolution	0.1°C
Accuracy	±0.1°C ±1 digit
Drift after 1 year	0.1°C/year

ORDERING CODES

HD3405.2: The kit is composed of: instrument HD3405.2 **datalogger**, for measurement of pH - redox - temperature, 3 1.5V alkaline batteries, operating manual and **DeltaLog9 version 2.0**.

pH/mV electrodes, temperature probes, standard reference solutions connection cables for pH electrodes with S7 connector, cables for data download to PC or printer have to be ordered separately.

Accessories

HD2110CSNM: 8-pole connection cable Mini Din - Sub D 9-pole female for RS232C, for connection to PC without USB input.

HD2101/USB: Connection cable USB 2.0 connector type A - 8-pole Mini Din for connection to PC with USB input.

SWD10: Stabilized power supply at 100-240Vac/12 Vdc-1A mains voltage.

HD40.1: Portable, serial input, 24 column thermal printer, 57mm paper width.

HD22.2: Laboratory electrode holder composed of basis plate with incorporated magnetic stirrer, staff and replaceable electrode holder. Height max. 380mm. Powered by bench-top meters of the series HD22... with cable HD22.2.1 (**optional**) or supplier SWD10 (**optional**).

HD22.3: Laboratory electrode holder with metal basis plate. Flexible electrode holder for free positioning. For Ø 12mm probes.

TP47: Module for the connection of Pt100 4-wire and Pt1000 2-wire probes to instrument series HD34..., without amplifying electronics.

pH Electrodes

KP 20: Gel pH combined electrode for general use, with S7 screw connector, EPOXY body.

KP 30: Gel pH combined electrode for general use, 1m cable with BNC, EPOXY body.

KP 50: Gel pH combined electrode, porous Teflon ring junction, suitable for emulsions, demineralised water and waste water with S7 screw connector, glass body.

KP 61: 3 diaphragm liquid filled pH combined electrode for wine, milk, cream, etc., S7 screw connector, liquid reference filling, glass body.

KP 62: 1 diaphragm gel pH combined electrode for general use, pure water, varnishes, gel filled, S7 screw connector, glass body.

KP 63: 1 liquid filled pH combined electrode for general use, varnishes, 1m cable with BNC, glass body.

KP 64: Liquid filled pH combined electrode, Teflon ring diaphragm, for wine, varnishes, emulsions, S7 screw connector, glass body.

KP 70: Pointed gel combined pH microelectrode diam. 6 x L=70 mm., with S7 screw connector, EPOXY body, glass tip, open junction for meat and cheese.

KP 80: Pointed gel pH combined electrode, with S7 screw connector, glass body, for cream, milk, viscous material, open junction.

KP100: Flat membrane gel combined pH electrode with S7 screw connector, glass body, for skin, leather, paper.

Characteristics and dimensions of the probes on page WA-76.

CP: 1.5m extension cable with BNC/S7 connector for electrode without cable, thread S7.

CP 5: 5m extension cable with BNC/S7 connector for electrode without cable, thread S7.

CP 10: 10m extension cable with BNC/S7 connector for electrode without cable, thread S7.

CP 15: 15m extension cable with BNC/S7 connector for electrode without cable, thread S7.

CE : S7 screw connector for pH electrode.

BNC: female BNC for extension cable

ORP Electrodes

KP90: Redox Platinum electrode, with screw connector S7, electrolyte KCl 3M, body in glass.

KP91: Redox Platinum electrode with 1m cable, GEL filled, body in glass.

Electrode dimensions and characteristics at page WA-76

pH buffer solutions

HD8642: Buffer solution 4.01pH - 200cc.

HD8672: Buffer solution 6.86pH - 200cc.

HD8692: Buffer solution 9.18pH - 200cc.

Redox buffer solutions

HDR220: Redox buffer solution 220mV 0,5 l.

HDR468: Redox buffer solution 468mV 0,5 l.

Electrolyte solutions

KCL 3M: 100cc ready for use solution for refilling of the electrodes.

Cleaning and maintenance

HD62PT: Diaphragm cleaning (tiourea in HCl) - 500ml.

HD62PP: Protein cleaning (pepsin in HCl) - 500ml.

HD62RF: Regeneration (fluorhydric acid) - 100ml.

HD62SC: Solution for electrode preservation - 500ml.

Temperature probes equipped with SICRAM module

TP472I: Wire wound Pt100 sensor, immersion probe. Stem Ø 3 mm, length 300 mm. Cable length 2 m.

TP472I.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 230 mm. Cable length 2 m.

TP473P.I: Wire wound Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP473P.0: Thin film Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP474C.I: Wire wound Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP474C.0: Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP475A.0: Thin film Pt100 sensor, air probe. Stem Ø 4mm, length 230mm. Cable length 2 m.

TP472I.5: Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 500 mm. Cable length 2 m.

TP472I.10: Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 1000mm. Cable length 2 m.

TP49A.0: Thin film Pt100 sensor, immersion probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AC.0: Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AP.0: Thin film Pt100 sensor, penetration probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP875.I: Wire wound Pt100 sensor, 150mm diameter globe-thermometer equipped with handle. Cable length 2 m.

TP876.I: Wire wound Pt100 sensor, 50mm diameter globe-thermometer equipped with handle. Cable length 2 m.

TP87.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 70 mm. Cable length 2 m.

TP878.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 2 m.

TP878.1.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 5 m.

TP879.0: Thin film Pt100 sensor, penetration probe for compost. Stem Ø 8 mm, length 1000 mm. Cable length 2 m.

Temperature probes without SICRAM module

TP47.100.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP47.1000.0: Thin film Pt1000 sensor, immersion probe. Probe's Stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP47: Connector for Pt100 4-wire and Pt1000 2-wire probes without SICRAM module.

TP87.100.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 70mm. 4-wires connection cable with connector, length 1 m.

TP87.1000.0: Thin film Pt1000 sensor, immersion probe. Stem Ø 3mm, length 70mm. 2-wires connection cable with connector, length 1 m.



pH



mV



HD 2205.2 BENCH-TOP pH METER

The **HD2205.2** is a bench top instrument for electrochemical measures: **pH** and **temperature**. It is fitted with a large backlit LCD display.

The **HD2205.2** is equipped with two BNC inputs for the measurement of **pH**, **mV**, **redox potential** (ORP) with pH or redox electrodes, or electrodes with separate reference, and one input for combined pH/temperature probes fitted with SICRAM module.

All models are fitted with input for the measurement of **temperature** with Pt100 or Pt1000 immersion, penetration or contact probes. The temperature probes are equipped with an automatic recognition module and factory calibration data are stored inside.

The pH electrode calibration can be carried out up to five points and the calibration sequence can be chosen from a list of 13 buffers. Temperature compensation can be automatic or manual.

The HD 2205.2 is a **datalogger**, it stores up to 2000 samples of data: pH, mV and temperature. The data can be transferred from the instrument connected to a PC via the RS232C and USB 2.0 serial ports. The storing parameters can be configured using the menu.

The RS232C serial port can be used to transfer the acquired measurements to a 24 column portable printer in real time (HD40.1 or HD40.2).

The instruments equipped with **HD22BT** (Bluetooth) option can transfer data, without any connection, to a PC or printer fitted with Bluetooth input or through Bluetooth/RS232C converter.

The software DeltaLog11 allows instrument management and configuration, and data processing on PC.

The instruments have IP66 protection degree.

Technical characteristics HD2205.2

pH - mV - °C - °F measurement

Instrument

Dimensions (Length x Width x Height)	265x185x70mm
Weight	490g
Materials	ABS, rubber
Display	Back lighted, matrix point display. 240x64 points, visible area: 128x35mm

Operating conditions

Working temperature	-5 ... 50°C
Storing temperature	-25 ... 65°C
Working relative humidity	0 ... 90% R.H. without condensate

Protection degree

IP66

Power

Mains adapter (cod. SWD10) 12Vdc/1A

Auxiliary socket

For supplying of electrode holder with built-in stirrer HD22.2

Security of memorized data

Unlimited

Time

Date and hour

Real time schedule with backup battery 3.6V - ½AA

Accuracy

1 min/month max drift

Measured values storing

Quantity	2000 screens
Storage interval	1s ... 999s


Calibration storage

Quantity	Last 8 calibrations of each physical quantity
----------	---

RS232C serial interface

Type

RS232C electrically isolated



Baud rate	Can be set from 1200 to 115200 baud
Data bit	8
Parity	None
Stop bit	1
Flow Control	Xon/Xoff
Length of serial cable	Max 15m

USB Interface

Type	1.1 - 2.0 electrically isolated
Bluetooth Interface	HD22BT optional

Connections

Input for temperature probes with SICRAM module⑤	8-pole male DIN45326 connector
Inputs pH/mV ① - ②	female BNC
Inputs for SICRAM module pH/temperature ③ - ④	8-pole male DIN45326 connector
Serial interface	DB9 connector (9-pole male)
USB interface	USB connector type B
Bluetooth	Optional
Mains adapter	2-pole connector (Ø5.5mm-2.1mm). Positive at centre.
Socket for power supply of electrode holder with built-in magnetic stirrer	2-pole connector (Ø5.5mm-2.1mm). Positive at centre (output 12Vdc/200mA max).

Measurement of pH by instrument

Measuring range	-9.999...+19.999pH
Resolution	0.01 o 0.001pH selectable from menu
Accuracy	±0.001pH ±1digit
Input impedance	>10 ¹² Ω
Calibration error @25°C	IOffsetl > 20mV Slope > 63mV/pH o Slope < 50mV/pH Sensitivity > 106.5% or Sensitivity < 85%
Calibration points	Up to 5 points from a list of 13 automatically detected buffers.
Automatically detected pH standard solutions (@25°C)	1.679pH - 2.000pH - 4.000pH - 4.008pH - 4.010pH 6.860pH - 6.865pH - 7.000pH - 7.413pH - 7.648pH 9.180pH - 9.210pH - 10.010pH

Measurement of mV by instrument

Measuring range	-1999.9...+1999.9mV
Resolution	0.1mV
Accuracy	±0.1mV ±1digit
Drift after 1 year	0.5mV/year

Measurement of temperature by instrument

Pt100 Measuring range	-50...+150°C
Pt1000 Measuring range	-50...+150°C
Resolution	0.1°C
Accuracy	±0.1°C ±1digit
Drift after 1 year	0.1°C/year

TECHNICAL DATA OF PROBES AND MODULES EQUIPPED WITH INSTRUMENT Temperature probes Pt100 sensor with SICRAM module

Model	Type	Application field	Accuracy
TP472I	Immersion	-196°C...+500°C	±0.25°C (-196°C...+300°C) ±0.5°C (+300°C...+500°C)
TP472I.0 1/3 DIN Thin Film	Immersion	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP473P.I	Penetration	-50°C...+400°C	±0.25°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP473P.0 1/3 DIN Thin Film	Penetration	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP474C.I	Contact	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP474C.0 1/3 DIN Thin Film	Contact	-50°C...+300°C	±0.3°C (-50°C...+300°C)
TP475A.0 1/3 DIN Thin Film	Air	-50°C...+250°C	±0.3°C (-50°C...+250°C)
TP472I.5	Penetration	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP472I.10	Penetration	-50°C...+400°C	±0.30°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP49A.0 Class A Thin Film	Immersion	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AC.0 Class A Thin Film	Contact	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AP.0 Class A Thin Film	Penetration	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP875.I	Globe-thermometer Ø150mm	-30°C...+120°C	±0.25°C
TP876.I	Globe-thermometer Ø50mm	-30°C...+120°C	±0.25°C
TP87.0 1/3 DIN Thin Film	Immersion	-50°C...+200°C	±0.25°C
TP878.0 1/3 DIN Thin Film TP878.1.0 1/3 DIN Thin Film	Photovoltaic	+4°C...+85°C	±0.25°C
TP879.0 1/3 DIN Thin Film	Compost	-20°C...+120°C	±0.25°C

Common characteristics

Temperature drift @ 20°C	0.003%/°C
--------------------------	-----------

4 wires Pt100 and 2 wires Pt1000 Probes

Model	Type	Application field	Accuracy
TP47.100.0 1/3 DIN Thin Film	4 wires Pt100	-50...+250°C	1/3 DIN
TP47.1000.0 1/3 DIN Thin Film	2 wires Pt1000	-50...+250°C	1/3 DIN
TP87.100.0 1/3 DIN Thin Film	4 wires Pt100	-50...+200°C	1/3 DIN
TP87.1000.0 1/3 DIN Thin Film	2 wires Pt1000	-50...+200°C	1/3 DIN

Common features

Temperature drift @20°C	
Pt100	0.003%/°C
Pt1000	0.005%/°C



pH



mV

Ordering codes

HD2205.2: The kit is composed of: instrument HD2205.2 for measurement of pH - redox - temperature, **data logger**, SWD10 stabilized power supply at mains voltage 100-240Vac/12Vdc-1A, instructions manual and software DeltaLog11.

pH/mV electrodes, temperature probes, standard reference solutions, connection cables for pH electrodes with S7 connector, cables for data download to PC or printer have to be ordered separately.

Accessories

9CPRS232: Connection cable SubD female 9- pole for serial output RS232C.

CP22: USB 2.0 connection cable - connector type A - connector type B.

DeltaLog11: Software for download and management of the data on PC using Windows operating systems.

SWD10: Stabilized power supply at 100-240Vac/12Vdc-1A mains voltage.

HD40.1: 24-column portable thermal printer, serial interface, 57mm paper width, four NiMH 1.2V rechargeable batteries, SWD10 power supply, instruction manual, 5 thermal paper rolls. It uses the cable 9CPRS 232 (optional).

HD40.2: 24-column portable thermal printer, **Bluetooth and serial interface**, 57mm paper width, four NiMH 1.2V rechargeable batteries, SWD10 power supply, instruction manual, 5 thermal paper rolls. Requires the module HD22BT (optional) or the cable 9CPRS 232 (optional).

HD22.2: Laboratory electrode holder composed of basis plate with incorporated magnetic stirrer, staff and replaceable electrode holder. Height max. 380mm. Powered by bench-top meters of the series HD22... with cable HD22.2.1 (optional) or supplier SWD10 (optional).

HD22.3: Laboratory electrode holder with metal basis plate. Flexible electrode holder for free positioning. For Ø 12mm probes.

HD22BT: Bluetooth module for wireless data transmission from instrument to PC. **The fitting of the module into the instrument is made exclusively by Delta Ohm, at the time of placing the order.**

TP47: Connector for Pt100 4-wire and Pt1000 2-wire probes without SICRAM module.

pH electrodes without SICRAM module (Inputs ① and ②)

pH Electrodes

KP 20: Gel pH combined electrode for general use, with S7 screw connector, EPOXY body.

KP 30: Gel pH combined electrode for general use, 1m cable with BNC, EPOXY body.

KP 50: Gel pH combined electrode, porous Teflon ring junction, suitable for emulsions, demineralised water and waste water with S7 screw connector, glass body.

KP 61: 3 diaphragm liquid filled pH combined electrode for wine, milk, cream, etc., S7 screw connector, liquid reference filling, glass body.

KP 62: 1 diaphragm gel pH combined electrode for general use, pure water, varnishes, gel filled, S7 screw connector, glass body.

KP 63: liquid filled pH combined electrode for general use, varnishes, 1m cable with BNC, glass body.

KP 64: Liquid filled pH combined electrode, Teflon ring diaphragm, for wine, varnishes, emulsions, S7 screw connector, glass body.

KP 70: Pointed gel combined pH microelectrode diam. 6 x L=70 mm., with S7 screw connector, EPOXY body, glass tip, open junction for meat and cheese.

KP 80: Pointed gel pH combined electrode, with S7 screw connector, glass body, for cream, milk, viscous material, open junction.

KP100: Flat membrane gel combined pH electrode with S7 screw connector, glass body, for skin, leather, paper.

Electrode dimensions and characteristics at page WA-76

CP: Extension cable 1.5m with BNC connector on one side and S7 on the other side for electrode with S7 connector.

CP5: Extension cable 5m with BNC connector on one side and S7 on the other side for electrode with S7 connector.

CP10: Extension cable 10m with BNC connector on one side and S7 on the other side, for electrode without cable.

CP15: Extension cable 15m with BNC connector on one side and S7 on the other side, for electrode without cable.

CE: S7 screw connector for pH electrode.

BNC: Female BNC for electrode extension.

pH electrodes with SICRAM module (Input ③)


KP63TS: Combined pH/temperature electrode with SICRAM module, body in Epoxy, Ag/AgCl sat KCl, one diaphragm, for general purpose.

SICRAM Module with S7 input for pH electrodes (Input ③)

pH 471.1: SICRAM module for pH electrodes with S7 standard connection, cable L=1m.

pH 471.2: SICRAM module for pH electrodes with S7 standard connection, cable L=2m.

pH 471.5: SICRAM module for pH electrodes with S7 standard connection, cable L=5m.



ORP Electrodes (Inputs ① and ②)

KP90: Redox Platinum electrode, with screw connector S7, electrolyte KCl 3M, body in glass.

KP91: Redox Platinum electrode with 1m cable, GEL filled, body in glass.

Characteristics and dimensions of the probes at page WA-76



HD40.1

pH buffer solutions

HD8642: Buffer solution 4.01pH - 200cc.

HD8672: Buffer solution 6.86pH - 200cc.

HD8692: Buffer solution 9.18pH - 200cc.

Redox buffer solutions

HDR220: Redox buffer solution 220mV 0,5 l.

HDR468: Redox buffer solution 468mV 0,5 l.

Electrolyte solutions

KCL 3M: 100cc ready for use solution for electrode refilling.

Cleaning and maintenance

HD62PT: Diaphragm cleaning (tiourea in HCl) - 500ml.

HD62PP: Protein cleaning (pepsin in HCl) - 500ml.

HD62RF: Regeneration (fluorhydric acid) - 100ml.

HD62SC: Solution for electrode preservation - 500ml.

Temperature probes complete with SICRAM module (Input ⑤)

TP472I: Wire wound Pt100 sensor, immersion probe. Stem Ø 3 mm, length 300 mm. Cable length 2 m.

TP472I.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 230 mm. Cable length 2 m.

TP473PI: Wire wound Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP473P.0: Thin film Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP474C.I: Wire wound Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP474C.0: Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP475A.0: Thin film Pt100 sensor, air probe. Stem Ø 4mm, length 230mm. Cable length 2 m.

TP472I.5: Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 500 mm. Cable length 2 m.

TP472I.10: Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 1000mm. Cable length 2 m.

TP49A.0: Thin film Pt100 sensor, immersion probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AC.0: Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AP.0: Thin film Pt100 sensor, penetration probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP875.I: Wire wound Pt100 sensor, 150mm diameter globe-thermometer equipped with handle. Cable length 2 m.

TP876.I: Wire wound Pt100 sensor, 50mm diameter globe-thermometer equipped with handle. Cable length 2 m.

TP87.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 70 mm. Cable length 2 m.

TP878.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 2 m.

TP878.1.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 5 m.

TP879.0: Thin film Pt100 sensor, penetration probe for compost. Stem Ø 8 mm, length 1000 mm. Cable length 2 m.

Temperature probes complete with TP47 module (input ⑤)

TP47.100.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP47.1000.0: Thin film Pt1000 sensor, immersion probe. Probe's Stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP47: Connector for Pt100 4-wire and Pt1000 2-wire probes without SICRAM module.

TP87.100.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 70mm. 4-wires connection cable with connector, length 1 m.

TP87.1000.0: Thin film Pt1000 sensor, immersion probe. Stem Ø 3mm, length 70mm. 2-wires connection cable with connector, length 1 m.



HD22.2



HD 3456.2 BENCH-TOP pH AND CONDUCTIVITY METER

The **HD3456.2** is a bench top instrument for electrochemical measures: **pH, conductivity and temperature**.

The displayed data can be stored (**datalogger**) and can be transferred to PC or serial printer thanks to the multi-standard serial port RS232C and USB2.0 and software DeltaLog9 (Vers.2.0 and subsequent ones). The storing and printing parameters can be set from menu.

The **HD3456.2** measures **pH, mV, redox potential (ORP), conductivity, resistivity in liquids, total dissolved solids (TDS), and salinity** using combined 4-ring and 2-ring conductivity/temperature probes. **Temperature** is measured by Pt100 or Pt1000 immersion, penetration or contact probes.

The pH electrode calibration, as well as manual, can be carried out automatically on one, two or three points and the calibration sequence can be chosen from a list of 13 buffers.

The conductivity probe calibration can be performed automatically in one or more of the 147µS, 1413µS, 12880µS or 111800µS/cm conductivity calibration solutions.



The display shows continually the temperature in °C or °F and one selectable parameter according to the connected probe type. Printing and storing data always include temperature in °C or °F and one selectable parameter for every kind of probe: i.e. in case of conductivity probe it is possible to select between χ or Ω or TDS or g/l. Other functions of this instrument include: Max, Min and Avg function, the Auto-HOLD function, the automatic turning off which can also be disabled.

The instruments have IP66 protection degree.

Technical characteristics HD3456.2

pH, mV, χ , Ω , TDS, Sal, °C/°F measurement

Instrument

Dimensions (Length x Width x Height)	220x120x55mm
Weight	460g (complete with batteries)
Materials	ABS, rubber
Display	2x4½ characters plus symbols visible area: 52x42mm

Operating conditions

Working temperature	-5 ... 50°C
Storage temperature	-25 ... 65°C
Working relative humidity	0 ... 90% RH without condensation

Protection degree

IP66

Power

Batteries	3 batteries 1.5V type AA
Autonomy (only batteries)	100 hours with 1800mAh alkaline batteries
Mains (cod. SWD10)	Output mains adapter 100-240Vac/ 12Vdc-1A

Security of memorized data

Unlimited

Storage of measured values

Type	2000 pages of 10 samples each
Quantity	20,000 terns of measures made up of [pH or mV], [χ or Ω or TDS or salinity] and temperature.

Selectable storage interval

1s, 5s, 10s, 15s, 30s, 1min, 2min, 5min, 10min, 15min, 20min, 30min and 1hour

Time

Date and hour	Schedule in real time
Accuracy	1min/month max drift

Serial interface RS232C


Type	RS232C electrically isolated
Baud rate	Can be set from 1200 to 38400 baud
Data bit	8
Parity	None
Stop bit	1
Flow Control	Xon/Xoff
Serial cable length	Max 15m
Selectable print interval	immediate or 1s, 5s, 10s, 15s, 30s, 1min, 2min, 5min, 10min, 15min, 20min, 30min and 1ora

USB Interface

Type	1.1 - 2.0 electrically isolated
------	---------------------------------

Connections

Serial interface and USB	8-pole MiniDin connector
Mains adapter (cod. SWD10)	2-pole connector (positive at centre) 12Vdc/1A



Connections
pH/mV input
Conductivity input
Input for temperature probes

Female BNC connector
8-pole male DIN45326 connector
8-pole male DIN45326 connector

Measurement of pH by Instrument

Measurement range -2.000...+19.999pH
Resolution 0.01 or 0.001pH selectable from menu
Accuracy $\pm 0.001\text{pH} \pm 1\text{digit}$
Input impedance $> 10^{12}\Omega$
Calibration error @25°C
Offset $> 20\text{mV}$
Slope $> 63\text{mV/pH}$ or Slope $< 50\text{mV/pH}$
Sensitivity $> 106.5\%$ or Sensitivity $< 85\%$
Automatic / manual temperature compensation -50...+150°C

Measurement of mV by Instrument

Measurement range -1999.9...+1999.9mV
Resolution 0.1mV
Accuracy $\pm 0.1\text{mV} \pm 1\text{digit}$
Drift after 1 year 0.5mV/year

Standard solutions automatically detected (@25°C)

1.679pH - 2.000pH - 4.000pH - 4.008pH
4.010pH - 6.860pH - 6.865pH - 7.000pH
7.413pH - 7.648pH - 9.180pH - 9.210pH
10.010pH

Measurement of conductivity by Instrument

Measurement range (Kcell=0.01)	0.000...1.999 $\mu\text{S/cm}$	Resolution 0.001 $\mu\text{S/cm}$
Measurement range (Kcell=0.1)	0.00...19.99 $\mu\text{S/cm}$	0.01 $\mu\text{S/cm}$
Measurement range (Kcell=1)	0.0...199.9 $\mu\text{S/cm}$	0.1 $\mu\text{S/cm}$
	200...1999 $\mu\text{S/cm}$	1 $\mu\text{S/cm}$
	2.00...19.99mS/cm	0.01mS/cm
	20.0...199.9mS/cm	0.1mS/cm
Measurement range (Kcell=10)	200...1999mS/cm	1mS/cm
Accuracy (conductivity)	$\pm 0.5\% \pm 1\text{digit}$	

Measurement of resistivity by Instrument

Measurement range (Kcell=0.01)	Up to 1G $\Omega\cdot\text{cm}$	Resolution (*)
Measurement range (Kcell=0.1)	Up to 100M $\Omega\cdot\text{cm}$	(*)
Measurement range (Kcell=1)	5.0...199.9M $\Omega\cdot\text{cm}$	0.1M $\Omega\cdot\text{cm}$
	200...999M $\Omega\cdot\text{cm}$	1M $\Omega\cdot\text{cm}$
	1.00k...19.99k $\Omega\cdot\text{cm}$	0.01k $\Omega\cdot\text{cm}$
	20.0k...99.9k $\Omega\cdot\text{cm}$	0.1k $\Omega\cdot\text{cm}$
	100k...999k $\Omega\cdot\text{cm}$	1k $\Omega\cdot\text{cm}$
	1...10M $\Omega\cdot\text{cm}$	1M $\Omega\cdot\text{cm}$
Measurement range (Kcell=10)	0.5...5.0M $\Omega\cdot\text{cm}$	0.1M $\Omega\cdot\text{cm}$
Accuracy (resistivity)	$\pm 0.5\% \pm 1\text{digit}$	

Measurement of total dissolved solids (with coefficient $\chi/\text{TDS}=0.5$)

Measurement range (Kcell=0.01)	0.00...1.999mg/l	0.005mg/l
Measurement range (Kcell=0.1)	0.00...19.99mg/l	0.05mg/l
Measurement range (Kcell=1)	0.0...199.9mg/l	0.5mg/l
	200...1999mg/l	1mg/l
	2.00...19.99g/l	0.01g/l
	20.0...99.9g/l	0.1g/l
Measurement range (Kcell=10)	100...999g/l	1g/l
Accuracy (total dissolved solids)	$\pm 0.5\% \pm 1\text{digit}$	

Measurement of salinity

Measurement range	0.000...1.999g/l	1mg/l
	2.00...19.99g/l	10mg/l
	20.0...199.9g/l	0.1g/l
Accuracy (salinity)	$\pm 0.5\% \pm 1\text{digit}$	

Automatic/manual temperature compensation

0...100°C with α_T that can be selected from 0.00 to 4.00%/°C

Reference temperature

20°C or 25°C selectable from menu

χ/TDS conversion factor

0.4...0.8

Cell constant K (cm⁻¹)

0.01 - 0.1 - 0.7 - 1.0 - 10.0

Standard solutions automatically detected (@25°C)

147 $\mu\text{S/cm}$
1413 $\mu\text{S/cm}$
12880 $\mu\text{S/cm}$
111800 $\mu\text{S/cm}$

Measurement of temperature by Instrument

Pt100 measurement range -50...+200°C
Pt1000 measurement range -50...+200°C
Resolution 0.1°C
Accuracy $\pm 0.25^\circ\text{C}$
Drift after 1 year 0.1°C/year

(*) The resistivity measurement is obtained from the reciprocal of conductivity measurement. Close to the bottom of the scale, the indication of resistivity appears like reported in the table below:

K cell = 0.01 cm ⁻¹		K cell = 0.1 cm ⁻¹	
Conductivity ($\mu\text{S/cm}$)	Resistivity (M $\Omega\cdot\text{cm}$)	Conductivity ($\mu\text{S/cm}$)	Resistivity (M $\Omega\cdot\text{cm}$)
0.001 $\mu\text{S/cm}$	1000 M $\Omega\cdot\text{cm}$	0.01 $\mu\text{S/cm}$	100 M $\Omega\cdot\text{cm}$
0.002 $\mu\text{S/cm}$	500 M $\Omega\cdot\text{cm}$	0.02 $\mu\text{S/cm}$	50 M $\Omega\cdot\text{cm}$
0.003 $\mu\text{S/cm}$	333 M $\Omega\cdot\text{cm}$	0.03 $\mu\text{S/cm}$	33 M $\Omega\cdot\text{cm}$
0.004 $\mu\text{S/cm}$	250 M $\Omega\cdot\text{cm}$	0.04 $\mu\text{S/cm}$	25 M $\Omega\cdot\text{cm}$

ORDERING CODES

HD3456.2: The kit is composed of: instrument HD3456.2 **data logger**, for the measurement of pH - redox - conductivity - resistivity - TDS - salinity - temperature, 3 1.5V alkaline batteries, operating manual and **DeltaLog9 version 2.0**.

pH/mV electrodes, conductivity probes, temperature probes, standard reference solutions for different measurement types, connection cables for pH electrodes with S7 connector, cables for data download to PC or printer have to be ordered separately.

Accessories

HD2110CSNM: 8-pole connection cable Mini Din - Sub D 9-pole female for RS232C, for connection to PC without USB input.

HD2101/USB: Connection cable USB 2.0 connector type A - 8-pole Mini Din for connection to PC with USB input.

SWD10: Stabilized power supply at 230Vac/9Vdc-300mA mains voltage.

HD40.1: Portable, serial input, 24 column thermal printer, 57mm paper width.

HD22.2: Laboratory electrode holder composed of basis plate with incorporated magnetic stirrer, staff and replaceable electrode holder. Height max. 380mm. Powered by bench-top meters of the series HD22... with cable HD22.2.1 (**optional**) or supplier SWD10 (**optional**).

HD22.3: Laboratory electrode holder with metal basis plate. Flexible electrode holder for free positioning. For \varnothing 12mm probes.

TP47: Connector for Pt100 4-wire and Pt1000 2-wire probes without SICRAM module.



pH



χ



Ω



TDS

pH Electrodes

- KP 20:** Gel pH combined electrode for general use, with S7 screw connector, EPOXY body.
- KP 30:** Gel pH combined electrode for general use, 1m cable with BNC, EPOXY body.
- KP 50:** Gel pH combined electrode, porous Teflon ring junction, suitable for emulsions, demineralised water and waste water with S7 screw connector, glass body.
- KP 61:** 3 diaphragm liquid filled pH combined electrode for wine, milk, cream, etc., S7 screw connector, liquid reference filling, glass body.
- KP 62:** 1 diaphragm gel pH combined electrode for general use, pure water, varnishes, gel filled, S7 screw connector, glass body.
- KP 63:** liquid filled pH combined electrode for general use, varnishes, 1m cable with BNC, glass body.
- KP 64:** Liquid filled pH combined electrode, Teflon ring diaphragm, for wine, varnishes, emulsions, S7 screw connector, glass body.
- KP 70:** Pointed gel combined pH microelectrode diam. $6 \times L=70$ mm., with S7 screw connector, EPOXY body, glass tip, open junction for meat and cheese.
- KP 80:** Pointed gel pH combined electrode, with S7 screw connector, glass body, for cream, milk, viscous material, open junction.
- KP100:** Flat membrane gel combined pH electrode with S7 screw connector, glass body, for skin, leather, paper.

Characteristics and dimensions of the probes on page WA-76.

- CP:** Extension cable 1.5m with BNC connectors on one side and S7 on the other side for electrode with S7 connector.
- CP5:** Extension cable 5m with BNC connectors on one side and S7 on the other side for electrode with S7 connector.
- CP10:** Extension cable 10m with BNC connector on one side and S7 on the other side, for electrode without cable.
- CP15:** Extension cable 15m with BNC connector on one side and S7 on the other side, for electrode without cable.
- CE:** S7 screw connector for pH electrode.
- BNC:** Female BNC for electrode extension.

ORP Electrodes

- KP90:** Redox Platinum electrode, with screw connector S7, electrolyte KCl 3M, body in glass.
- KP91:** Redox Platinum electrode with 1m cable, GEL filled, body in glass.
- Electrode dimensions and characteristics at page WA-76**

pH buffer solutions

- HD8642:** Buffer solution 4.01pH - 200cc.
- HD8672:** Buffer solution 6.86pH - 200cc.
- HD8692:** Buffer solution 9.18pH - 200cc.

Redox buffer solutions

- HDR220:** Redox buffer solution 220mV 0,5 l.
- HDR468:** Redox buffer solution 468mV 0,5 l.

Electrolyte solutions

- KCL 3M:** 100cc ready for use solution for refilling of the electrodes.

Cleaning and maintenance

- HD62PT:** Diaphragm cleaning (tiourea in HCl) - 500ml.
- HD62PP:** Protein cleaning (pepsin in HCl) - 500ml.
- HD62RF:** Regeneration (fluorhydric acid) - 100ml.
- HD62SC:** Solution for electrode preservation - 500ml.

Combined conductivity and temperature probes

- SP06T:** Combined conductivity and temperature 4-electrode cell in Platinum, body in Pocan. Cell constant $K = 0.7$. Measurement range $5\mu\text{S}/\text{cm} \dots 200\text{mS}/\text{cm}$, $0 \dots 90^\circ\text{C}$. max. pressure 5bar.
- SPT401.001:** Combined conductivity and temperature 2- electrode cell in stainless steel AISI 316. Cell constant $K = 0.01$. Measurement range $0.04\mu\text{S}/\text{cm} \dots 20\mu\text{S}/\text{cm}$, $0 \dots 120^\circ\text{C}$. **Measurement in closed-cell.** Max. pressure 5bar.
- SPT016:** Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant $K = 0.1$. Measurement range $0.1\mu\text{S}/\text{cm} \dots 500\mu\text{S}/\text{cm}$, $0 \dots 80^\circ\text{C}$. Max. pressure 5bar.
- SPT16:** Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant $K = 1$. Measurement range $10\mu\text{S}/\text{cm} \dots 10\text{mS}/\text{cm}$, $0 \dots 80^\circ\text{C}$. Max. pressure 5bar.
- SPT106:** Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant $K = 10$. Measurement range $500\mu\text{S}/\text{cm} \dots 200\text{mS}/\text{cm}$, $0 \dots 80^\circ\text{C}$. Max. pressure 5bar.

Electrode dimensions and characteristics at page WA-77

Standard conductivity calibration solutions

- HD8747:** Standard calibration solution 0.001mol/l equal to $147\mu\text{S}/\text{cm}$ @ 25°C - 200cc.
- HD8714:** Standard calibration solution 0.01mol/l equal to $1413\mu\text{S}/\text{cm}$ @ 25°C - 200cc.
- HD8712:** Standard calibration solution 0.1mol/l equal to $12880\mu\text{S}/\text{cm}$ @ 25°C - 200cc.
- HD87111:** Standard calibration solution 1mol/l equal to $111800\mu\text{S}/\text{cm}$ @ 25°C - 200cc.

Temperature probes complete with SICRAM module

- TP472I:** Wire wound Pt100 sensor, immersion probe. Stem \varnothing 3 mm, length 300 mm. Cable length 2 m.
- TP472I.0:** Thin film Pt100 sensor, immersion probe. Stem \varnothing 3 mm, length 230 mm. Cable length 2 m.
- TP473P.I:** Wire wound Pt100 sensor, penetration probe. Stem \varnothing 4mm, length 150 mm. Cable length 2 m.
- TP473P.0:** Thin film Pt100 sensor, penetration probe. Stem \varnothing 4mm, length 150 mm. Cable length 2 m.
- TP474C.I:** Wire wound Pt100 sensor, contact probe. Stem \varnothing 4mm, length 230mm, contact surface \varnothing 5mm. Cable length 2 m.
- TP474C.0:** Thin film Pt100 sensor, contact probe. Stem \varnothing 4mm, length 230mm, contact surface \varnothing 5mm. Cable length 2 m.
- TP475A.0:** Thin film Pt100 sensor, air probe. Stem \varnothing 4mm, length 230mm. Cable length 2 m.
- TP472I.5:** Thin film Pt100 sensor, penetration probe. Stem \varnothing 6mm, length 500 mm. Cable length 2 m.
- TP472I.10:** Thin film Pt100 sensor, penetration probe. Stem \varnothing 6mm, length 1000mm. Cable length 2 m.
- TP49A.0:** Thin film Pt100 sensor, immersion probe. Stem \varnothing 2,7mm, length 150mm. Cable length 2 m. Aluminium handle
- TP49AC.0:** Thin film Pt100 sensor, contact probe. Stem \varnothing 4mm, length 150mm. Cable length 2 m. Aluminium handle
- TP49AP.0:** Thin film Pt100 sensor, penetration probe. Stem \varnothing 2,7mm, length 150mm. Cable length 2 m. Aluminium handle
- TP875.I:** Wire wound Pt100 sensor, 150mm diameter globe-thermometer equipped with handle. Cable length 2 m.
- TP876.I:** Wire wound Pt100 sensor, 50mm diameter globe-thermometer equipped with handle. Cable length 2 m.
- TP870:** Thin film Pt100 sensor, immersion probe. Stem \varnothing 3 mm, length 70 mm. Cable length 2 m.
- TP878.0:** Thin film Pt100 sensor, contact probe for solar panels. Cable length 2 m.
- TP878.1.0:** Thin film Pt100 sensor, contact probe for solar panels. Cable length 5 m.
- TP879.0:** Thin film Pt100 sensor, penetration probe for compost. Stem \varnothing 8 mm, length 1000 mm. Cable length 2 m.

Temperature probes complete with TP47 module

- TP47.100.0:** Thin film Pt100 sensor, immersion probe. Stem \varnothing 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.
- TP47.1000.0:** Thin film Pt1000 sensor, immersion probe. Probe's Stem \varnothing 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.
- TP47:** Connector for Pt100 4-wire and Pt1000 2-wire probes without SICRAM module.
- TP87.100.0:** Thin film Pt100 sensor, immersion probe. Stem \varnothing 3mm, length 70mm. 4-wires connection cable with connector, length 1 m.
- TP87.1000.0:** Thin film Pt1000 sensor, immersion probe. Stem \varnothing 3mm, length 70mm. 2-wires connection cable with connector, length 1 m.





HD 2256.2 BENCH-TOP pH AND CONDUCTIVITY METER

The **HD2256.2** is a bench top instrument for electrochemical measures: **pH, conductivity and temperature**. It is fitted with a large backlit LCD display.

The **HD2256.2** measures **pH, mV, redox potential (ORP)** with pH, redox electrodes or electrodes with separate reference. **Conductivity and resistivity** in liquids, **total dissolved solids (TDS)** and **salinity** with combined 4-ring and 2-ring conductivity/temperature probes. The conductivity probes can have a direct input or with SICRAM module. The inputs are separate.

The instrument is fitted with input for the measurement of **temperature** with Pt100 or Pt1000 immersion, penetration or contact probes. The temperature probes are equipped with an automatic recognition module and factory calibration data are stored inside.

- The pH electrode calibration can be carried out on one or more (up to five) points and the calibration sequence can be chosen from a list of 13 buffers Temperature compensation can be automatic or manual.

- The conductivity probe calibration can be performed automatically with automatically detected conductivity calibration solutions: 147µS/cm, 1413µS/cm, 12880µS/cm or 111800µS/cm or manually with calibration solutions having different values.

- Conductivity, pH and temperature probes fitted with SICRAM module can store factory and calibration data inside.

The HD 2256.2 is a datalogger, it stores up to 2000 samples of: pH or mV, conductivity or resistivity or TDS or salinity and temperature.

The data can be transferred from the instrument connected to a PC via the multi-standard RS232C serial port and USB 2.0. The storing parameters can be configured using the menu.

The RS232C serial port can be used to transfer the acquired measurements to a 24 column portable printer in real time (HD40.1 or HD40.2).

The instruments equipped with **HD22BT** (Bluetooth) option can transfer data without any connection, to a PC or printer fitted with Bluetooth input (HD40.2) or through Bluetooth/RS232C converter.

The software DeltaLog11 allows instrument management and configuration, and data processing on PC.

The instruments have IP66 protection degree.

Technical characteristics HD2256.2

pH - mV - χ - Ω - TDS - NaCl - °C - °F

Instrument

Dimensions

(Length x Width x Height)

265x185x70mm

Weight

490g

Materials

ABS, rubber

Display

Back lighted, matrix point display.

240x64 points, visible area: 128x35mm

Operating conditions

Working temperature

-5 ... 50°C

Storage temperature

-25 ... 65°C

Working relative humidity

0 ... 90% R.H. without condensate

Protection degree

IP66

Power supply

Mains adapter (cod. SWD10) 12Vdc/1A

Auxiliary socket

For supplying of electrode holder with built-in stirrer HD22.2

Security of memorized data

Unlimited

Time

Date and hour

Real time schedule with backup battery 3.6V - 1/2AA

Accuracy

1min/month max drift


Measured values storing

Quantity

2000 screens

Storage interval

1s ... 999s



Calibration storage

Quantity Last 8 calibrations of each physical quantity

RS232C serial interface

Type RS232C electrically isolated
Baud rate Can be set from 1200 to 115200 baud
Data bit 8
Parity None
Stop bit 1
Flow Control Xon/Xoff
Length of serial cable Max 15m

Digital data Interface

USB 1.1 - 2.0 electrically isolated

Bluetooth Interface optional

Connections

Input for temperature probes with SICRAM modules⑤ 8-pole male DIN45326 connector
pH/mV input ① BNC female
Input for SICRAM module pH/ temperature ③ 8-pole male DIN45326 connector
2/ 4- electrode direct conductivity input ⑥ 8-pole male DIN45326 connector
Input conductivity electrodes with SICRAM module⑦ 8-pole male DIN45326 connector
Serial interface DB9 connector (9- pole male)
USB interface USB connector type B
Bluetooth Optional
Mains adapter 2-pole connector (Ø5.5mm-2.1mm).
Positive at centre
Outlet for power supply of electrode holder with built-in magnetic stirrer 2-pole connector (Ø5.5mm-2.1mm).
Positive at centre
(output 12Vdc/200mA max).

pH measurement by instrument

Measuring range -9.999...+19.999pH
Resolution 0.01 or 0.001pH selectable from menu
Accuracy $\pm 0.001\text{pH} \pm 1\text{digit}$
Input impedance $> 10^{12}\Omega$
Calibration error @25°C Offset $> 20\text{mV}$
Slope $> 63\text{mV/pH}$ or Slope $< 50\text{mV/pH}$
Sensitivity $> 106.5\%$ or Sensitivity $< 85\%$
Calibration points Up to 5 points with 13 automatically detected buffer solutions
Standard solutions automatically detected (@25°C) 1.679pH - 2.000pH - 4.000pH - 4.008pH - 4.010pH
6.860pH - 6.865pH - 7.000pH - 7.413pH - 7.648pH
9.180pH - 9.210pH - 10.010pH

mV measurement by instrument

Measuring range -1999.9...+1999.9mV
Resolution 0.1mV
Accuracy $\pm 0.1\text{mV} \pm 1\text{digit}$
Drift after 1 year 0.5mV/year

Conductivity measurement by instrument

Measuring range (Kcell=0.01) 0.000...1.999 $\mu\text{S/cm}$ Resolution 0.001 $\mu\text{S/cm}$
Measuring range (Kcell=0.1) 0.00...19.99 $\mu\text{S/cm}$ 0.01 $\mu\text{S/cm}$
Measuring range (K cell=1) 0.0...199.9 $\mu\text{S/cm}$ 0.1 $\mu\text{S/cm}$
200...1999 $\mu\text{S/cm}$ 1 $\mu\text{S/cm}$
2.00...19.99mS/cm 0.01mS/cm
20.0...199.9mS/cm 0.1mS/cm
Measuring range (Kcell=10) 200...1999mS/cm 1mS/cm
Accuracy (conductivity) $\pm 0.5\% \pm 1\text{digit}$

Measurement of resistivity by instrument

Measuring range (Kcell=0.01) Up to 1G $\Omega\cdot\text{cm}$ Resolution (*)
Measuring range (Kcell=0.1) Up to 100M $\Omega\cdot\text{cm}$ (*)
Measuring range (K cell=1) 5.0...199.9 $\Omega\cdot\text{cm}$ 0.1 $\Omega\cdot\text{cm}$
200...999 $\Omega\cdot\text{cm}$ 1 $\Omega\cdot\text{cm}$
1.00k...19.99k $\Omega\cdot\text{cm}$ 0.01k $\Omega\cdot\text{cm}$
20.0k...99.9k $\Omega\cdot\text{cm}$ 0.1k $\Omega\cdot\text{cm}$
100k...999k $\Omega\cdot\text{cm}$ 1k $\Omega\cdot\text{cm}$
1...10M $\Omega\cdot\text{cm}$ 1M $\Omega\cdot\text{cm}$
Measuring range (Kcell=10) 0.5...5.0 $\Omega\cdot\text{cm}$ 0.1 $\Omega\cdot\text{cm}$
Accuracy (resistivity) $\pm 0.5\% \pm 1\text{digit}$

Measurement of total dissolved solids (with coefficient $\chi/\text{TDS}=0.5$)

Measuring range (Kcell=0.01) 0.00...1.999mg/l 0.005mg/l
Measuring range (Kcell=0.1) 0.00...19.99mg/l 0.05mg/l
Measuring range (K cell=1) 0.0...199.9 mg/l 0.5 mg/l
200...1999 mg/l 1 mg/l
2.00...19.99 g/l 0.01 g/l
20.0...199.9 g/l 0.1 g/l
Measuring range (Kcell=10) 100...999 g/l 1 g/l
Accuracy (total dissolved solids) $\pm 0.5\% \pm 1\text{digit}$

Measurement of salinity by instrument

Measuring range 0.000...1.999g/l 1mg/l
2.00...19.99g/l 10mg/l
20.0...199.9 g/l 0.1 g/l
Accuracy (salinity) $\pm 0.5\% \pm 1\text{digit}$

Automatic/manual temperature compensation

0...100°C with $\alpha_T = 0.00...4.00\%/^{\circ}\text{C}$
Reference temperature 0...50°C
Conversion factor χ/TDS 0.4...0.8
Cell constant K (cm⁻¹) already set on instrument 0.01 - 0.1 - 0.5 - 0.7 - 1.0 - 10.0
Cell constants K(cm⁻¹) that can be set by user 0.01...20.00

Standard solutions automatically detected (@25°C)


147 $\mu\text{S/cm}$
1413 $\mu\text{S/cm}$
12880 $\mu\text{S/cm}$
111800 $\mu\text{S/cm}$

Measurement of temperature by instrument

Pt100 measuring range -50...+150°C
Pt1000 measuring range -50...+150°C
Resolution 0.1°C
Accuracy $\pm 0.1^{\circ}\text{C} \pm 1\text{digit}$
Drift after 1 year 0.1°C/year

(*) The resistivity measurement is obtained from the reciprocal of conductivity measurement. Close to the bottom of the scale, the indication of resistivity appears like reported in the table below:

K cell = 0.01 cm ⁻¹		K cell = 0.1 cm ⁻¹	
Conductivity ($\mu\text{S/cm}$)	Resistivity (M $\Omega\cdot\text{cm}$)	Conductivity ($\mu\text{S/cm}$)	Resistivity (M $\Omega\cdot\text{cm}$)
0.001 $\mu\text{S/cm}$	1000 M $\Omega\cdot\text{cm}$	0.01 $\mu\text{S/cm}$	100 M $\Omega\cdot\text{cm}$
0.002 $\mu\text{S/cm}$	500 M $\Omega\cdot\text{cm}$	0.02 $\mu\text{S/cm}$	50 M $\Omega\cdot\text{cm}$
0.003 $\mu\text{S/cm}$	333 M $\Omega\cdot\text{cm}$	0.03 $\mu\text{S/cm}$	33 M $\Omega\cdot\text{cm}$
0.004 $\mu\text{S/cm}$	250 M $\Omega\cdot\text{cm}$	0.04 $\mu\text{S/cm}$	25 M $\Omega\cdot\text{cm}$



TECHNICAL DATA OF PROBES AND MODULES EQUIPPED WITH INSTRUMENT **Temperature probes Pt100 sensor with SICRAM module**

Model	Type	Application field	Accuracy
TP472I	Immersion	-196°C...+500°C	±0.25°C (-196°C...+300°C) ±0.5°C (+300°C...+500°C)
TP472I.0 1/3 DIN Thin Film	Immersion	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP473P.I	Penetration	-50°C...+400°C	±0.25°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP473P.0 1/3 DIN Thin Film	Penetration	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP474C.I	Contact	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP474C.0 1/3 DIN Thin Film	Contact	-50°C...+300°C	±0.3°C (-50°C...+300°C)
TP475A.0 1/3 DIN Thin Film	Air	-50°C...+250°C	±0.3°C (-50°C...+250°C)
TP472I.5	Penetration	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP472I.10	Penetration	-50°C...+400°C	±0.30°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP49A.0 Class A Thin Film	Immersion	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AC.0 Class A Thin Film	Contact	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AP.0 Class A Thin Film	Penetration	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP875.I	Globe-thermometer Ø150mm	-30°C...+120°C	±0.25°C
TP876.I	Globe-thermometer Ø50mm	-30°C...+120°C	±0.25°C
TP87.0 1/3 DIN Thin Film	Immersion	-50°C...+200°C	±0.25°C
TP878.0 1/3 DIN Thin Film TP878.1.0 1/3 DIN Thin Film	Photovoltaic	+4°C...+85°C	±0.25°C
TP879.0 1/3 DIN Thin Film	Compost	-20°C...+120°C	±0.25°C

Common characteristics

Temperature drift @ 20°C 0.003%/°C

4 wires Pt100 and 2 wires Pt1000 Probes

Model	Type	Application field	Accuracy
TP47.100.0 1/3 DIN Thin Film	4 wires Pt100	-50...+250°C	1/3 DIN
TP47.1000.0 1/3 DIN Thin Film	2 wires Pt1000	-50...+250°C	1/3 DIN
TP87.100.0 1/3 DIN Thin Film	4 wires Pt100	-50...+200°C	1/3 DIN
TP87.1000.0 1/3 DIN Thin Film	2 wires Pt1000	-50...+200°C	1/3 DIN

Common features

Temperature drift @20°C

Pt100 0.003%/°C
Pt1000 0.005%/°C

ORDERING CODES

HD2256.2: The kit is composed of: instrument HD2256.2 for the measurement of pH - redox - conductivity - resistivity - TDS - salinity - temperature, **datalogger**, stabilized power supply at mains voltage 100-240Vac/12Vdc-1A., instructions manual and software DeltaLog11.

pH/mV electrodes, conductivity probes, probes, temperature probes, standard reference solutions for different measurement types, connection cables for pH electrodes with S7 connector, cables for data download to PC or printer have to be ordered separately.

Accessories

9CPRS232: Connection cable SubD female 9- pole for serial output RS232C.

CP22: USB 2.0 connection cable - connector type A - connector type B.

DeltaLog11: Software for download and management of the data on PC using Windows operating systems.

SWD10: Stabilized power supply at 100-240Vac/12Vdc -1A mains voltage.

HD40.1: Portable, serial input, 24 column thermal printer, 57mm paper width.

HD40.2: 24-column portable thermal printer, **Bluetooth and serial interface**, 57mm paper width, four NiMH 1.2V rechargeable batteries, SWD10 power supply, instruction manual, 5 thermal paper rolls. Requires the module HD22BT (**optional**) or the cable 9CPRS232 (**optional**).

HD22.2: Laboratory electrode holder composed of basis plate with incorporated magnetic stirrer, staff and replaceable electrode holder. Height max. 380mm. Powered by bench-top meters of the series HD22... with cable HD22.2.1 (**optional**) or supplier SWD10 (**optional**).

HD22.3: Laboratory electrode holder with metal basis plate. Flexible electrode holder for free positioning. For Ø 12mm probes.

HD22BT: Bluetooth module for wireless data transmission from instrument to PC. **The fitting of the module into the instrument is made exclusively by Delta Ohm, at the time of placing the order.**

TP47: Connector for Pt100 4-wire and Pt1000 2-wire probes without SICRAM module.

pH electrodes (Input ①)

KP 20: Gel pH combined electrode for general use, with S7 screw connector, EPOXY body.

KP 30: Gel pH combined electrode for general use, 1m cable with BNC, EPOXY body.

KP 50: Gel pH combined electrode, porous Teflon ring junction, suitable for emulsions, demineralised water and waste water with S7 screw connector, glass body.

KP 61: 3 diaphragm liquid filled pH combined electrode for wine, milk, cream, etc., S7 screw connector, liquid reference filling, glass body.

KP 62: 1 diaphragm gel pH combined electrode for general use, pure water, varnishes, gel filled, S7 screw connector, glass body.

KP 63: liquid filled pH combined electrode for general use, varnishes, 1m cable with BNC, glass body.

KP 64: Liquid filled pH combined electrode, Teflon ring diaphragm, for wine, varnishes, emulsions, S7 screw connector, glass body.

KP 70: Pointed gel combined pH microelectrode diam. 6 x L=70 mm., with S7 screw connector, EPOXY body, glass tip, open junction for meat and cheese.

KP 80: Pointed gel pH combined electrode, with S7 screw connector, glass body, for cream, milk, viscous material, open junction.

KP100: Flat membrane gel combined pH electrode with S7 screw connector, glass body, for skin, leather, paper.

Characteristics and dimensions of the probes on page WA-76.

pH electrodes with SICRAM module (Input ③)



KP63TS: Combined pH/temperature electrode with SICRAM module, body in Epoxy, Ag/AgCl sat KCl, one diaphragm, for general purpose.

SICRAM Module with S7 input for pH electrodes (Input ③)

pH 471.1: SICRAM module for pH electrodes with S7 standard connection, cable L=1m.

pH 471.2: SICRAM module for pH electrodes with S7 standard connection, cable L=2m.

pH 471.5: SICRAM module for pH electrodes with S7 standard connection, cable L=5m.



mV, Ω, T

ORP Electrodes (Inputs ① and ②)

KP90: Redox Platinum electrode, with screw connector S7, electrolyte KCl 3M, body in glass.

KP91: Redox Platinum electrode with 1m cable, GEL filled, body in glass.

CP: Extension cable 1.5m with BNC connectors on one side and S7 on the other side for electrode with S7 connector.

CP5: Extension cable 5m with BNC connectors on one side and S7 on the other side for electrode with S7 connector.

CP10: Extension cable 10m with BNC connector on one side and S7 on the other side, for electrode without cable.

CP15: Extension cable 15m with BNC connector on one side and S7 on the other side, for electrode without cable.

CE: S7 screw connector for pH electrode.

BNC: Female BNC for electrode extension.

pH buffer solutions

HD8642: Buffer solution 4.01pH - 200cc.

HD8672: Buffer solution 6.86pH - 200cc.

HD8692: Buffer solution 9.18pH - 200cc.

Redox buffer solutions

HDR220: Redox buffer solution 220mV 0,5 l.

HDR468: Redox buffer solution 468mV 0,5 l.

Electrolyte solutions

KCL 3M: 100cc ready for use solution for electrode refilling.

Cleaning and maintenance

HD62PT: Diaphragm cleaning (tiourea in HCl) - 500ml.

HD62PP: Protein cleaning (pepsin in HCl) - 500ml.

HD62RF: Regeneration (fluorhydric acid) - 100ml.

HD62SC: Solution for electrode preservation - 500ml.

Conductivity probes and combined conductivity and temperature probes

without SICRAM module (Input ②)

SP06T: Combined conductivity and temperature 4-electrode cell in Platinum, body in Pocan. Cell constant $K = 0.7$. Measurement range $5\mu\text{S}/\text{cm} \dots 200\text{mS}/\text{cm}$, $0 \dots 90^\circ\text{C}$. Max. working pressure 5bar.

SPT401.001: Combined conductivity and temperature 2- electrode cell in stainless steel AISI 316. Cell constant $K = 0.01$. Measurement range $0.04\mu\text{S}/\text{cm} \dots 20\mu\text{S}/\text{cm}$, $0 \dots 120^\circ\text{C}$. **Measurement in closed-cell.** Max. working pressure 5bar.

SPT01G: Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant $K = 0.1$. Measurement range $0.1\mu\text{S}/\text{cm} \dots 500\mu\text{S}/\text{cm}$, $0 \dots 80^\circ\text{C}$. Max. working pressure 5bar.

SPT1G: Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant $K = 1$. Measurement range $10\mu\text{S}/\text{cm} \dots 10\text{mS}/\text{cm}$, $0 \dots 80^\circ\text{C}$. Max. working pressure 5bar.

SPT10G: Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant $K = 10$. Measurement range $500\mu\text{S}/\text{cm} \dots 200\text{mS}/\text{cm}$, $0 \dots 80^\circ\text{C}$. Max. working pressure 5bar.

Combined conductivity / temperature probes with SICRAM module (Input ②)

SPT1GS: Combined conductivity /temperature 2-electrode Platinum- wire cell, body in glass with SICRAM module. Cell constant $K = 1$. Measuring range $10\mu\text{S}/\text{cm} \dots 10\text{mS}/\text{cm}$, $0 \dots 80^\circ\text{C}$. Max. working pressure 5bar.

pH electrode and conductivity probes characteristics at page WA-76

Standard conductivity calibration solutions

HD8747: Standard calibration solution 0.001mol/l equal to $147\mu\text{S}/\text{cm}$ @ 25°C - 200cc.

HD8714: Standard calibration solution 0.01mol/l equal to $1413\mu\text{S}/\text{cm}$ @ 25°C - 200cc.

HD8712: Standard calibration solution 0.1mol/l equal to $12880\mu\text{S}/\text{cm}$ @ 25°C - 200cc.

HD87111: Standard calibration solution 1mol/l equal to $111800\mu\text{S}/\text{cm}$ @ 25°C - 200cc.

Temperature probes complete with SICRAM module (Input ⑤)

TP472I: Wire wound Pt100 sensor, immersion probe. Stem \varnothing 3 mm, length 300 mm. Cable length 2 m.

TP472I.0: Thin film Pt100 sensor, immersion probe. Stem \varnothing 3 mm, length 230 mm. Cable length 2 m.

TP473P.I: Wire wound Pt100 sensor, penetration probe. Stem \varnothing 4mm, length 150 mm. Cable length 2 m.

TP473P.0: Thin film Pt100 sensor, penetration probe. Stem \varnothing 4mm, length 150 mm. Cable length 2 m.

TP474C.I: Wire wound Pt100 sensor, contact probe. Stem \varnothing 4mm, length 230mm, contact surface \varnothing 5mm. Cable length 2 m.

TP474C.0: Thin film Pt100 sensor, contact probe. Stem \varnothing 4mm, length 230mm, contact surface \varnothing 5mm. Cable length 2 m.

TP475A.0: Thin film Pt100 sensor, air probe. Stem \varnothing 4mm, length 230mm. Cable length 2 m.

TP472I.5: Thin film Pt100 sensor, penetration probe. Stem \varnothing 6mm, length 500 mm. Cable length 2 m.

TP472I.10: Thin film Pt100 sensor, penetration probe. Stem \varnothing 6mm, length 1000mm. Cable length 2 m.

TP49A.0: Thin film Pt100 sensor, immersion probe. Stem \varnothing 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AC.0: Thin film Pt100 sensor, contact probe. Stem \varnothing 4mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AP.0: Thin film Pt100 sensor, penetration probe. Stem \varnothing 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP875.I: Wire wound Pt100 sensor, 150mm diameter globe-thermometer equipped with handle. Cable length 2 m.

TP876.I: Wire wound Pt100 sensor, 50mm diameter globe-thermometer equipped with handle. Cable length 2 m.

TP87.0: Thin film Pt100 sensor, immersion probe. Stem \varnothing 3 mm, length 70 mm. Cable length 2 m.

TP878.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 2 m.

TP878.1.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 5 m.

TP879.0: Thin film Pt100 sensor, penetration probe for compost. Stem \varnothing 8 mm, length 1000 mm. Cable length 2 m.

Temperature probes complete with TP47 module

TP47.100.0: Thin film Pt100 sensor, immersion probe. Stem \varnothing 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP47.1000.0: Thin film Pt1000 sensor, immersion probe. Probe's Stem \varnothing 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP47: Connector for Pt100 4-wire and Pt1000 2-wire probes without SICRAM module.

TP87.100.0: Thin film Pt100 sensor, immersion probe. Stem \varnothing 3mm, length 70mm. 4-wires connection cable with connector, length 1 m.

TP87.1000.0: Thin film Pt1000 sensor, immersion probe. Stem \varnothing 3mm, length 70mm. 2-wires connection cable with connector, length 1 m.

Instrument

Dimensions (Length x Width x Height)
Weight
Materials
Display

220x120x55mm
460g (complete with batteries)
ABS, rubber
2x4½ characters plus symbols
visible area: 52x42mm

Operating conditions

Working temperature
Storage temperature
Working relative humidity
Protection degree

-5 ... 50°C
-25 ... 65°C
0 ... 90% RH without condensation
IP66

Power

Batteries
Autonomy (only batteries)
Mains (cod. **SWD10**)

3 batteries 1.5V type AA
100 hours with 1800mAh alkaline batteries
Output mains adapter 100-240Vac/ 12Vdc-1A

Security of memorized data

Unlimited

Storage of measured values

Type
Quantity

2000 pages of 18 samples each
36,000 sets of measures made up of [X - Ω or TDS or NaCl] and [°C- °F]

Selectable storage interval

1s, 5s, 10s, 15s, 30s, 1min, 2min, 5min, 10min, 15min, 20min, 30min and 1hour

Time

Date and hour
Accuracy

Schedule in real time
1min/month max drift

Serial interface RS232C

Type
Baud rate
Data bit
Parity
Stop bit
Flow Control
Serial cable length
Selectable print interval

RS232C electrically isolated
Can be set from 1200 to 38400 baud
8
None
1
Xon/Xoff
Max 15m
immediate or 1s, 5s, 10s, 15s, 30s, 1min, 2min, 5min, 10min, 15min, 20min, 30min and 1hour

USB Interface

Type

1.1 - 2.0 electrically isolated

Common connections to all models

Serial interface and USB
Mains adapter (cod. SWD10)

8-pole MiniDin connector
2-pole connector (positive at centre) 12Vdc/1A

Measurement connections

Input conductivity
Input for temperature probes
complete with TP47 modules

8-pole male DIN45326 connector
8-pole male DIN45326 connector

Measurement of conductivity by instrument

Measurement range (Kcell=0.01)
Measurement range (Kcell=0.1)
Measurement range (Kcell=1)

Measurement range	Resolution
0.000...1.999 μ S/cm	0.001 μ S/cm
0.00...19.99 μ S/cm	0.01 μ S/cm
0.0...199.9 μ S/cm	0.1 μ S/cm
200...1999 μ S/cm	1 μ S/cm
2.00...19.99mS/cm	0.01mS/cm
20.0...199.9mS/cm	0.1mS/cm
200...1999mS/cm	1mS/cm
Accuracy (conductivity)	$\pm 0.5\% \pm 1$ digit



HD 3406.2

BENCH-TOP CONDUCTIVITY METER

The **HD3406.2** is a bench top instrument for electrochemical measures: **conductivity and temperature**.

The displayed data can be stored (**datalogger**) and can be transferred to PC or serial printer thanks to the multi-standard serial port RS232C and USB2.0 and software DeltaLog9 (Vers.2.0 and subsequent ones). The storing and printing parameters can be set from menu.


The **HD3406.2** measures **conductivity, resistivity in liquids, total dissolved solids (TDS) and salinity** using combined 4-ring and 2-ring conductivity/temperature probes. Temperature is measured by Pt100 or Pt1000 immersion, penetration or contact probes.

The probe calibration can be performed automatically in one or more of the 147 μ S, 1413 μ S, 12880 μ S or 111800 μ S/cm conductivity calibration solutions.

The display shows continually the temperature in °C or °F and one selectable parameter according to the connected probe type, i.e. in case of conductivity probe it is possible to select between X, Ω or TDS or NaCl.

Other functions of this instrument include: Max, Min and Avg function, the Auto-HOLD function, the automatic turning off which can also be excluded.

The instruments have IP66 protection degree.



Measurement of resistivity by instrument

Measurement range (Kcell=0.01)	Up to 16Ω-cm
Measurement range (Kcell=0.1)	Up to 100MΩ-cm
Measurement range (Kcell=1)	5.0...199.9Ω-cm 200...999Ω-cm 1.00k...19.99kΩ-cm 20.0k...99.9kΩ-cm 100k...999kΩ-cm 1...10MΩ-cm
Measurement range (Kcell=10)	0.5...5.0Ω-cm
Accuracy (resistivity)	±0.5% ±1digit

Measurement of total dissolved solids (with coefficient χ /TDS=0.5)

Measurement range (Kcell=0.01)	0.00...1.999mg/l	0.005mg/l
Measurement range (Kcell=0.1)	0.00...19.99mg/l	0.05mg/l
Measurement range (Kcell=1)	0.0...199.9 mg/l 200...1999 mg/l 2.00...19.99 g/l 20.0...99.9 g/l	0.5 mg/l 1 mg/l 0.01 g/l 0.1 g/l
Measurement range (Kcell=10)	100...999 g/l	1 g/l
Accuracy (total dissolved solids)	±0.5% ±1digit	

Measurement of salinity

Measurement range	0.000...1.999g/l 2.00...19.99g/l 20.0...199.9g/l	1mg/l 10mg/l 0.1g/l
Accuracy (salinity)	±0.5% ±1digit	

Temperature measurement by instrument

Measurement range Pt100	-50...+200°C
Measurement range Pt1000	-50...+200°C
Resolution	0.1°C
Accuracy	±0.25°C
Drift after 1 year	0.1°C/year

Automatic/manual temperature compensation

Reference temperature	0...100°C with $\alpha_T = 0.00...4.00\%/^{\circ}\text{C}$
Conversion factor χ /TDS	20°C or 25°C selectable from menu
Cell constant K (cm ⁻¹)	0.4...0.8
	0.01 - 0.1 - 0.7 - 1.0 - 10.0

Standard solutions automatically detected (@25°C)

147μS/cm
1413μS/cm
12880μS/cm
111800μS/cm

(*) The resistivity measurement is obtained from the reciprocal of conductivity measurement. Close to the bottom of the scale, the indication of resistivity appears like reported in the table below:

K cell = 0.01 cm ⁻¹		K cell = 0.1 cm ⁻¹	
Conductivity (μS/cm)	Resistivity (MΩ-cm)	Conductivity (μS/cm)	Resistivity (MΩ-cm)
0.001 μS/cm	1000 MΩ-cm	0.01 μS/cm	100 MΩ-cm
0.002 μS/cm	500 MΩ-cm	0.02 μS/cm	50 MΩ-cm
0.003 μS/cm	333 MΩ-cm	0.03 μS/cm	33 MΩ-cm
0.004 μS/cm	250 MΩ-cm	0.04 μS/cm	25 MΩ-cm

ORDERING CODES

HD3406.2: The kit is composed of: instrument HD3406.2 **data logger**, for measurement of conductivity - resistivity - TDS - salinity - temperature, 3 1.5V alkaline batteries, operating manual and **DeltaLog9 version 2.0**.

Conductivity probes, temperature probes, standard reference solutions, cables for data download to PC or printer have to be ordered separately.

Accessories

HD2110CSNM: 8-pole connection cable Mini Din - Sub D 9-pole female for RS232C, for connection to PC without USB input.

HD2101/USB: Connection cable USB 2.0 connector type A - 8-pole Mini Din for connection to PC with USB input.

SWD10: Stabilized power supply at 100-240Vac/12Vdc-1A mains voltage.

HD40.1: Portable, serial input, 24 column thermal printer, 57mm paper width.

HD22.2: Laboratory electrode holder composed of basis plate with incorporated magnetic stirrer, staff and replaceable electrode holder. Height max. 380mm.

HD22.3: Laboratory electrode holder with metal basis plate. Flexible electrode holder for free positioning. For Ø 12mm probes.

TP47: Connector for Pt100 4-wire and Pt1000 2-wire probes without SICRAM module.

Combined conductivity and temperature probes

SP06T: Combined conductivity and temperature 4-electrode cell in Platinum, body in POCAN. Cell constant K = 0.7. Measurement range 5μS/cm...200mS/cm, 0...90°C. Max. working pressure 5bar.

SPT401.001: Combined conductivity and temperature 2-electrode cell in stainless steel AISI 316. Cell constant K = 0.01. Measurement range 0.04μS/cm...20μS/cm, 0...120°C. **Measurement in closed-cell.** Max. working pressure 5bar.

SPT01G: Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant K = 0.1. Measurement range 0.1μS/cm...500μS/cm, 0...80°C. Max. working pressure 5bar.

SPT1G: Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant K = 1. Measurement range 10μS/cm...10mS/cm, 0...80°C. Max. working pressure 5bar.

SPT10G: Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant K = 10. Measurement range 500μS/cm...200mS/cm, 0...80°C. Max. working pressure 5bar.

Conductivity probes characteristics at pag. WA-77

Standard conductivity calibration solutions

HD8747: Standard calibration solution 0.001mol/l equal to 147μS/cm @25°C - 200cc.

HD8714: Standard calibration solution 0.01mol/l equal to 1413μS/cm @25°C - 200cc.

HD8712: Standard calibration solution 0.1mol/l equal to 12880μS/cm @25°C - 200cc.

HD87111: Standard calibration solution 1mol/l equal to 111800μS/cm @25°C - 200cc.

Temperature probes complete with SICRAM module

TP472I: Wire wound Pt100 sensor, immersion probe. Stem Ø 3 mm, length 300 mm. Cable length 2 m.

TP472I.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 230 mm. Cable length 2 m.

TP473P.I: Wire wound Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP473P.0: Thin film Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP474C.I: Wire wound Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP474C.0: Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP475A.0: Thin film Pt100 sensor, air probe. Stem Ø 4mm, length 230mm. Cable length 2 m.

TP472I.5: Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 500 mm. Cable length 2 m.

TP472I.10: Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 1000mm. Cable length 2 m.

TP49A.0: Thin film Pt100 sensor, immersion probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AC.0: Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AP.0: Thin film Pt100 sensor, penetration probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP875.I: Wire wound Pt100 sensor, 150mm diameter globe-thermometer equipped with handle. Cable length 2 m.

TP876.I: Wire wound Pt100 sensor, 50mm diameter globe-thermometer equipped with handle. Cable length 2 m.

TP87.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 70 mm. Cable length 2 m.

TP878.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 2 m.

TP878.1.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 5 m.

TP879.0: Thin film Pt100 sensor, penetration probe for compost. Stem Ø 8 mm, length 1000 mm. Cable length 2 m.

Temperature probes without SICRAM module


TP47.100.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP47.1000.0: Thin film Pt1000 sensor, immersion probe. Probe's stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP47: Connector for Pt100 4-wire and Pt1000 2-wire probes without SICRAM module.

TP87.100.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 70mm. 4-wires connection cable with connector, length 1 m.

TP87.1000.0: Thin film Pt1000 sensor, immersion probe. Stem Ø 3mm, length 70mm. 2-wires connection cable with connector, length 1 m.





HD 2206.2 BENCH-TOP CONDUCTIVITY METER

The **HD2206.2** is a bench top instrument for electrochemical measures: **conductivity** and **temperature**. It is fitted with a large backlit LCD display. The **HD2206.2** measures **conductivity**, **resistivity** in liquids, **total dissolved solids** (TDS), and **salinity** with combined 4-ring and 2-ring conductivity/temperature probes. The conductivity probes can have a direct input or with SICRAM module. The inputs are separate. The instrument is fitted with input for the measurement of **temperature** with Pt100 or Pt1000 immersion, penetration or contact probes. The temperature probes are equipped with an automatic recognition module and factory calibration data are stored inside.

- The conductivity probe calibration can be performed automatically with automatically detected conductivity calibration solutions: 147µS/cm, 1413µS/cm, 12880µS/cm or 111800µS/cm or manually with calibration solutions having different values.
- Conductivity and temperature probes fitted with SICRAM module can store factory and calibration data inside.

The instrument HD2206.2 is a **data logger**, it can memorize up to 2,000 samples of data of conductivity (or resistivity or TDS or salinity) and temperature.

The data can be transferred from the instrument connected to a PC via the RS232C and USB 2.0 serial ports. The storing parameters can be configured using the menu.

The RS232C serial port can be used to transfer the acquired measurements to a 24 column portable printer in real time (HD40.1 or HD40.2). The instruments equipped with **HD22BT** (Bluetooth) option can transfer data without any connection to a PC or printer fitted with Bluetooth input or through Bluetooth/RS232C converter. The software DeltaLog11 allows instrument management and configuration, and data processing on PC.

The instruments have IP66 protection degree.

Technical characteristics HD2206.2

X - Ω - TDS - NaCl - °C - °F measurement

Instrument

Dimensions (Length x Width x Height)	265x185x70mm
Weight	490g
Materials	ABS, rubber
Display	Back lighted, matrix point display. 240x64 points, visible area: 128x35mm

Operating conditions

Working temperature	-5 ... 50°C
Storage temperature	-25 ... 65°C
Working relative humidity	0 ... 90% R.H. without condensate

Protection degree

IP66

Power

Mains adapter (cod. SWD10) 12Vdc/1A

Auxiliary socket

For supplying of electrode holder with built-in stirrer HD22.2

Security of memorized data

Unlimited

Time

Date and hour

Real time schedule with backup battery 3.6V - ½AA

Accuracy

1min/month max drift

Measured values storing

Quantity	2000 screens
Storage interval	1s ... 999s

Calibration storage


Quantity
Last 8 calibrations of each physical quantity

RS232C serial interface

Type	RS232C electrically isolated
Baud rate	Can be set from 1200 to 115200 baud
Data bit	8
Parity	None
Stop bit	1
Flow Control	Xon/Xoff
Length of serial cable	Max 15m

Data Interface

USB	1.1 - 2.0 electrically isolated
Bluetooth Interface	optional



Connections

Input for temperature probes with SICRAM module ⑤	8-pole male DIN45326 connector
2/4 ring direct conductivity input	8-pole male DIN45326 connector
Conductivity probe with SICRAM module input ⑦	8-pole male DIN45326 connector
Serial interface	DB9 connector (9- pole male)
USB interface	USB connector type B
Bluetooth	Optional
Mains adapter	2 -pole (Ø5.5mm-2.1mm). Positive at centre
Socket for power supply of electrode holder with built-in magnetic stirrer	2- pole connector (Ø5.5mm-2.1mm). Positive at centre (output 12Vdc/200mA max).

Measurement of conductivity by instrument

Measuring range (Kcell=0.01)	0.000...1.999µS/cm	Resolution 0.001µS/cm
Measuring range (Kcell=0.1)	0.00...19.99µS/cm	0.01µS/cm
Measuring range (K cell=1)	0.0...199.9µS/cm	0.1µS/cm
	200...1999µS/cm	1µS/cm
	2.00...19.99mS/cm	0.01mS/cm
	20.0...199.9mS/cm	0.1mS/cm
	200...1999mS/cm	1mS/cm
Measuring range (Kcell=10)		
Accuracy (conductivity)	±0.5% ±1digit	

Measurement of resistivity by instrument

Measuring range (Kcell=0.01)	Up to 1GΩ·cm	(*)
Measuring range (Kcell=0.1)	Up to 100MΩ·cm	(*)
Measuring range (K cell=1)	5.0...199.9Ω·cm	0.1Ω·cm
	200...999Ω·cm	1Ω·cm
	1.00k...19.99kΩ·cm	0.01kΩ·cm
	20.0k...99.9kΩ·cm	0.1kΩ·cm
	100k...999kΩ·cm	1kΩ·cm
	1...10MΩ·cm	1MΩ·cm
Measuring range (Kcell=10)	0.5...5.0Ω·cm	0.1Ω·cm
Accuracy (resistivity)	±0.5% ±1digit	

Measurement of total dissolved solids (with coefficient χ /TDS=0.5)

Measuring range (Kcell=0.01)	0.00...1.999mg/l	0.005mg/l
Measuring range (Kcell=0.1)	0.00...19.99mg/l	0.05mg/l
Measuring range (K cell=1)	0.0...199.9 mg/l	0.5 mg/l
	200...1999 mg/l	1 mg/l
	2.00...19.99 g/l	0.01 g/l
	20.0...199.9 g/l	0.1 g/l
Measurement range (Kcell=10)	100...999 g/l	1 g/l
Accuracy (total dissolved solids)	±0.5% ±1digit	

Measurement of salinity

Measuring range	0.000...1.999g/l	Resolution 1mg/l
	2.00...19.99g/l	10mg/l
	20.0...199.9 g/l	0.1 g/l
Accuracy (salinity)	±0.5% ±1digit	

Automatic/manual temperature compensation

Reference temperature	0...100°C with $\alpha_T = 0.00...4.00\%/^{\circ}\text{C}$
χ /TDS conversion factor	0...50°C
Cell constants K (cm^{-1}) already set on the instrument	0.4...0.8
Cell constants K (cm^{-1}) that can be set by user	0.01 - 0.1 - 0.5 - 0.7 - 1.0 - 10.0
	0.01...20.00

Standard solutions automatically detected (@25°C)

147µS/cm
1413µS/cm
12880µS/cm
111800µS/cm

Measurement of temperature by instrument

Pt100 measuring range	-50...+150°C
Pt1000 measuring range	-50...+150°C
Resolution	0.1°C
Accuracy	±0.1°C ±1digit
Drift after 1 year	0.1°C/year

(*) The resistivity measurement is obtained from the reciprocal of conductivity measurement. Close to the bottom of the scale, the indication of resistivity appears like reported in the table below:


K cell = 0.01 cm^{-1}		K cell = 0.1 cm^{-1}	
Conductivity (µS/cm)	Resistivity (MΩ·cm)	Conductivity (µS/cm)	Resistivity (MΩ·cm)
0.001 µS/cm	1000 MΩ·cm	0.01 µS/cm	100 MΩ·cm
0.002 µS/cm	500 MΩ·cm	0.02 µS/cm	50 MΩ·cm
0.003 µS/cm	333 MΩ·cm	0.03 µS/cm	33 MΩ·cm
0.004 µS/cm	250 MΩ·cm	0.04 µS/cm	25 MΩ·cm

TECHNICAL DATA OF PROBES AND MODULES EQUIPPED WITH INSTRUMENT Temperature probes Pt100 sensor with SICRAM module


Model	Type	Application field	Accuracy
TP472I	Immersion	-196°C...+500°C	±0.25°C (-196°C...+300°C) ±0.5°C (+300°C...+500°C)
TP472I.0 1/3 DIN Thin Film	Immersion	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP473PI	Penetration	-50°C...+400°C	±0.25°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP473PO 1/3 DIN Thin Film	Penetration	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP474CI	Contact	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP474C.0 1/3 DIN Thin Film	Contact	-50°C...+300°C	±0.3°C (-50°C...+300°C)
TP475A.0 1/3 DIN Thin Film	Air	-50°C...+250°C	±0.3°C (-50°C...+250°C)
TP472I.5	Penetration	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP472I.10	Penetration	-50°C...+400°C	±0.30°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP49A.0 Class A Thin Film	Immersion	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AC.0 Class A Thin Film	Contact	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AP.0 Class A Thin Film	Penetration	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP875.I	Globe-thermometer Ø150mm	-30°C...+120°C	±0.25°C
TP876.I	Globe-thermometer Ø50mm	-30°C...+120°C	±0.25°C
TP87.0 1/3 DIN Thin Film	Immersion	-50°C...+200°C	±0.25°C
TP878.0 1/3 DIN Thin Film TP878.1.0 1/3 DIN Thin Film	Photovoltaic	+4°C...+85°C	±0.25°C
TP879.0 1/3 DIN Thin Film	Compost	-20°C...+120°C	±0.25°C

Common characteristics

Temperature drift @ 20°C	0.003%/°C
--------------------------	-----------



χ



TDS

4 wires Pt100 and 2 wires Pt1000 Probes

Model	Type	Application field	Accuracy
TP47.100.0 1/3 DIN Thin Film	4 wires Pt100	-50...+250°C	1/3 DIN
TP47.1000.0 1/3 DIN Thin Film	2 wires Pt1000	-50...+250°C	1/3 DIN
TP87.100.0 1/3 DIN Thin Film	4 wires Pt100	-50...+200°C	1/3 DIN
TP87.1000.0 1/3 DIN Thin Film	2 wires Pt1000	-50...+200°C	1/3 DIN

Common features

Temperature drift @20°C

Pt100 0.003%/°C

Pt1000 0.005%/°C

ORDERING CODES

HD2206.2: The kit is composed of: instrument HD2206.2 for the measurement of conductivity - resistivity - TDS - salinity - temperature, **datalogger**, stabilized power supply at mains voltage 100-240Vac/12Vdc-1A., instructions manual and software DeltaLog11.

Conductivity probes, temperature probes, standard reference solutions, cables for data download to PC or printer have to be ordered separately.

ACCESSORIES

9CPRS232: Connection cable SubD female 9- pole for serial output RS232C.

CP22: USB 2.0 connection cable - connector type A - connector type B.

DeltaLog11: Software for download and management of the data on PC using Windows operating systems.

SWD10: Stabilized power supply at 100...240Vac/12Vdc-1A mains voltage.

HD40.1: 24-column portable thermal printer, **serial interface**, 57mm paper width, four NiMH 1.2V rechargeable batteries, SWD10 power supply, instruction manual, 5 thermal paper rolls. Requires the cable 9CPRS232 (**optional**).

HD40.2: 24-column portable thermal printer, **Bluetooth and serial interface**, 57mm paper width, four NiMH 1.2V rechargeable batteries, SWD10 power supply, instruction manual, 5 thermal paper rolls. Requires the module HD22BT (**optional**) or the cable 9CPRS232 (**optional**).

HD22.2: Laboratory electrode holder composed of basis plate with incorporated magnetic stirrer, staff and replaceable electrode holder. Height max. 380mm. Powered by bench-top meters of the series HD22... with cable HD22.2.1 (**optional**) or supplier SWD10 (**optional**).

HD22.3: Laboratory electrode holder with metal basis plate. Flexible electrode holder for free positioning. For Ø 12mm probes.

HD22BT: Bluetooth module for wireless data transmission from instrument to PC. **The fitting of the module into the instrument is made exclusively by Delta Ohm, at the time of placing the order.**

Conductivity probes and combined conductivity and temperature probes without SICRAM module (Input ②)

SP06T: Combined conductivity and temperature 4-electrode cell in Platinum, body in Pocan. Cell constant K = 0.7. Measurement range 5µS/cm ...200mS/cm, 0...90°C. Max. pressure 5bar.

SPT401.001: Combined conductivity and temperature 2- electrode cell in stainless steel AISI 316. Cell constant K = 0.01. Measurement range 0.04µS/cm ...20µS/cm, 0...120°C. Measurement in closed-cell. Max.pressure 5bar.

SPT01G: Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant K = 0.1. Measurement range 0.1µS/cm ...500µS/cm, 0...80°C. Max.pressure 5bar.

SPT1G: Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant K = 1. Measurement range 10µS/cm ...10mS/cm, 0...80°C. Max. pressure 5bar.

SPT10G: Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant K = 10. Measurement range 500µS/cm ...200mS/cm, 0...80°C. Max.pressure 5bar.

Probe dimensions and characteristics at page WA-77

Combined conductivity / temperature probes with SICRAM module (Input ③)

SPT1GS: Combined conductivity /temperature 2-electrode Platinum- wire cell, body in glass with SICRAM module. Cell constant K = 1. Measuring range 10µS/cm ...10mS/cm, 0...80°C. Max.pressure 5bar.

Probe dimensions and characteristics at page WA-77

Standard conductivity calibration solutions

HD8747: Standard calibration solution 0.001mol/l equal to 147µS/cm @25°C - 200cc.

HD8714: Standard calibration solution 0.01mol/l equal to 1413µS/cm @25°C - 200cc.

HD8712: Standard calibration solution 0.1mol/l equal to 12880µS/cm @25°C - 200cc.

HD87111: Standard calibration solution 1mol/l equal to 111800µS/cm @25°C - 200cc.

Temperature probes complete with SICRAM module

TP472I: Wire wound Pt100 sensor, immersion probe. Stem Ø 3 mm, length 300 mm. Cable length 2 m.

TP472I.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 230 mm. Cable length 2 m.

TP473P.I: Wire wound Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP473P.0: Thin film Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP474C.I: Wire wound Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP474C.0: Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP475A.0: Thin film Pt100 sensor, air probe. Stem Ø 4mm, length 230mm. Cable length 2 m.

TP472I.5: Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 500 mm. Cable length 2 m.

TP472I.10: Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 1000mm. Cable length 2 m.

TP49A.0: Thin film Pt100 sensor, immersion probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AC.0: Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AP.0: Thin film Pt100 sensor, penetration probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP875.I: Wire wound Pt100 sensor, 150mm diameter globe-thermometer equipped with handle. Cable length 2 m.

TP876.I: Wire wound Pt100 sensor, 50mm diameter globe-thermometer equipped with handle. Cable length 2 m.

TP87I.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 70 mm. Cable length 2 m.

TP878.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 2 m.

TP878.1.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 5 m.

TP879.0: Thin film Pt100 sensor, penetration probe for compost. Stem Ø 8 mm, length 1000 mm. Cable length 2 m.

Temperature probes without SICRAM module

TP47100.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP471000.0: Thin film Pt1000 sensor, immersion probe. Probe's Stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP47: Connector for Pt100 4-wire and Pt1000 2-wire probes without SICRAM module.

TP87100.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 70mm. 4-wires connection cable with connector, length 1 m.

TP871000.0: Thin film Pt1000 sensor, immersion probe. Stem Ø 3mm, length 70mm. 2-wires connection cable with connector, length 1 m.



HD22.3

HD40.1



HD 3409.2 BENCH-TOP DISSOLVED OXYGEN METER

The **HD3409.2** is a bench top instrument for electrochemical measures: **dissolved oxygen** and **temperature**.

The displayed data can be stored (**datalogger**) and can be transferred to PC or serial printer thanks to the multi-standard serial ports RS232C and USB2.0 and software DeltaLog9 (Vers.2.0 and subsequent ones). The storing and printing parameters can be set from menu.

The **HD3409.2** measures the **concentration** (in mg/l) of **dissolved Oxygen in liquids**, the **saturation index** (in %) and the **temperature** using SICRAM combined probes of polarographic type with two or three electrodes or galvanic type, and integrated temperature sensor. **Temperature** is measured by Pt100-SICRAM or direct 4 wire-immersion, penetration or contact probes.

Thanks to an internal pressure sensor, the instruments automatically compensate for barometric pressure. The instrument anticipates automatic compensation of the Oxygen probe membrane permeability and of the salinity of the liquid being examined. The dissolved Oxygen probe's quick calibration function guarantees timely correctness of the performed measurements.



The display shows continually the temperature in °C or °F and one selectable parameter according to the connected probe type. Printing and storage always include the temperature in °C or °F and one selectable parameter for each probe type.

Other common function of this instrument series include: Max, Min and Avg function, the Auto-HOLD function, the automatic turning off which can also be excluded.

The instruments have IP66 protection degree.

Technical characteristics HD3409.2

mg/l O₂, %O₂, mbar, °C/°F measurement

Instrument

Dimensions (Length x Width x Height)	220x120x55mm
Weight	460g (complete with batteries)
Materials	ABS, rubber
Display	2x4½ characters plus symbols visible area: 52x42mm

Operating conditions

Working temperature	-5 ... 50°C
Storage temperature	-25 ... 65°C
Working relative humidity	0 ... 90% RH without condensation
Protection degree	IP66

Power

Batteries	3 batteries 1.5V type AA
Autonomy (only batteries)	100 hours with 1800mAh alkaline batteries
Mains (cod. SWD10)	Output mains adapter 100-240Vac/ 12Vdc-1A

Security of memorized data

Unlimited

Storage of the measured values

Type	2000 pages of 9 samples each
Quantity	18,000 measures made up of the four parameters mg/l O ₂ , %O ₂ , mbar, [°C or °F]

Selectable storage interval

1s, 5s, 10s, 15s, 30s, 1min, 2min, 5min, 10min, 15min, 20min, 30min and 1hour

Time

Date and hour	Schedule in real time
Accuracy	1min/month max drift

Serial interface RS232C


Type	RS232C electrically isolated
Baud rate	Can be set from 1200 to 38400 baud
Data bit	8
Parity	None
Stop bit	1
Flow Control	Xon/Xoff
Serial cable length	Max 15m
Selectable print interval	immediate or 1s, 5s, 10s, 15s, 30s, 1min, 2min, 5min, 10min, 15min, 20min, 30min and 1hour

USB Interface

Type	1.1 - 2.0 electrically isolated
------	---------------------------------

Connections

Serial interface and USB	8-pole MiniDin connector
Mains adapter (cod. SWD10)	2-pole connector (positive at centre) 12Vdc/1A



Power absorbed with instrument off

Without dissolved oxygen probe	20µA
With dissolved oxygen probe	40µA

Measurement connections

Input for Oxygen probes	8-pole male DIN45326 connector
Input for temperature probes with SICRAM module or TP47 module	8-pole male DIN45326 connector

Measurement of the concentration of dissolved Oxygen

Measurement range	0.00...90.00mg/l
Resolution	0.01mg/l
Accuracy	±0.03mg/l±1digit (0...90%RH, 1013mbar, 20...25°C)

Measurement of the saturation index of dissolved Oxygen

Measurement range	0.0...600.0%
Resolution	0.1%
Accuracy	±0.3% ±1digit (in the range 0.0...199.9%) ±1% ±1digit (in the range 200.0...600.0%)

Automatic/manual temperature compensation

0...50°C

Measurement of barometric pressure

Measurement range	0.0...1100.0mbar
Resolution	0.1mbar
Accuracy	±2mbar±1digit between 18 and 25°C ±(2mbar+0.1mbar/°C) in the remaining range

Salinity setting

Setting range	0.0...70.0g/l
Resolution	0.1g/l

Temperature measurement with the sensor inside the dissolved Oxygen probe

Measurement range	0.0...+45.0°C
Resolution	0.1°C
Accuracy	±0.1°C ±1digit
Drift after 1 year	0.1°C/year

Temperature measurement by Instrument with Pt100 probe

Pt100 Measurement range	-200...+650°C
Resolution	0.1°C
Accuracy	±0.1°C ±1digit
Drift after 1 year	0.1°C/year

Ordering codes

HD3409.2: The kit is composed of: instrument HD3409.2 **data logger**, for the measurement of dissolved oxygen concentration - saturation index - temperature, calibrator HD9709/20 (for polarographic probe) or DO9709/21 (for galvanic probe), 3 1.5V alkaline batteries, operating manual and **DeltaLog9**.

Dissolved oxygen probes, temperature probes, standard reference solutions, connection cables, cables for data download to PC or printer have to be ordered separately.

Accessories

HD2110CSNM: 8-pole connection cable Mini Din - Sub D 9-pole female for RS232C, for connection to PC without USB input.

HD2101/USB: Connection cable USB 2.0 connector type A - 8-pole Mini Din for connection to PC with USB input.

SWD10: Stabilized power supply at 100-240Vac/12Vdc/1A mains voltage.

HD40.1: Portable, serial input, 24 column thermal printer, 57mm paper width, four NiMH 1.2V rechargeable batteries, SWD10 power supply, instruction manual, 5 thermal paper rolls. Requires the cable HD2110CSNM (**optional**).

HD22.2: Laboratory electrode holder composed of basis plate with incorporated magnetic stirrer, staff and replaceable electrode holder. Height max. 380mm. Powered by bench-top meters of the series HD22... with cable HD22.2.1 (**optional**) or supplier SWD10 (**optional**).

HD22.3: Laboratory electrode holder with metal basis plate. Flexible electrode holder for free positioning. For Ø 12mm probes.

TP47: Connector for Pt100 4-wire probes without SICRAM module.

Solutions

D09700: zero oxygen solution.

D09701: electrolyte solution for polarographic probes D09709 SS and D09709 SS.5.

D09701.1: electrolyte solution for galvanic probes D09709 SS.1 and D09709 SS.5.1.

Combined dissolved Oxygen/temperature probes

DO 9709 SS Polarographic combined oxygen and temperature probe with possibility of membrane replacement. Ø12mm x 120mm. 2m cable. The code includes: probe, 2 membranes, electrolyte solution and zero point solution.

DO 9709 SS.5 Polarographic combined oxygen and temperature probe with possibility of membrane replacement. Ø12mm x 120mm. 5m cable. The code includes: probe, 2 membranes, electrolyte solution and zero point solution.

DO 9709 SS.1 Galvanic combined oxygen and temperature probe with possibility of membrane replacement. Ø12mm x 120mm. Ø16mm tip with membrane. 2m cable. The code includes: probe, 2 membranes in total, electrolyte solution and zero point solution.

DO 9709 SS.5.1 Galvanic combined oxygen and temperature probe with possibility of membrane replacement. Ø12mm x 120mm. Ø16mm tip with membrane. 5m cable. The code includes: probe, 2 membranes in total, electrolyte solution and zero point solution.

Probes' specifications at page WA-79.

Accessories

DO 9709/20: Calibrator for polarographic probes DO 9709SS and DO 9709SS.5

DO 9709/21: Calibrator for galvanic probes DO 9709SS.1 and DO 9709SS.5.1

DO 9709 SSK: Kit of accessories for probes DO 9709SS and DO 9709SS.5: 3 membranes, zero point solution and electrolyte.

DO 9709/21K: Kit of accessories for probes DO 9709SS.1 and DO 9709SS.5.1: 3 membranes, zero point solution and electrolyte.

Temperature probes equipped with SICRAM module

TP472I: Wire wound Pt100 sensor, immersion probe. Stem Ø 3 mm, length 300 mm. Cable length 2 m.

TP472I.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 230 mm. Cable length 2 m.

TP473P.I: Wire wound Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP473P.0: Thin film Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP474C.I: Wire wound Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP474C.0: Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP475A.0: Thin film Pt100 sensor, air probe. Stem Ø 4mm, length 230mm. Cable length 2 m.

TP472I.5: Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 500 mm. Cable length 2 m.

TP472I.10: Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 1000mm. Cable length 2 m.

TP49A.0: Thin film Pt100 sensor, immersion probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AC.0: Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AP.0: Thin film Pt100 sensor, penetration probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP875.I: Wire wound Pt100 sensor, 150mm diameter globe-thermometer equipped with handle. Cable length 2 m.

TP876.I: Wire wound Pt100 sensor, 50mm diameter globe-thermometer equipped with handle. Cable length 2 m.

TP87.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 70 mm. Cable length 2 m.

TP878.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 2 m.

TP878.1.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 5 m.

TP879.0: Thin film Pt100 sensor, penetration probe for compost. Stem Ø 8 mm, length 1000 mm. Cable length 2 m.

Temperature probes without SICRAM module

TP47.100.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP47: Connector for Pt100 4-wire probes without SICRAM module.

TP87.100.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 70mm. 4-wires connection cable with connector, length 1 m.



mg/l

%sat

mbar

°C



HD 2259.2 BENCH-TOP DISSOLVED OXYGEN AND pH METER

The **HD2259.2** a bench top instrument for electrochemical measures: **pH**, **dissolved oxygen**, and **temperature**. It is fitted with a large backlit LCD display.

The **HD2259.2** measures **pH**, **mV**, **redox potential** (ORP) with pH, redox electrodes or electrodes with separate reference; the **concentration of dissolved oxygen** in liquids (in mg/l), and **saturation index** (in %), using SICRAM combined probes of polarographic type with two or three electrodes or galvanic type, and integrated temperature sensor.

The instrument fitted with an input for the measurement of **temperature** with Pt100 or Pt1000 immersion, penetration or contact probes. The temperature probes are equipped with an automatic recognition module and factory calibration data are stored inside.

- The pH electrode calibration can be carried out on up to five points and the calibration sequence can be chosen from a list of 13 buffers Temperature compensation can be automatic or manual.
- The dissolved Oxygen probe's quick calibration function guarantees timely correctness of the performed measurements.

- Conductivity, dissolved oxygen and temperature probes fitted with SICRAM module can store factory and calibration data inside.

The instrument HD2259.2 is a **datalogger**, it can memorize up to 2,000 samples of data: pH or mV, concentration of dissolved oxygen or saturation index and saturation index and temperature:

The data can be transferred from the instrument connected to a PC via the RS232C and USB 2.0 serial ports. The storing parameters can be configured using the menu.

The RS232C serial port can be used to transfer the acquired measurements to a 24 column portable printer in real time (HD 40.1, HD 40.2).

The instruments equipped with **HD22BT** (Bluetooth) option can transfer data without any connection to a PC or printer fitted with Bluetooth input (HD40.2) or through Bluetooth/RS232C converter.

The software DeltaLog11 allows instrument management and configuration, and data processing on PC.

The instruments have IP66 protection degree.

Technical characteristics HD2259.2

pH - mV - mg/l O₂ - %O₂ - mbar - °C - °F measurement

Instrument

Dimensions (Length x Width x Height)	265x185x70mm
Weight	490g
Materials	ABS, rubber
Display	Back lighted, matrix point display. 240x64 points, visible area: 128x35mm

Operating conditions

Working temperature	-5 ... 50°C
Storage temperature	-25 ... 65°C
Working relative humidity	0 ... 90% R.H. without condensate
Protection degree	IP66

Power

Mains adapter (cod. SWD10) 12Vdc/1A

Auxiliary output socket

For supplying of electrode holder with built-in stirrer HD22.2

Security of memorized data

Unlimited

Time

Date and hour	Real time schedule with backup battery 3.6V - ½AA
Accuracy	1min/month max drift

Measured values storing


Quantity	2000 screens
Storage interval	1s ... 999s

Calibration storage

Quantity	Last 8 calibrations of each physical quantity
----------	---

RS232C serial interface

Type	RS232C electrically isolated
------	------------------------------



Baud rate	Can be set from 1200 to 115200 baud
Data bit	8
Parity	None
Stop bit	1
Flow Control	Xon/Xoff
Length of serial cable	Max 15m

USB Interface

Type	1.1 - 2.0 electrically isolated
Bluetooth Interface	optional

Connections

Input for temperature probes with SICRAM modules ⑤	8-pole male DIN45326 connector
pH/mV inputs ①	BNC female
Input SICRAM module	8-pole male DIN45326 connector
pH/ temperature probes ③	
Input dissolved oxygen ⑥	8-pole male DIN45326 connector
Serial interface	DB9 connector (9- pole male)
USB interface	USB connector type B
Bluetooth	Optional
Mains adapter	2-pole connector (Ø5.5mm-2.1mm). Positive at centre
Outlet for power supply of electrode holder with built-in magnetic stirrer	2- pole connector (Ø5.5mm-2.1mm). Positive at centre (output 12Vdc/200mA max).

Measurement of pH by instrument

Measuring range	-9.999...+19.999pH
Resolution	0.01 o 0.001pH selectable from menu
Accuracy	±.001pH ±1digit
Input impedance	>10 ¹² Ω
Calibration error @25°C	Offset > 20mV Slope > 63mV/pH o Slope < 50mV/pH Sensitivity < 85% or sensitivity < 85%
Calibration points	Up to 5 points with 13 automatically detected buffer solutions
Automatically detected pH standard solutions (@25°C)	1.679pH - 2.000pH - 4.000pH - 4.008pH 4.010pH - 6.860pH - 6.865pH - 7.000pH 7.413pH - 7.648pH - 9.180pH - 9.210pH 10.010pH

mV measurement by instrument

Measuring range	-1999.9...+1999.9mV
Resolution	0.1mV
Accuracy	±0.1mV ±1digit
Drift after 1 year	0.5mV/year

Measurement of dissolved oxygen by instrument

Measuring range	0.00...90.00mg/l
Resolution	0.01mg/l
Accuracy	±0.03mg/l±1digit (0...90%RH,±1013mbar, 20...25°C)

Measurement of saturation index of dissolved oxygen

Measuring range	0.0...600.0%
Resolution	0.1%
Accuracy	±0.3% ±1digit (in the range 0.0...199.9%) ±1% ±1digit (in the range 200.0...600.0%)

Automatic temperature compensation

	0...50°C
--	----------

Measurement of barometric pressure

Measuring range	0.0...1100.0mbar
Resolution	0.1mbar
Accuracy	±2mbar±1digit between 18 and 25°C ±(2mbar+0.1mbar/°C) in the remaining range

Salinity setting

Setting range	0.0...70.0g/l
Resolution	0.1g/l

Temperature measurement with the sensor inside the dissolved oxygen probe

Measuring range	0.0...50.0°C
Resolution	0.1°C
Accuracy	±0.1°C ±1digit
Drift after 1 year	0.1°C/year

Measurement of temperature by instrument

Pt100 measuring range	-50...+150°C
Pt1000 measuring range	-50...+150°C
Resolution	0.1°C
Accuracy	±0.1°C ±1digit
Drift after 1 year	0.1°C/year

Ordering codes

HD2259.2: The kit is composed of: instrument HD2259.2 for the measurement of pH - redox - concentration of dissolved oxygen, saturation index - temperature, **datalogger**, stabilized power supply at mains voltage 100-240Vac/12Vdc-1A, calibrator HD9709/20 (for polarographic probe) or D09709/21 (for galvanic probe), instructions manual and software DeltaLog11.

pH/mV electrodes, dissolved oxygen probes, temperature probes, standard reference solutions for different measurement types, connection cables for pH electrodes with S7 connector, cables for data download to PC or printer have to be ordered separately.

Accessories

9CPRS232: Connection cable SubD female 9- pole for serial output RS232C.

CP22: USB 2.0 connection cable - connector type A - connector type B.

DeltaLog11: Software for download and management of the data on PC using Windows operating systems.

SWD10: Stabilized power supply at 100-240Vac/12Vdc/1A mains voltage.

HD40.1: 24-column portable thermal printer, serial interface, 57mm paper width, four NiMH 1.2V rechargeable batteries, SWD10 power supply, instructions manual, 5 thermal paper rolls. Requires the cable 9CPRS232 (optional).

HD40.2: 24-column portable thermal printer, **Bluetooth and serial interface**, 57mm paper width, four NiMH 1.2V rechargeable batteries, SWD10 power supply, instructions manual, 5 thermal paper rolls. Requires the module HD22BT (**optional**) or the cable 9CPRS232 (**optional**).

HD22.2: Laboratory electrode holder composed of basis plate with incorporated magnetic stirrer, staff and replaceable electrode holder. Height max. 380mm. Powered by bench-top meters of the series HD22... with cable HD22.2.1 (**optional**) or supplier SWD10 (**optional**).

HD22.3: Laboratory electrode holder with metal basis plate. Flexible electrode holder for free positioning. For Ø 12mm probes.

HD22BT: Bluetooth module for wireless data transmission from instrument to PC. **The fitting of the module into the instrument is made exclusively by Delta Ohm, at the time of placing the order.**

TP47: Connector for Pt100 4-wire and Pt1000 2-wire probes without SICRAM module.

Accessories

pH electrodes without SICRAM module (Inputs ①)

KP 20: Gel pH combined electrode for general use, with S7 screw connector, EPOXY body.

KP 30: Gel pH combined electrode for general use, 1m cable with BNC, EPOXY body.

KP 50: Gel pH combined electrode, porous Teflon ring junction, suitable for emulsions, demineralised water and waste water with S7 screw connector, glass body.

KP 61: 3 diaphragm liquid filled pH combined electrode for wine, milk, cream, etc., S7 screw connector, liquid reference filling, glass body.


KP 62: 1 diaphragm gel pH combined electrode for general use, pure water, varnishes, gel filled, S7 screw connector, glass body.

KP 63: liquid filled pH combined electrode for general use, varnishes, 1m cable with BNC, glass body.


KP 64: Liquid filled pH combined electrode, Teflon ring diaphragm, for wine, varnishes, emulsions, S7 screw connector, glass body.

KP 70: Pointed gel combined pH microelectrode diam. 6 x L=70 mm., with S7 screw connector, EPOXY body, glass tip, open junction for meat and cheese.

KP 80: Pointed gel pH combined electrode, with S7 screw connector, glass body, for cream, milk, viscous material, open junction.



pH



mg/l

KP100: Flat membrane gel combined pH electrode with S7 screw connector, glass body, for skin, leather, paper.

Characteristics and dimensions of the probes on page WA-76.

pH electrodes with SICRAM module (Input ③)


KP63TS: Combined pH/temperature electrode with SICRAM module, body in glass, 1m cable, 1 diaphragm, for general use, internal liquid reference.

SICRAM Module with S7 input for pH electrodes (Input ③)

pH 471.1: SICRAM module for pH electrodes with S7 standard connection, cable L=1m.

pH 471.2: SICRAM module for pH electrodes with S7 standard connection, cable L=2m.

pH 471.5: SICRAM module for pH electrodes with S7 standard connection, cable L=5m.



ORP Electrodes (Inputs ① and ②)

KP90: Redox Platinum electrode, with screw connector S7, electrolyte KCl 3M, body in glass.

KP91: Redox Platinum electrode with 1m cable and BNC connector, GEL filled, body in Epoxy.

Electrode dimensions and characteristics at page WA-76

CP: Extension cable 1.5m with BNC connectors on one side and S7 on the other side for electrode with S7 connector.

CP5: Extension cable 5m with BNC connectors on one side and S7 on the other side for electrode with S7 connector.

CP 10: 10m extension cable with BNC/S7 connector for electrode without cable, thread S7.

CP 15: 15m extension cable with BNC/S7 connector for electrode without cable, thread S7.

CE: S7 screw connector for pH electrode.

BNC: Female BNC for cable extension.

pH buffer solutions

HD8642: Buffer solution 4.01pH - 200cc.

HD8672: Buffer solution 6.86pH - 200cc.

HD8692: Buffer solution 9.18pH - 200cc.

Redox buffer solutions

HDR220: Redox buffer solution 220mV 0,5 l.

HDR468: Redox buffer solution 468mV 0,5 l.

Electrolyte solutions

KCL 3M: Ready for use solution for refilling of the electrodes.

Cleaning and maintenance

HD62PT: Diaphragm cleaning (tiourea in HCl) - 500ml.

HD62PP: Protein cleaning (pepsin in HCl) - 500ml.

HD62RF: Regeneration (fluorhydric acid) - 100ml.

HD62SC: Solution for electrode preservation - 500ml.

Solutions

D09700: zero oxygen solution.

D09701: electrolyte solution for polarographic probes D09709 SS and D09709 SS.5.

D09701.1: electrolyte solution for galvanic probes D09709 SS.1 and D09709 SS.5.1.

Combined dissolved Oxygen/temperature probes

D0 9709 SS Polarographic combined oxygen and temperature probe with possibility of membrane replacement. Ø12mm x 120mm. 2m cable. The code includes: probe, 2 membranes, electrolyte solution and zero point solution.

D0 9709 SS.5 Polarographic combined oxygen and temperature probe with possibility of membrane replacement. Ø12mm x 120mm. 5m cable. The code includes: probe, 2 membranes, electrolyte solution and zero point solution.

D0 9709 SS.1 Galvanic combined oxygen and temperature probe with possibility of membrane replacement. Ø12mm x 120mm. Ø16mm tip with membrane. 2m cable. The code includes: probe, 2 membranes in total, electrolyte solution and zero point solution.

D0 9709 SS.5.1 Galvanic combined oxygen and temperature probe with possibility of membrane replacement. Ø12mm x 120mm. Ø16mm tip with membrane. 5m cable. The code includes: probe, 2 membranes in total, electrolyte solution and zero point solution.

Probes' specifications at page WA-79.

Accessories

D0 9709/20: Calibrator for polarographic probes D0 9709SS and D0 9709SS.5

D0 9709/21: Calibrator for galvanic probes D0 9709SS.1 and D0 9709SS.5.1

D0 9709 SSK: Kit of accessories for probes D0 9709SS and D0 9709SS.5: 3 membranes, zero

point solution and electrolyte.

D0 9709/21K: Kit of accessories for probes D0 9709SS.1 and D0 9709SS.5.1: 3 membranes, zero point solution and electrolyte.

Temperature probes equipped with SICRAM module

TP472I: Wire wound Pt100 sensor, immersion probe. Stem Ø 3 mm, length 300 mm. Cable length 2 m.

TP472I.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 230 mm. Cable length 2 m.

TP473P.I: Wire wound Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP473P.0: Thin film Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP474C.I: Wire wound Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP474C.0: Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP475A.0: Thin film Pt100 sensor, air probe. Stem Ø 4mm, length 230mm. Cable length 2 m.

TP472I.5: Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 500 mm. Cable length 2 m.

TP472I.10: Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 1000mm. Cable length 2 m.

TP49A.0: Thin film Pt100 sensor, immersion probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AC.0: Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AP.0: Thin film Pt100 sensor, penetration probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP875.I: Wire wound Pt100 sensor, 150mm diameter globe-thermometer equipped with handle. Cable length 2 m.

TP876.I: Wire wound Pt100 sensor, 50mm diameter globe-thermometer equipped with handle. Cable length 2 m.

TP87I.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 70 mm. Cable length 2 m.

TP878.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 2 m.

TP878.1.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 5 m.

TP879.0: Thin film Pt100 sensor, penetration probe for compost. Stem Ø 8 mm, length 1000 mm. Cable length 2 m.

Temperature probes without SICRAM module

TP47.100.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP47.1000.0: Thin film Pt1000 sensor, immersion probe. Probe's Stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP47: Connector for Pt100 4-wire and Pt1000 2-wire probes without SICRAM module.

TP87.100.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 70mm. 4-wires connection cable with connector, length 1 m.

TP87.1000.0: Thin film Pt1000 sensor, immersion probe. Stem Ø 3mm, length 70mm. 2-wires connection cable with connector, length 1 m.



The instrument is fitted with an input for the measurement of **temperature** with Pt100 or Pt1000 immersion, penetration or contact probes. The temperature probes are equipped with an automatic recognition module and factory calibration data are stored inside.

- The pH electrode calibration can be carried out on up to five points and the calibration sequence can be chosen from a list of 13 buffers. Temperature compensation can be automatic or manual.
- The conductivity probe calibration can be performed automatically with automatically detected conductivity calibration solutions: 147µS/cm, 1413µS/cm, 12880µS/cm or 111800µS/cm or manually with calibration solutions having different values.
- The dissolved Oxygen probe's quick calibration function guarantees timely correctness of the performed measurements.
- Conductivity, dissolved oxygen and temperature probes fitted with SICRAM module can store factory and calibration data inside.

The instrument HD22569.2 is a **datalogger**, it can memorize up to 2,000 samples of data: pH or mV, conductivity or resistivity or TDS or salinity, concentration of dissolved oxygen and temperature.

The data can be transferred from the instrument connected to a PC via the RS232C or USB 2.0 serial ports. The storing parameters can be configured using the menu.

The RS232C serial port can be used to transfer the acquired measurements to a 24 column portable printer in real time (HD40.1 or HD40.2).

The instrument equipped with **HD22BT** (Bluetooth) option can transfer data without any connection to a PC or printer fitted with Bluetooth input or through Bluetooth/RS232C converter.

The software DeltaLog11 allows instrument management and configuration, and data processing on PC.

The instruments have IP66 protection degree.

Technical characteristics of HD22569.2

pH - mV - χ - Ω - TDS - NaCl - mg/l O_2 - % O_2 - mbar - $^{\circ}C$ - $^{\circ}F$ measurement

Instrument

Dimensions (Length x Width x Height)	265x185x70mm
Weight	490g
Materials	ABS, rubber
Display	Back lighted, matrix point display. 240x64 points, visible area: 128x35mm

Operating conditions

Working temperature	-5 ... 50 $^{\circ}C$
Storage temperature	-25 ... 65 $^{\circ}C$
Working relative humidity	0 ... 90% R.H. without condensate
Protection degree	IP66

Power

Mains adapter (cod. SWD10) 12Vdc/1A

Auxiliary socket

For supplying of electrode holder with built-in stirrer HD22.2

Security of memorized data

Unlimited




HD 22569.2

BENCH-TOP METER FOR pH - CONDUCTIVITY - DISSOLVED OXYGEN

The instrument **HD22569.2** is a bench top instrument for electrochemical measures: **pH, conductivity, dissolved oxygen**, and **temperature**. It is fitted with a large backlit LCD display.

The **HD22569.2** measures **pH, mV, redox potential (ORP)** with pH, redox electrodes or electrodes with separate reference; **conductivity, resistivity** in liquids, **total dissolved solids (TDS)** and **salinity** with combined 4-ring and 2-ring conductivity/temperature probes with direct input or SICRAM module; **concentration of dissolved oxygen** in liquids (in mg/l) and **saturation index** (in %), using SICRAM combined probes of polarographic type with two or three electrodes or galvanic type, and integrated temperature sensor.



Time		Measurement of resistivity by instrument		Resolution
Date and hour	Real time schedule with backup battery	Measurement range (K cell=0.01)	Up to 1GΩcm	
	3.6V - ½AA	Measurement range (K cell=0.1)	Up to 100MΩ·cm	(*)
Accuracy	1min/month max drift	Measurement range (K cell=1)	5.0...199.9Ω·cm	0.1Ωcm
			200...999Ω·cm	1Ω·cm
Measured values storing			1.00k...19.99kΩ·cm	0.01kΩcm
Quantity	2000 screens		20.0k...99.9kΩ·cm	0.1kΩcm
Storage interval	1s ... 999s		100k...999kΩ·cm	1kΩcm
			1...10MΩ·cm	1MΩ·cm
Calibration storage		Measurement range (K cell=10)	0.5...5.0Ω·cm	0.1Ωcm
Quantity	Last 8 calibrations of each physical quantity	Accuracy (resistivity) instrument	±0.5% ±1digit	
RS232C serial interface		Measurement of total dissolved solids		Resolution
Type	RS232C electrically isolated	<i>(with coefficient χ/TDS=0.5)</i>		
Baud rate	Can be set from 1200 to 115200 baud	Measurement range (K cell=0.01)	0.00...1.999mg/l	0.005mg/l
Data bit	8	Measurement range (K cell=0.1)	0.00...19.99mg/l	0.05mg/l
Parity	None	Measurement range (K cell=1)	0.0...199.9 mg/l	0.5 mg/l
Stop bit	1		200...1999 mg/l	1 mg/l
Flow Control	Xon/Xoff		2.00...19.99 g/l	0.01 g/l
Length of serial cable	Max 15m		20.0...199.9 g/l	0.1 g/l
			100...999 g/l	1 g/l
Serial data Interface		Measurement range (K cell=10)		
USB	1.1 - 2.0 electrically isolated	Accuracy (total dissolved solids)		
Bluetooth	optional	instrument	±0.5% ±1digit	
Connections		Measurement of salinity		Resolution
Input for temperature probes with SICRAM modules⑤	8-pole male DIN45326 connector	Measurement range	0.000...1.999g/l	1mg/l
pH/mV input ①	BNC female		2.00...19.99g/l	10mg/l
Input for SICRAM module pH/ temperature ③	8-pole male DIN45326 connector		20.0...199.9 g/l	0.1 g/l
2/ 4- electrode direct conductivity input ⑧	8-pole male DIN45326 connector	Accuracy (salinity) instrument	±0.5% ±1digit	
Conductivity SICRAM module input ⑦	8-pole male DIN45326 connector			
Dissolved Oxygen input ⑥	8-pole male DIN45326 connector	Automatic/manual temperature compensation		
Serial interface	DB9 connector (9- pole male)		0...100°C with $\alpha_T = 0.00...4.00\%/^{\circ}\text{C}$	
USB interface	USB connector type B	Reference temperature	0...50°C	
Bluetooth	Optional	Conversion factor χ/TDS	0.4...0.8	
Mains adapter	2- pole connector (Ø5.5mm-2.1mm). Positive at centre	Admitted cell constants K (cm⁻¹)	0.01- 0.1 - 0.5 - 0.7 - 1.0 - 10.0	
		Cell constants K (cm⁻¹)	0.01...20.00	
Outlet for power supply of electrode holder with built-in magnetic stirrer	2-pole connector (Ø5.5mm-2.1mm). Positive at centre (output 12Vdc/200mA max).	that can be set by user		
		Automatically detected standard solutions (@25°C)		
			147µS/cm	
			1413µS/cm	
			12880µS/cm	
			111800µS/cm	
Measurement of pH by instrument		Measurement of concentration of dissolved oxygen		
Measuring range	-9.999...+19.999pH	Measurement range	0.00...90.00mg/l	
Resolution	0.01 or 0.001pH selectable from menu	Resolution	0.01mg/l	
Accuracy	±0.001pH ±1digit	Accuracy instrument	±0.03mg/l ±1digit (0...90%,1013mbar, 20...25°C)	
Input impedance	>10 ¹² Ω			
Calibration error @25°C	Offset > 20mV Slope > 63mV/pH or Slope < 50mV/pH Sensitivity > 106.5% or Sensitivity < 85%	Measurement of saturation index of dissolved oxygen		
Calibration points	Up to 5 points from a list of 8 automatically detected buffers	Measurement range	0.0...600.0%	
Temperature compensation	-50...150°C	Resolution	0.1%	
Automatically detected standard solutions @25°C	1.679pH - 4.000pH - 4.010pH 6.860pH - 7.000pH - 7.648pH 9.180pH - 10.010pH	Accuracy instrument	±0.3% ±1digit (in the range 0.0...199.9%) ±1% ±1digit (in the range 200.0...600.0%)	
Measurement of mV by instrument		Measurement of barometric pressure		
Measuring range	-1999.9...+1999.9mV	Measuring range	0.0...1100.0mbar	
Resolution	0.1mV	Resolution	0.1mbar	
Accuracy	±0.1mV ±1digit	Accuracy	±2mbar±1digit between 18 and 25°C ±(2mbar+0.1mbar/°C) in the remaining range	
Drift after 1 year	0.5mV/year	Salinity setting		
		Setting	directly from menu or automatically by conductivity measurement	
Measurement of conductivity by instrument		Setting range	0.0...70.0g/l	
Measurement range (K cell=0.01)	0.000...1.999µS/cm	Resolution	0.1g/l	
Measurement range (K cell=0.1)	0.00...19.99µS/cm			
Measurement range (K cell=1)	0.0...199.9µS/cm	Temperature measurement with the sensor inside the O₂ probe		
	200...1999µS/cm	Measurement range	0.0...50.0°C	
	2.00...19.99mS/cm	Resolution	0.1°C	
	20.0...199.9mS/cm	Accuracy instrument	±0.1°C ±1digit	
Measurement range (K cell=10)	200...1999mS/cm	Drift after 1 year	0.1°C/year	
Accuracy (conductivity) instrument	±0.5% ±1digit	Automatic temperature compensation	0...50°C	

Measurement of temperature by instrument

Pt100 Measurement range	-50...+150°C
Pt1000 Measurement range	-50...+150°C
Resolution	0.1°C
Accuracy instrument	±0.1°C ±1digit
Drift after 1 year	0.1°C/year

(*) The resistivity measurement is obtained from the reciprocal of conductivity measurement. Close to the bottom of the scale, the indication of resistivity appears like reported in the table below:

K cell = 0.01 cm ⁻¹		K cell = 0.1 cm ⁻¹	
Conductivity (μS/cm)	Resistivity (MΩ·cm)	Conductivity (μS/cm)	Resistivity (MΩ·cm)
0.001 μS/cm	1000 MΩ·cm	0.01 μS/cm	100 MΩ·cm
0.002 μS/cm	500 MΩ·cm	0.02 μS/cm	50 MΩ·cm
0.003 μS/cm	333 MΩ·cm	0.03 μS/cm	33 MΩ·cm
0.004 μS/cm	250 MΩ·cm	0.04 μS/cm	25 MΩ·cm

ORDERING CODES

HD22569.2: The kit is composed of: instrument HD22569.2 for the measurement of pH - redox - conductivity - resistivity - TDS - salinity - concentration of dissolved oxygen, saturation index - temperature, **datalogger**, stabilized power supply at mains voltage 100-240Vac/12Vdc-1A, calibrator HD9709/20 (for polarographic probe) or DO9709/21 (for galvanic probe), instructions manual and software DeltaLog11.

pH/mV electrodes, conductivity probes, dissolved oxygen probes, temperature probes, standard reference solutions for different measurement types, connection cables for pH electrodes with S7 connector, cables for data download to PC or printer have to be ordered separately.

Accessories

9CPRS232: Connection cable SubD female 9- pole for serial output RS232C.

CP22: USB 2.0 connection cable - connector type A - connector type B.

DeltaLog11: Software for download and management of the data on PC using Windows operating systems.

SWD10: Stabilized power supply at 100-240Vac/12Vdc-1A mains voltage.

HD40.1: 24-column portable thermal printer, serial interface, 57mm paper width, four NiMH 1.2V rechargeable batteries, SWD10 power supply, instruction manual, 5 thermal paper rolls. Requires the cable 9CPRS232 (optional).

HD40.2: 24-column portable thermal printer, **Bluetooth and serial interface**, 57mm paper width, four NiMH 1.2V rechargeable batteries, SWD10 power supply, instruction manual, 5 thermal paper rolls. Requires the module HD22BT (optional) or the cable 9CPRS232 (optional).

HD22.2: Laboratory electrode holder composed of basis plate with incorporated magnetic stirrer, staff and replaceable electrode holder. Height max. 380mm. Powered by bench-top meters of the series HD22... with cable HD22.2.1 (optional) or supplier SWD10 (optional).

HD22.3: Laboratory electrode holder with metal basis plate. Flexible electrode holder for free positioning. For Ø 12mm probes.

HD22BT: Bluetooth module for wireless data transmission from instrument to PC. **The fitting of the module into the instrument is made exclusively by Delta Ohm, at the time of placing the order.**

TP47: Connector for Pt100 4-wire and Pt1000 2-wire probes without SICRAM module.

pH electrodes without SICRAM module (Inputs ①)

KP 20: Gel pH combined electrode for general use, with S7 screw connector, EPOXY body.

KP 30: Gel pH combined electrode for general use, 1m cable with BNC, EPOXY body.

KP 50: Gel pH combined electrode, porous Teflon ring junction, suitable for emulsions, demineralised water and waste water with S7 screw connector, glass body.

KP 61: 3 diaphragm liquid filled pH combined electrode for wine, milk, cream, etc., S7 screw connector, liquid reference filling, glass body.

KP 62: 1 diaphragm gel pH combined electrode for general use, pure water, varnishes, gel filled, S7 screw connector, glass body.

KP 63: liquid filled pH combined electrode for general use, varnishes, 1m cable with BNC, glass body.

KP 64: Liquid filled pH combined electrode, Teflon ring diaphragm, for wine, varnishes, emulsions, S7 screw connector, glass body.

KP 70: Pointed gel combined pH microelectrode diam. 6 x L=70 mm., with S7 screw connector, EPOXY body, glass tip, open junction for meat and cheese.

KP 80: Pointed gel pH combined electrode, with S7 screw connector, glass body, for cream, milk, viscous material, open junction.

KP100: Flat membrane gel combined pH electrode with S7 screw connector, glass body, for skin, leather, paper.

Characteristics and dimensions of the probes on page WA-76.

CP: Extension cable 1.5m with BNC connectors on one side and S7 on the other side for electrode with S7 connector.

CP5: Extension cable 5m with BNC connectors on one side and S7 on the other side for electrode with S7 connector.

CP 10: 10m extension cable with BNC/S7 connector for electrode without cable, thread S7.

CP 15: 15m extension cable with BNC/S7 connector for electrode without cable, thread S7.

CE: S7 screw connector for pH electrode.

BNC: Female BNC for cable extension.

pH electrodes with SICRAM module (Input ③)


KP63TS: Combined pH/temperature electrode with SICRAM module, body in glass, 1m cable, 1 diaphragm, for general use, internal liquid reference.

SICRAM Module with S7 input for pH electrodes (Input ③)

pH 471.1: SICRAM module for pH electrodes with S7 standard connection, cable L=1m.

pH 471.2: SICRAM module for pH electrodes with S7 standard connection, cable L=2m.

pH 471.5: SICRAM module for pH electrodes with S7 standard connection, cable L=5m.



ORP Electrodes (Inputs ① and ②)

KP90: Redox Platinum electrode, with screw connector S7, liquid electrolyte, body in glass.

KP91: Redox Platinum electrode with 1m cable with BNC, GEL filled, body in Epoxy.

Electrode dimensions and characteristics at page WA-76

pH buffer solutions

HD8642: Buffer solution 4.01pH - 200cc.

HD8672: Buffer solution 6.86pH - 200cc.

HD8692: Buffer solution 9.18pH - 200cc.

Redox buffer solutions

HDR220: Redox buffer solution 220mV 0,5 l.

HDR468: Redox buffer solution 468mV 0,5 l.

Electrolyte solutions

KCL 3M: Ready for use solution for refilling of the electrodes.

Cleaning and maintenance

HD62PT: Diaphragm cleaning (thiourea in HCl) - 500ml.

HD62PP: Protein cleaning (pepsin in HCl) - 500ml.

HD62RF: Regeneration (fluorhydric acid) - 100ml.

HD62SC: Solution for electrode preservation - 500ml.

Conductivity probes and combined conductivity and temperature probes without SICRAM module (Input ①)

SP06T: Combined conductivity and temperature 4-electrode cell, body in POCAN. Cell constant K=0.7. Measurement range 5μS/cm...200mS/cm, 0...90°C. Max pressure 5bar.

SPT401.001: Combined conductivity and temperature 2-electrode cell in stainless steel AISI 316. Cell constant K=0.01. Cable 2m. Measurement range 0.04μS/cm ...20μS/cm, 0...120°C. Measurement in closed-cell. Max pressure 5bar.

SPT01G: Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant K=0.1. Measurement range 0.1μS/cm ...500μS/cm, 0...80°C. Max pressure 5bar.

SPT1G: Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant K=1. Measurement range 10μS/cm ...10mS/cm, 0...80°C. Max pressure 5bar.

SPT10G: Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant K=10. Measurement range 500μS/cm ...200mS/cm, 0...80°C. Max pressure 5bar.

Probe dimensions and characteristics at page WA-77

Combined conductivity / temperature probes with SICRAM module (Input ©)

SPT16S: Combined conductivity /temperature 2-electrode Platinum- wire cell, body in glass with SICRAM module. Cell constant K = 1. Measuring range 10µS/cm ... 10mS/cm, 0...80°C.

Probe characteristics at page WA-77

Standard calibration solutions

HD8747: Standard calibration solution 0.001mol/l equal to 147µS/cm @25°C - 200cc.

HD8714: Standard calibration solution 0.01mol/l equal to 1413µS/cm @25°C - 200cc.

HD8712: Standard calibration solution 0.1mol/l equal to 12880µS/cm @25°C - 200cc.

HD87111: Standard calibration solution 1mol/l equal to 111800µS/cm @25°C - 200cc.

Combined dissolved Oxygen/temperature probes

DO 9709 SS Polarographic combined oxygen and temperature probe with possibility of membrane replacement. Ø12mm x 120mm. 2m cable. The code includes: probe, 2 membranes, electrolyte solution and zero point solution.

DO 9709 SS.5 Polarographic combined oxygen and temperature probe with possibility of membrane replacement. Ø12mm x 120mm. 5m cable. The code includes: probe, 2 membranes, electrolyte solution and zero point solution.

DO 9709 SS.1 Galvanic combined oxygen and temperature probe with possibility of membrane replacement. Ø12mm x 120mm. Ø16mm tip with membrane. 2m cable. The code includes: probe, 2 membranes in total, electrolyte solution and zero point solution.

DO 9709 SS.5.1 Galvanic combined oxygen and temperature probe with possibility of membrane replacement. Ø12mm x 120mm. Ø16mm tip with membrane. 5m cable. The code includes: probe, 2 membranes in total, electrolyte solution and zero point solution.

Probe dimensions and characteristics at page WA-79

Accessories

DO 9709/20: Calibrator for polarographic probes DO 9709SS and DO 9709SS.5

DO 9709/21: Calibrator for galvanic probes DO 9709SS.1 and DO 9709SS.5.1

DO 9709 SSK: Kit of accessories for probes DO 9709SS and DO 9709SS.5: 3 membranes, zero point solution and electrolyte.

DO 9709/21K: Kit of accessories for probes DO 9709SS.1 and DO 9709SS.5.1: 3 membranes, zero point solution and electrolyte.

DO9700: zero oxygen solution.

DO9701: electrolyte solution for polarographic probes DO9709 SS and DO9709 SS.5.

DO9701.1: electrolyte solution for galvanic probes DO9709 SS.1 and DO9709 SS.5.1.

Temperature probes equipped with SICRAM module

TP472I: Wire wound Pt100 sensor, immersion probe. Stem Ø 3 mm, length 300 mm. Cable length 2 m.

TP472I.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 230 mm. Cable length 2 m.

TP473PI: Wire wound Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP473P.0: Thin film Pt100 sensor, penetration probe. Stem Ø 4mm, length 150 mm. Cable length 2 m.

TP474C.I: Wire wound Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP474C.0: Thin film Pt100 sensor, contact probe. Stem Ø 4mm, length 230mm, contact surface Ø 5mm. Cable length 2 m.

TP475A.0: Thin film Pt100 sensor, air probe. Stem Ø 4mm, length 230mm. Cable length 2 m.

TP472I.5: Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 500 mm. Cable length 2 m.

TP472I.10: Thin film Pt100 sensor, penetration probe. Stem Ø 6mm, length 1000mm. Cable length 2 m.

TP49A.0: Thin film Pt100 sensor, immersion probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AC.0: Thin film Pt100 sensor, contact probe. Stem Ø 1mm, length 150mm. Cable length 2 m. Aluminium handle

TP49AP.0: Thin film Pt100 sensor, penetration probe. Stem Ø 2,7mm, length 150mm. Cable length 2 m. Aluminium handle

TP875.I: Wire wound Pt100 sensor, 150mm diameter globe-thermometer equipped with handle. Cable length 2 m.

TP876.I: Wire wound Pt100 sensor, 50mm diameter globe-thermometer equipped with handle. Cable length 2 m.

TP87.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3 mm, length 70 mm. Cable length 2 m.

TP878.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 2 m.

TP878.1.0: Thin film Pt100 sensor, contact probe for solar panels. Cable length 5 m.

TP879.0: Thin film Pt100 sensor, penetration probe for compost. Stem Ø 8 mm, length 1000 mm. Cable length 2 m.

Temperature probes without SICRAM module

TP47.100.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.











TP47.1000.0: Thin film Pt1000 sensor, immersion probe. Probe's Stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP47: Connector for Pt100 4-wire and Pt1000 2-wire probes without SICRAM module.



TP87.100.0: Thin film Pt100 sensor, immersion probe. Stem Ø 3mm, length 70mm. 4-wires connection cable with connector, length 1 m.

TP87.1000.0: Thin film Pt1000 sensor, immersion probe. Stem Ø 3mm, length 70mm. 2-wires connection cable with connector, length 1 m.



pH electrodes without SICRAM module for HD2205.2, HD2256.2, HD2259.2, HD22569.2 ① ②, HD2305.0, HD2105.X, HD3405.2, HD3456.2

ORDERING CODE	MEASUREMENT RANGE AND USE	DIMENSIONS
KP 20	0...14pH / 0...80°C / 2bar Body in Epoxy - GEL filled 1 ceramic diaphragm Waste water, drinking water, paints, water emulsions, galvanic baths, fruit juices, water suspensions, titration, varnishes.	
KP 30	0...14pH / 0...80°C / 2bar Body in Epoxy - GEL filled 1 ceramic diaphragm Cable L=1m with BNC Waste water, drinking water, water emulsions, galvanic baths, paints, varnishes, water suspensions, fruit juices, titration.	
KP 50	0...14pH / -5...+100°C / 2bar Body in glass - GEL filled 1 Teflon ring diaphragm Varnishes, cosmetics, water emulsions, galvanic baths, creams, deionised water, TRIS solutions, drinking water, fruit juices, low-ion-content solutions, mayonnaise, preserved food, paints, titration, titration in non-water solutions, water suspensions, detergents, waste water, viscous samples.	
KP 61	2...14pH / 0...80°C / 2bar Body in glass Liquid reference filling Triple ceramic diaphragm Waste water, paste, bread, fruit juices, varnishes, cosmetics, creams, deionised water, drinking water, water emulsions, galvanic baths, detergents, yoghurt, milk, titration, preserved food, titration in non-water solutions, water suspensions, mayonnaise, wine, low ion-content solution, butter, proteic substances, paints, viscous samples	
KP 62	0...14pH / 0...80°C / 2bar Body in glass - gel filled 1 ceramic diaphragm Paints, varnishes, drinking water, water emulsions, fruit juices, galvanic baths, water suspensions, titration, waste water.	
KP 63	0...14pH / 0...80°C / 1bar Body in glass Reference filling solution KCl 3M 1 ceramic diaphragm Cable L=1m with BNC Paints, varnishes, drinking water, water solutions, fruit juices, galvanic baths, water suspensions, titrations, waste water.	
KP 64	0...14pH / 0...80°C / 0.1bar Body in glass Liquid reference KCl 3M Teflon collar diaphragm Paints, varnishes, cosmetics, creams, deionised water, drinking water, water emulsions, fruit juices, detergents, low ion-content solutions, preserved food, water suspensions, titration, titration in non-water solutions, TRIS solutions, waste water, viscous samples, wine.	
KP 70	2...14pH / 0...50°C / 0.1bar Body in Epoxy - gel filled 1 open junction Paste, bread, paints, varnishes, cosmetics, creams, drinking water, water emulsions, fruit juices, galvanic baths, detergents, mayonnaise, preserved foods, cheese, milk, water suspensions, viscous samples, waste water, butter, yoghurt.	
KP 80	2...14pH / 0...60°C / 1bar Body in glass - gel filled 1 open junction Paste, bread, paints, varnishes, cosmetics, creams, drinking water, water emulsions, fruit juices, galvanic baths, detergents, mayonnaise, preserved food, water suspensions, titration, titration in non-water solutions, viscous samples, waste water, yoghurt, milk, butter.	
KP 100	2...14pH / 0...80°C / 1bar Body in glass Liquid reference KCl 3M Teflon ring diaphragm Flat membrane gel combined pH electrode, S7 connector, for skin, leather, paper.	




pH electrodes complete with SICRAM module for HD2205.2, HD2256.2, HD2259.2, HD22569.2 ③ ④ and HD98569

ORDERING CODE	MEASUREMENT RANGE AND USE	DIMENSIONS
KP 63TS	0...14pH / 0...80°C / 1bar Pt100 sensor Body in glass Reference filling solution KCl 3M 1 ceramic diaphragm Cable L=1m with SICRAM module Paints, varnishes, drinking water, water solutions, fruit juices, galvanic baths, water suspensions, titrations, waste water.	
pH 471.1 pH 471.2 pH 471.5	SICRAM module for pH electrodes with S7 connector. Cable L= 1, 2 or 5m. See characteristics of the connected electrode. Stores the calibration data of the connected Electrode.	



Redox Electrodes for HD2205.2, HD2256.2, HD2259.2, HD22569.2 ① ②, HD2305.0, HD2105.X, HD2156.X, HD3405.2, HD3456.2

ORDERING CODE	MEASUREMENT RANGE AND USE	DIMENSIONS
KP 90	±2000mV 0...80°C 5bar Body in glass Reference filling solution KCl 3M General use	
KP 91	±1000mV 0...60°C 1bar Body in Epoxy - GEL Cable L=1m with BNC General use No heavy tasks	






Combined 2-ring or 4-ring conductivity probes without SICRAM module for HD2206.2, HD2256.2, and HD22569.2 ⑦, HD2306.0, HD2106.X, HD2156.X, HD3406.2, HD3456.2

ORDERING CODE	MEASUREMENT RANGE AND USE	DIMENSIONS
SP06T	K=0.7 5µS/cm ...200mS/cm 0...90°C 4-electrode cell in Pocan/Platinum Probe material Pocan General use No heavy tasks Max. pressure 5bar	
SPT401.001 Not suitable for HD2306.0	K=0.01 0.04µS/cm ...20µS/cm 0...120°C 2-electrode cell in AISI 316 Ultrapure water Measurement in closed-cell Max. pressure 5bar	
SPT01G	K=0.1 0.1µS/cm ...500µS/cm 0...80°C 2-electrode cell in Platinum wire Probe material glass Pure water Max. pressure 5bar	

2-ring or 4-ring conductivity probes without SICRAM module for HD2206.2, HD2256.2 and HD22569.2 ⑦, HD2306.0, HD2106.X, HD2156.X, HD3406.2, HD3456.2

ORDERING CODE	MEASUREMENT RANGE AND USE	DIMENSIONS
SPT1G	K=1 10µS/cm ...10mS/cm 0...80°C 2-electrode cell in Platinum wire Probe material glass General heavy tasks, average conductivity Max. pressure 5bar	
SPT10G	K=10 500µS/cm ...200mS/cm 0...80°C 2-electrode cell in Platinum wire Probe material glass General heavy tasks, high conductivity Max. pressure 5bar	

Conductivity probes with SICRAM module for HD2206.2, HD2256.2, HD22569.2 ⑦ and HD98569


ORDERING CODE	MEASUREMENT RANGE AND USE	DIMENSIONS
SP06TS	K=0.7 5µS...200mS/cm 0...90°C 4-electrode cell in Pocan/Platinum General use No heavy tasks Max pressure 5bar	
SPT401.001S Not suitable for HD2306.0	K=0.01 0,04...20µS/cm 0...120°C 2-electrode cell AISI 316 - Teflon Ultrapure water Measurement in closed-cell Max pressure 5bar	
SPT01GS	K=0.1 0.1µS...500µS/cm 0...80°C 2-electrode cell in Glass/Platinum Pure water Max pressure 5bar	
SPT1GS	K=1 10µS...10mS/cm 0...80°C 2-electrode cell in Glass/Platinum General heavy tasks, average conductivity Max pressure 5bar	
SPT10GS	K=10 500µS...200mS/cm 0...80°C 2-electrode cell in Glass/Platinum General heavy tasks, high conductivity Max pressure 5bar	

Dissolved oxygen probe for HD2259.2, HD22569.2 ©, HD2109.X, HD3409.2


Model	D09709 SS	D09709 SS.5	D09709 SS.1	D09709 SS.5.1
Type	Polarographic probe, Silver anode, Platinum cathode		Galvanic probe, Zinc anode, Silver cathode	
O ₂ measuring range	0.00...60.00mg/l		0.00...20.00mg/l	
Functioning temperature	0...45°C		0...50°C	
Accuracy instrument with probe	±1% FS		±2% FS.	
Membrane	Replaceable		Replaceable	
Cable length	2m	5m(*)	2m	5m(*)

(*) Cable with connector


D09709 SS




D09709 SS.5



D09709 SS.1



D09709 SS.5.1



TECHNICAL DATA OF PROBES AND MODULES EQUIPPED WITH INSTRUMENT
Temperature probes Pt100 sensor with SICRAM module

Model	Type	Application field	Accuracy
TP472I	Immersion	-196°C...+500°C	±0.25°C (-196°C...+300°C) ±0.5°C (+300°C...+500°C)
TP472I.0 1/3 DIN Thin Film	Immersion	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP473P.I	Penetration	-50°C...+400°C	±0.25°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP473P.0 1/3 DIN Thin Film	Penetration	-50°C...+300°C	±0.25°C (-50°C...+300°C)
TP474C.I	Contact	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.5°C (+300°C...+400°C)
TP474C.0 1/3 DIN Thin Film	Contact	-50°C...+300°C	±0.3°C (-50°C...+300°C)
TP475A.0 1/3 DIN Thin Film	Air	-50°C...+250°C	±0.3°C (-50°C...+250°C)
TP472I.5	Penetration	-50°C...+400°C	±0.3°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP472I.10	Penetration	-50°C...+400°C	±0.30°C (-50°C...+300°C) ±0.6°C (+300°C...+400°C)
TP49A.0 Class A Thin Film	Immersion	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AC.0 Class A Thin Film	Contact	-70°C...+250°C	±0.3°C (-50°C...-50°C) ±0.25°C (-50°C...+250°C)
TP49AP.0 Class A Thin Film	Penetration	-70°C...+250°C	±0.3°C (-70°C...-50°C) ±0.25°C (-50°C...+250°C)
TP875.I	Globe-thermometer Ø150mm	-30°C...+120°C	±0.25°C
TP876.I	Globe-thermometer Ø50mm	-30°C...+120°C	±0.25°C
TP87.0 1/3 DIN Thin Film	Immersion	-50°C...+200°C	±0.25°C
TP878.0 1/3 DIN Thin Film TP878.1.0 1/3 DIN Thin Film	Photovoltaic	+4°C...+85°C	±0.25°C
TP879.0 1/3 DIN Thin Film	Compost	-20°C...+120°C	±0.25°C

Common characteristics

Temperature drift @ 20°C 0.003%/°C

4 wires Pt100 and 2 wires Pt1000 Probes

Model	Type	Application field	Accuracy
TP47.100.0 1/3 DIN Thin Film	4 wires Pt100	-50...+250°C	1/3 DIN
TP47.1000.0 1/3 DIN Thin Film	2 wires Pt1000	-50...+250°C	1/3 DIN
TP87.100.0 1/3 DIN Thin Film	4 wires Pt100	-50...+200°C	1/3 DIN
TP87.1000.0 1/3 DIN Thin Film	2 wires Pt1000	-50...+200°C	1/3 DIN

Common features

Temperature drift @20°C

Pt100 0.003%/°C

Pt1000 0.005%/°C



HD 22.2 - HD 22.3 ELECTRODE HOLDER

HD 22.2: Laboratory electrode holder composed of base plate with built-in magnetic stirrer, shaft and replaceable electrode holder. Suitable diameter 12mm. Powered by bench top instruments of the series **HD22...** with cable HD22.2.1 (**optional**), or with power supplier SWD10 (**optional**).

HD 22.3: Laboratory electrode holder with base plate. Flexible arm for free positioning. Suitable for electrodes with diameter 12mm.



HD22.2

HD22.3



HD 40.1, HD 40.2 PORTABLE THERMAL PRINTER

The **HD40.1** and **HD40.2** are lightweight, compact, portable thermal printers.

The **HD40.1** is connected to instruments or PC through the **RS232** serial input. The **HD40.2** features a dual mode data reception system - **RS232** serial and **Bluetooth**.

The Bluetooth wireless connection makes the HD40.2 printer very useful "in the field", since it does not require any connection to the instrument. A careful design allows you to replace the thermal paper roll in a few seconds. A four NiMH rechargeable battery pack provides power supply and ensures long autonomy: you can print up to 3000 lines at full charge.

Standard thermal paper roll width: 57mm.

Print resolution: 203 dpi

Print characters (each line): 24

Protection degree: IP40.

Specifications

Printing method	Thermal
Resolution	203 DPI (8 dot/mm)
Printing width	48mm centered in the paper roll
Paper roll width	57mm ... 58mm
Max. paper roll diameter	32mm
Number of columns	24
Printing speed	Up to 90 mm/sec (depending on battery charge and ambient conditions)

Sensors	Paper detection
Character set	IBM II 858 table
Printing formats	Normal or extended
Character font	1 (16 x 24 dot – 2mm x 3mm)

Thermal head durability	100 million pulses (temperature: 20...25°C)
Mechanism life	50km of paper (temperature: 20...25°C)
Abrasion resistance	2000 opening/closing cycles or more
Cover group durability	RS232
Communication interfaces	Bluetooth (for HD40.2)
RS232 Baud rate	9600, 19200 and 38400 baud (the factory parameter is 38400 baud)
Bluetooth Baud rate	38400 baud (for HD40.2)
Bluetooth operating distance	Up to 10m without hindrance (for HD40.2)

Mains power supply (cod. SWD10)	100-240Vac/12Vdc-1A mains battery charger
Batteries	Four 1.2V AA rechargeable batteries (NiMH)
Printing autonomy	3000 lines 24 characters each printing, one line every 10 seconds
Switch-off function	0, 5, 10 or 15 minutes - excludable.
Dimensions	105mm x 165mm x 53mm
Weight	380g (with batteries and paper roll)
Material	ABS

Operating conditions

Operating temperature	0°C ... 50°C
Operating relative humidity	20%RH ... 85%RH not condensing
Storage Temperature / Relative humidity	-25°C ... +70°C / 10%RH ... 90%RH not condensing
Protection degree	IP40

Connections

Serial interface	9-pole D sub male connector
Battery charger power supply (cod. SWD10)	2-pole connector (positive in the centre)

ORDERING CODES

HD40.1: The kit includes: 24-column portable thermal printer, **serial interface RS232**, 57mm paper width, four NiMH 1.2V rechargeable batteries, SWD10 power supply, instructions manual, 5 thermal paper rolls.

HD40.2: The kit includes: 24-column portable thermal printer, **Bluetooth and serial interface RS232**, 57mm paper width, four NiMH 1.2V rechargeable batteries, SWD10 power supply, instructions manual, 5 thermal paper rolls.

The serial cable for PC/instrument connection must be ordered separately.

HD2110CSNM: RS232C 8-pole MiniDin - 9-pole D Sub female null-modem cable for connecting the printer to instruments with MiniDIN connector (HD21xx.1 and HD21xx.2 series, HD34xx.2, HD98569, etc.).

9CPRS232: RS232C 9-pole D Sub female null-modem cable for connecting the printer to instrument with 9-pole D Sub connectors (Delta Ohm instruments: HD22xx.2 series, HD25.2, etc.).

HD 2110RS: M12 - 9-pole D Sub connectors cable for connecting the printer to instruments with M12 connector (Delta Ohm instruments: HD 2010UC, HD 2010UC/A, HD 2110L).

SWD10: 100-240Vac/12Vdc-1A Mains battery charger.

BAT40: Spare battery pack for HD40.1 and HD40.2 printers with in-built temperature sensor.

RCT: The kit includes 4 thermal paper rolls 57mm wide and 32mm diameter.





HD 25.2 BENCH-TOP TURBIDITY METER

The **HD25.2** is a digital turbidity meter for laboratory and mobile use, suitable for measurements in drinking water, waste water and process liquids. The working principle is based on the nephelometric (90° scattered light sensor) method.

It is equipped with three light detectors and two LED light sources (white and infrared) which are permanently kept under control in order to guarantee long-term stability. The instrument performs measurements according to the standards **EPA 180.1, ISO-NEPH (ISO 7027), EBC and ASBC**. It is also able to carry out measures of transmission factor percentage of white and infrared light.

The initial factory calibration is based on Formazin primary standard. For routine calibration a set of stabilized secondary standard solutions is available: **STCAL** (Turbidity standards for calibration):

- **STCAL 1** less than 0,05 NTU
- **STCAL 2** equal to 8 NTU
- **STCAL 3** equal to 80 NTU
- **STCAL 4** equal to 800 NTU



User Calibration is automatic on one or four points, depending on the measuring variable. Stabilized power supply and advanced electronics guarantees optimal performances over time.

The HD25.2 is a **datalogger** that stores up to 999 samples.

The data can be transferred from the instrument connected to a PC via the RS232C and USB 2.0 serial ports.

The RS232C serial port can be used to transfer the acquired measurements to a 24 column printer (i.e. HD40.1)

The Print function allows to print labels with progressive and automatically incrementing numeration, with all data related to the sample being examined. The dedicated software **DeltaLog11** allows instrument management and data processing on PC.

The use of the HD25.2 by more users is facilitated by the "User Management" function, which allows, according to the case, to enable or disable some advanced functions of the instrument through password.

The protection degree is IP66.

Technical characteristics

Instrument

Dimensions (Length x Width x Height)	220x120x55mm
Weight	400g (batteries included)
Materials	ABS, rubber
Display LCD	4½ characters plus symbols
	Visible area: 52x42mm

Operating conditions

Instrument working temperature	0 ... 50°C
Storage temperature instrument	-25 ... 65°C
Working relative humidity	0 ... 90% R.H. without condensation
Storing of Calibration standards	5...25°C (temperature should not exceed, protect from light)

Protection degree

IP66

Power supply

Batteries	3 1,5 V AA type batteries
Autonomy	100 hours with 1800mAh alkaline
Mains	Mains adapter (cod. SWD10) 100-240Vac/12Vdc-1A

Measuring methods

Standard	EPA180.1, ISO-NEPH (ISO 7027), EBC, ASBC, WHITE %T e IR %T
Light source	LED IR (850nm) and white LED (470nm)
Receiver	Silicium photodiode
Sample cell	Ø24mm - height 68mm, 20cc



Measurement of turbidity

Method / Measuring range

EPA180.1	(0...1000 NTU)
ISO-NEPH	(0...1000 FNU)
EBC	(0...250 EBC)
ASBC	(0...9999 ASBC)
WHITE %T	(0...100 %T)
IR %T	(0...100 %T)
0.01 NTU	(0...9.99 NTU)
0.1 NTU	(10.0...99.9 NTU)
1 NTU	(100...1000 NTU)
Accuracy	±2% reading + 0.01 NTU (0...500 NTU)
	±3% reading (500...1000 NTU)
Repeatability	±2% reading or 0.01 NTU (the major one)

Security of memorized data

Unlimited

Time

Date and hour

real time schedule

Accuracy

1min/month max error

Measured values storing

Quantity

999 samples

Serial interface RS232C

Type

RS232C electrically isolated

Baud rate

Can be set from 1200 to 38400 baud

Data bit

8

Parity

None

Stop bit

1

Flow Control

Xon/Xoff

Serial cable length

Max 15m

USB interface

Type

1.1 - 2.0 electrically isolated

Connections

Serial interface

DB9 connector (9- pole male)

USB interface

USB connector type B

Mains adapter

2- pole connector (Ø5.5mm-2.1mm). Positive at centre.

Ordering codes

HD 25.2: The kit is composed of: instrument HD25.2, 4 empty cells, 4 calibration standards STCAL, 3 1.5Vdc alkaline batteries, lubricant rag, 25cc Silicon oil, instructions manual, carrying case and software DeltaLog11 for PCs running Windows operating systems.

Accessories

9CPRS232: Connection cable SubD female 9- pole for serial output RS232C

CP22: Connection cable USB 2.0 connector type A - type B

SWD10: Stabilized power supply at 230Vac/9Vdc-300mA mains voltage.

HD40.1: 24-column portable thermal printer, serial interface, 57mm paper width, four NiMH 1.2V rechargeable batteries, SWD10 power supply, instruction manual, 5 thermal paper rolls. Requires the cable 9CPRS232 (optional).

PL: Lubricant rag

OS1: Silicon oil - 25cc.

KCV: 4 empty sample cells Ø24x68mm

Turbidity calibration standards

STCAL 1: Calibration standard with low turbidity Formazin reference less than 0,05 NTU.

STCAL 2: Calibration standard with Formazin reference 8 NTU - 20cc.

STCAL 3: Calibration standard with Formazin reference 80 NTU - 20cc.

STCAL 4: Calibration standard with Formazin reference 800 NTU - 20cc.

KS: Kit 4 calibration standard with Formazin reference STCAL 1, STCAL 2, STCAL 3, STCAL 4.



HD40.1



ASBC



FNU



EBC



NTU